MAPS

Instruction Manual

Manual No.: 14597001

Date of Release 10.12.2009

MAPS
Instruction Manual

Table of Contents	Page
User's Guide	3
Quick Reference Guide	99

Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the *Date* and *Article No* of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2009.

Struers A/S
Pederstrupvej 84
DK-2750 Ballerup
Denmark
Telephone +45 44 600 800
Fax +45 44 600 801

MAPS Safety Precaution Sheet

To be read carefully before use

1. The operator should be fully aware of the use of the module according to the Instruction Manual.
2. The module must be installed in compliance with local safety regulations.
3. The actual voltage must correspond to the voltage stated on the side of the module. The module must be earthed.
4. The water connections must be leakproof. Cut the water off if the module stands idle for a long period.
5. Be careful that the emergency stop is in order.
6. Always clamp the specimens securely in the specimen holder disc before starting the preparation process.
7. Never touch any moving parts of MAPS, i.e. specimen holder discs, arm or preparation discs, until they stop completely.
8. If any malfunction or unusual noise is detected, the apparatus should be stopped immediately and technical service called.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.)

Table of contents	Page
1. Introduction	
MAPS - Automated Preparation	5
Development or production?	6
MAPS Module design	6
Work stations	6
Cleaning station	8
Conveyor belts	8
Single Unit without Conveyors	8
Preparation methods	8
Method steps	8
Step parameters	9
Preparation Module Software	9
Production line	9
Menu structure	10
2. Operation	
Controls	11
Main controls	12
PC keys/function keys	12
Struers keys	13
Dosing controls	13
Security system	13
Numerical keys	13
Manual predosing	14
Software	15
Start screens	15
Main menu	16
Process mode	17
F2 P-CTRL, Process control	20
Development mode	22
F1 P-EDIT (development mode)	25
Option menus	26
F6 M-FUNC, Manual functions	28
Method database	30
F1 EDIT (Method database)	32
F3 COPY	33
F7 DELETE	34
Consumables	35
Refill/dispose liquids	35
Changing preparation discs	38
Changing grinding stone	39

3. Preparation of specimens

Preparation methods	41
Struers Metalog Guide™	41

4. Maintenance

Mechanical maintenance	42
Daily	42
Every evening or before breaks of more than 8 hours	43
Every second day	44
Every week or before breaks of more than 24 hours	44
Every month	44
Software maintenance	46
Configuration menu	46
Discs/Cloths, configuration	47
Lubricants, configuration	48
Suspensions, configuration	49
Cleaning programs, configuration	50
Module configuration	51

5. Accessories

Preparation modules	53
Conveyor belt	53
Water recirculation	53
Liquids supply	54
Specimen handling	54
Specimen holder discs	55
Preparation stones and discs	55
Software	56
Tools	56

6. Installation

External hose connections	57
Connection tube collar for optional exhaustion	57
MAPS data communication and control	58

7. Trouble-shooting

Error code listing	59
Function Errors (Fatal errors)	59


8. Spare Parts	69
9. Technical data	
Preparation module	70
Electrical data	70
Connections	71
Main data	72
Specimen holder discs	72
Transportation arm	72
Plane grinding work station	72
Fine grinding and polishing work station	72
Cleaning station	73
Built-in computer facilities	73
External computer facilities (accessory)	73
Method programming	73
Conveyor belt	74
Belt connections	74
Belt main dimensions	74
General data, Conveyor belt	74
Belt capacity	74
Recirculation cooling unit	74
Cooling unit capacity	74
Dosing system	74
Bottle unit capacity	74
Alcohol separation unit capacity	74
MAPS Master, PC kit	75
Code numbering of MAPS modules	75
Work station codes	75
Dosing system/Bottle unit codes	75
Bottle unit numbering	76
Construction descriptions	77
Specimen holder disc	78
Transportation arm	79
Plane grinding station	79
Dresser for stone	80
Fine grinding/polishing station	80
Cleaning station	81
Recirculation cooling Unit	82
Alcohol separation unit	83
Dosing system	83
Dosing system level indicator	83
Conveyor belt	84
Drying fan	85
Safety	85

10. Method parameter values

Method parameters	86
Method number in the module database	86
Method description.....	86
Method number in the MAPS Master PC database	86
Step parameters.....	88
Work position	88
Grinding/polishing disc.....	88
Grain/grit size of disc	88
Suspension	88
Lubricant	89
Rotational speed of the disc.....	89
Rotational direction of the specimen holder	89
Cleaning program number	90
Dressing depth.....	90
Phase parameters.....	91
Process time	91
Stock removal	91
Force.....	91
Dosing level for suspension (Susp) and lubricant (Lub.)	92
Positioning the transportation arm during the process.....	94
Cleaning programs.....	95
Disc cooling.....	96
11. Index.....	97

1. Introduction

MAPS - Automated Preparation MAPS preparation module is not only powerful and efficient - it is easy to use. All the important functions are available on the front panel keyboard equipped with modern touch buttons and LC display. Apart from feeding the preparation module with specimen holders and changing consumables it is seldom necessary to handle the specimens manually.

- ① Front panel with LC display and PC keys (standard)
- ② Conveyor delivery module with drying fan (optional)
- ③ Conveyor feed module (optional)
- ④ Dosing systems with bottle units (optional)
- ⑤ Recirculation cooling system and separation unit (optional)

Development or production?

You may use MAPS in two ways - or in a combination of both:

Development mode

Before you can run MAPS as an automated system you need to have a set of thoroughly tested preparation methods. You design and test the methods in the DEVELOPMENT MODE (see Section 2, Operation).

Process mode

To run the preparation methods you need a queue system capable of handling a large number of different materials. This is achieved in the PROCESS MODE (see Section 2, Operation).

Which mode to use largely depends on the type of preparation work you are doing. A great number of identical preparations may be carried out in the PROCESS MODE, whereas experiments or test preparations are more likely to be done in DEVELOPMENT MODE.

MAPS Module design

Work stations

A MAPS preparation module is either equipped with 1 plane grinding and 1 fine grinding/polishing work station (MAPSE/MAPON) or 2 fine grinding/polishing work stations (MAPFI/MAPVA). There are 2 different basic work station types:


Stone/Diamond Grinding: a heavy duty work station equipped with dresser and separate safety lid. Version available with grinding stone (PLAN/POL) or diamond disc (DIA/POL). Used for plane grinding (**PG**).

Fine Grinding/Polishing: an open work station with room for many types of easily exchangeable discs (POL/POL). Used for fine grinding (**FG**), diamond polishing (**DP**) or oxide polishing (**OP**).

All the above types of work station can be equipped with an optional Recirculation cooling unit (MAPRE): see Section 5, Accessories.

Basic MAPS module configurations			
Type	1st Work station	2nd Work station	Code word
PLAN/POL	Plane grinding	Fine grinding/polishing	MAPSE/ MAPON
DIA/POL	Plane grinding	Fine grinding/polishing	
POL/POL	Fine grinding/polishing	Fine grinding/polishing	MAPFI/ MAPVA

Full MAPS configuration with 5 modules (example)

Cleaning station

All MAPS module configurations have a programmable built-in cleaning station. The automatic cleaning stations work by high-pressure water rinsing, alcohol cleaning and air-drying. The cleaning station is placed between the two work stations. All cleaning of specimen holder discs take place in the cleaning station.

Conveyor belts

If you have more than one MAPS preparation module you need a number of conveyor belts to transport the specimen holder discs between the preparation modules. The belts work automatically, communicating with each other and the MAPS modules through optical sensors.

Single Unit without Conveyors

The special version of MAPS for use as a single unit without conveyors has the following specifications which are different from the standard specifications of a MAPS preparation module:

- Instead of conveyor belts you have a combined feed and delivery station placed opposite to the front control panel side.
- The single unit without conveyors is provided with 2 extra safety guards compared to the standard MAPS module. The 2 extra safety guards are placed in the back of the machine, facilitating the access to the feed and delivery station.

Preparation methods

Method steps

A preparation method is divided into a number of steps, usually between 5 and 10 steps (See Section 3, Preparation according to Struers Metalog Guide™). Each step belongs to one of the following types:

PG = Plane Grinding

For plane grinding (PG) you need a MAPSE or MAPON preparation module equipped with a grinding stone or a Diamond Pad. As the MAPSE/MAPON preparation module is also equipped with a fine grinding/polishing work station, the module may be used as a single unit without conveyors.

FG = Fine Grinding

For fine grinding (FG) you can use both a MAPSE/MAPON preparation module (one fine grinding/polishing work station) or a MAPFI/MAPVA preparation module (two fine grinding/polishing work stations).

DP = Diamond Polishing

You may use all types of MAPS preparation modules for both diamond polishing (DP) or in case of soft and ductile materials:

OP = Oxide Polishing

Oxide polishing (OP) or diamond polishing (DP) are done on a fine grinding/polishing work station.

Step parameters

The MAPS preparation module database enables you to fine-tune all details of the method, including the parameters of the individual method steps (see also Section 10, Method parameter values).

PG: stone/Diamond Pad including grain size, stock removal depth or timer- based grinding, force on specimens, lubrication, dressing depth (only for stone) and disc speed/rotational direction.

FG: grinding disc including grain size, timer, force on specimens, lubrication and disc speed/rotational direction.

DP: polishing disc, timer, force on specimens, suspension type, lubricant and disc speed/rotational direction.

OP: as above but without lubricant.

Cleaning: all step types have a choice of 10 user defined cleaning programs.

Preparation Module Software

The MAPS preparation modules are each equipped with a built-in preparation method database and editor. The MAPS preparation module database contains 30 methods in all, each with 1-10 steps. The database is designed to work both within the preparation module or connected to the optional MAPS Master PC database program.

In the method description in the database, you appoint each grinding/polishing step to a *position*, i.e. a disc on a specific module. When the conveyor belt delivers the specimen holder, the receiving module automatically "knows" which work station to use. If the method does not use any of the module's work stations (position numbers), the specimen holder disc is transported to the next preparation module.

To edit the methods, MAPS provides you with a *Method editor*. The editor allows you to change the data in a method and send or receive a method between the Method editor and the optional MAPS Master PC program.

Production line


MAPS is a production line system with a number of grinding and polishing discs. Each module may be identified by a specific *module number* which includes 2 work station identification numbers, called "*position numbers*". See also Section 10, Method parameter values.

On the LC display you can monitor the preparation process. You can even manipulate the methods during the process if you work in Development mode.

Menu structure


The MAPS preparation module software supports processing, monitoring and editing of preparation methods. The program is based upon a database with room for detailed information on up to 30 methods.

You can use the MAPS front panel almost like a PC, including easy LC display menus and well-known computer keys like ENTER, ESC and PGDN/PGUP. The menu structure reflects the main functions of MAPS; Process mode, Development mode, Database editing and configuration of consumables.


2. Operation

Controls


- ① PC keys/function keys
- ② Struers keys
- ③ Numerical keys

Main controls

Main switch: blue/grey switch placed on the front panel below the keyboard.

Emergency stop: red push button placed on the front below the keyboard (on certain models two emergency buttons are placed on the left and right panels of the module).

PC keys/function keys

PC keys

ENTER accepts a parameter or shows a pick list of consumables to choose from ("Option menu"). In the option menus ENTER is used to accept the option parameter.

ESC (escape) is used to return to the previous screen.

DEL (delete) deletes the previous digit, like the backspace key on a PC.

ARROW keys are used for selecting the desired items. Other uses of the ARROW keys are explained in the function tables for each screen.

**PGDN/
PGUP** In some of the option menus, you can scroll up and down through several pages of information. In the Method edit screen, the PGDN/PGUP keys are used for moving between method steps.

Function keys

F1 to F8 are used to activate various functions and to change screens. The function keys are explained in a table below each screen.

Struers keys

Start: re-starts the selected method.

Stop: stops any active method or mechanical function.

Continue: continues a method process that has been temporarily stopped.

Dosing controls

Predose: starts the polishing disc rotating. Used while predosing on the polishing cloth.

Lubricant: forced manual dosing of lubricant on the chosen polishing disc.

Abrasive: forced manual dosing of abrasive on the chosen grinding disc.

Security system

Lock MAPS: locks the MAPS system to prevent unauthorized access to the module. A special Security Manual Sheet has been delivered with the MAPS module. If you need an extra copy, please call your local Struers dealer.

Numerical keys

To type method number and other parameters (data), press the desired number on the keypad and then press ENTER to accept the parameter.

Manual predosing

You can perform manually controlled predosing by using the following Struers keys:

Predose: rotate the disc selected by the actual method step. The key does not perform the actual predosing by itself.

Lubricant: dose with lubricant, if available in the actual method step.

Abrasive: dose with abrasive, if available in the actual method step.

Disc position and dosing type depend on the parameters of the present method step.

Software

Start screens

When you turn on the MAPS preparation module with the main switch, a short message: "MAPS CONSOLE READY" and the MAPS screen appear. Press ENTER to activate the program. If you have more than one preparation module you must start each module individually.

Start screen for PLAN/POL (Stone plane grinding/polishing) module

```
          ///  Struers MAPS
=====

          Version: 1.08.00
          Grind/Polish

          Press ENTER to continue
```


Start screen for POL/POL (fine grinding/polishing) module


```
          ///  Struers MAPS
=====

          Version: 1.08.00
          Polish/Polish

          Press ENTER to continue
```

Start screen for DIA/POL (Diamond plane grinding/polishing) module

Main menu

Having pressed ENTER, the MAIN MENU screen appears.

Keys	Explanation
ARROW UP/DOWN	Select a menu
ENTER	Activate a menu

Process mode

Each MAPS preparation module is equipped with a "process queue" function, enabling you to process materialographic samples as batch jobs with a number of batches in the range of 1 to 30. A *batch* may consist of a *quantity* of specimen holder discs in the range of 1 to 99. The optional MAPS Master program allows you to supervise and control the process queue. In the case of 2 or more MAPS preparation modules working together, we strongly recommend that you operate from a PC with the MAPS Master program.

The special single unit without conveyors cannot, however, use the Process mode as the module is not equipped with conveyor belts.

The most straightforward way to prepare multiple specimens is by setting up a batch program for the specimen holder discs.

Programming a batch queue

A batch queue is very simple to program in PROCESS MODE, almost like writing a shopping list:

- Type the method number and quantity of each batch. If you need to change the specimen holder size this can be done in the list too.
- Press \diamond to start the actual process.

Editing a batch queue

You can insert a batch (item) between two existing batches (F3=INSERT) or you can remove a batch (F7=REMOVE). To delete the total batch queue program press F5=CLEAR.

Process queue screen

PROCESS QUEUE				
ITEM	METHOD	DESCRIPTION	SIZE	QUA
1	0		0	0
2	0		0	0
3	0		0	0
4	0		0	0
↓ 5	0		0	0

Type new parameter value

F1: F2:P-CTRL F3:INSERT F4:
F5: CLEAR F6: F7: REMOVE F8:

Keys	Explanation
ARROW UP/DOWN	Choose method number, specimen holder disc diameter size and number of specimen holder discs in the batch. The arrows scroll the screen while still keeping the first batch line in view (4 batch lines scroll).
◊	Start the MAPS preparation module and proceed to the Process control screen.
⊖	Stop the MAPS preparation module.
F2 P-CTRL	Move direct to the Process control screen.
F3 INSERT	Insert an empty batch line in the process queue above the marked line.
F5 CLEAR	Clear all batch lines in the queue.
F7 REMOVE	Remove a batch line (empty or with a method) in the process queue.

Method number format

Methods can be transmitted to the process queue in two ways: either externally from the MAPS Master PC program or internally by typing the method number on the PC keys of the front panel.

Method number from the PC: format [X (main group no.), Y (sub group no.), Z (method no.)]. Ex. (1,1,1). The numbers refer to the MAPS Master database in the PC.

Method number in the preparation module: a number between 1 and 30. The number refers to the internal database of the MAPS module.

F2 P-CTRL

View the contents of the active method (batch line 1).

F5 CLEAR

F5 clears all batch lines in the queue screen. This may take a little while as MAPS resets all parameters of the queue system.

F3/F7 INSERT/REMOVE

F3=INSERT and F7=REMOVE change the order of the 30 possible batches in the process queue. When inserting a new batch job (method), you delete the last batch line (number 30). When removing a batch line, you move all the batch lines in the queue, below the removed line, 1 line up.

IMPORTANT

Take care when using INSERT and REMOVE

F2 P-CTRL, Process control

You can view the contents of the active method (batch line 1) by pressing F2=P-CTRL. In the Process Control screen, you can omit one or more steps by pressing F8=OMIT (all steps less than the marked step are omitted) or vice versa by pressing F4=REPEAT. By pressing F6=M-FUNC (Manual functions) you can perform a number of mechanical functions.

```

Method:  1*                               PROCESS CONTROL
Description:
Holder size: 200   Remaining time: 0:00

STEP  POS  DISC      SUSP/#  LUB      TIME/μm
1    0
2     0              0:00
                                0:00

Use cursor keys to select parameter

F1:           F2:           F3:           F4: REPEAT
F5:           F6: M-FUNC  F7:           F8: OMIT
    
```

Keys	Explanation
	Starts the MAPS preparation module with the selected method. The first executable step will be the first step to start.
	Stops the MAPS preparation module.
	Continues processing of the active method. The selected step will be the first step to start.
F4 REPEAT	Activates the previous step in the method
F6 M-FUNC	Activates the Manual functions screen
F8 OMIT	Omits the active step and moves onto the next step

Method number format*)

- If the active method (batch line 1) has been transmitted from the MAPS Master PC program, the method no. will be of the format [X (main group no.), Y (sub group no.), Z (method no.)]. Ex. (1,1,1). The numbers refer to the MAPS Master database in the PC.
- If the active method (batch line 1) is loaded from the internal database in the module, the method number will be a single number between 1 and 30.

REPEAT/OMIT

With the F4 and F8 keys, you can quickly re-activate (F4=REPEAT) or cancel (F8=OMIT) one or more method steps. This function does not delete the steps in the database.

M-FUNC

See description under [Development mode](#).

START/CONTINUE

- The START  button processes the method from *the first executable step*.
- The CONTINUE  button starts ("continues") the method with *the method step already active*.

Development mode

Development mode is the way to process samples on a MAPS preparation module without using the Queue. Unlike in the Process mode, you can edit the method while the process is running. You may also repeat the method easily without having to define it again.

In Development mode all specimen holder discs will be treated by the same method until you change to another method.

Development mode is the correct mode for the special single unit without conveyors version.

Single methods

If you want to run only one single method at a time, you should use the DEVELOPMENT MODE.

- Choose a method in the Process Control screen.
- Choose holder size.
- Press \diamond to start the actual process.

Development and testing

You can change all data while using the Process Edit screen, but only the following parameters will change immediately during the process:

- Rotational speed (RPM), rotational direction, cleaning¹, dressing¹, force, lubricant and suspension dosing.

All other changes will be ignored until the process is restarted.

Development screen

Method: 1 (1,1,1)*	PROCESS CONTROL				
Description:					
Holder size: 200	Remaining time: 9:30				
STEP	POS	DISC	SUSP/#	LUB	TIME/ μ m
1	0				1:30
2	0				0:00
Type new parameter value					
F1:EDIT	F2:	F3:	F4:REPEAT		
F5:A-STOP	F6:M-FUNC	F7:	F8:OMIT		

¹ While the actual cleaning or dressing process is running changes will be ignored, while changes made before the method process reaches the cleaning/dressing step will be performed.

Keys	Explanation
ARROW UP/DOWN	Select a parameter.
ENTER	Accept a parameter value
	Starts the MAPS preparation module with the selected method. The first executable step will be the first step to start.
	Stops the MAPS preparation module.
	Continues processing of the active method. The selected step will be the first step to start.
F1 EDIT	Edits the active method while the process is running.
F4 REPEAT	Activates the previous step in the method
F5 A-STOP	Stops the preparation method process after the present method step number (Auto stop)
F6 M-FUNC	Activates the Manual functions screen
F8 OMIT	Omits the active step and moves onto the next step

- Method no. In the upper left field of the screen, you can choose between 30 methods. Type a method number between 1 and 30*) on the keypad and then press ENTER to accept the method number.
- Method number format*) The method number without parenthesis always refers to the method number of the internal database in the MAPS module. If the method has been transmitted from the MAPS Master PC program a number of the format [**X** (main group no.), **Y** (sub group no.), **Z** (method no.)] appears beside the internal method number. The PC database number will only appear if the following conditions apply:
- The method has been transmitted from the MAPS Master PC program. The data communication between the module and the PC (MAPS Master) is active.
- Specimen Holder diameter You can choose between 50 and 200 mm specimen holder discs.
- F4/F8 REPEAT/OMIT With the F4 and F8 keys, you can quickly re-activate (F4=REPEAT) or cancel (F8=OMIT) one or more method steps. This function does not delete the steps in the database.
- F5 A-STOP F5 (Auto stop) selects the actual method step as the last processed step in the method. If you need another step to be the last processed step, type in the step number and press ENTER. Pressing F5 again cancels the function.

F6 M-FUNC

Perform a number of manual functions; delivery of specimen holder, cleaning and others.

START/CONTINUE

- The START ◊ button processes the method from *the first executable step.*
- The CONTINUE ◊ button starts ("continues") the method with *the method step already active.*

F1 P-EDIT (development mode)

To change (edit) the data in a materialographic method, Struers has provided each MAPS preparation module with a method editor in the software. In the method editor, you can add or delete steps and change all method parameters. Both Development mode and the Method database screens work with the same edit function.

See also Section 10, [Method parameter values](#).

```

Method:  1      Step: 1      PROCESS EDIT
Stone   #150      Rec. wate 100
1 Stone   #150      Rec. water 100
   ①         ②         ③         ④
Time   :  --:--:--  --:--  Disc RPM :1500
S. rem:    25    75    0  Direction:→→
Force   :    100  300    0  Cleaning  :  0
Susp.   :     0    0    0  Dressing  :  40
Lub.    :    ON    ON    ON
Type new parameter value

F1:          F2:          F3:INSERT F4:
F5:          F6:          F7:REMOVE F8:
    
```

Keys	Explanation
ARROW UP/DOWN/LEFT/RIGHT	Select a parameter.
ENTER	<ul style="list-style-type: none"> Pressing ENTER on ①, ②, ③, and ④ activates an option menu. Otherwise ENTER accepts a parameter. If you know the parameter numbers by heart you may specify the number directly and press ENTER to accept without using the Option menu system.
PGUP/PGDN	Move between the steps. The list of steps is circular, meaning that step 1 follows step 10.
F3 INSERT	Insert an empty step. WARNING! This will delete step number 10.
F7 REMOVE	Remove a step. This function moves <i>all the following steps 1 step back</i> and creates an empty step at the end of the list (step 10). WARNING! Lost (removed or deleted) data cannot be restored!

WARNING!
You cannot insert/remove a step while a process is running!

Option menus

The Option menus for parameters marked with ①, ②, ③ and ④ on the Development mode screen (Process edit) are selected with the ENTER key. Use the ARROWS to select the desired consumable/setting and press ENTER to accept the choice. (See also Section 10, Method parameter values.)

Disc/cloth parameters

OPTION: DISC/CLOTH			
0:	1: Stone	2: Dia.PAD	
3:	4: `DiaPAD	5: Dia.PAD	
6: `Dia.PAD	7:	8: PG-Paper	
9:	10: SiC-Paper	11: `SiC-Paper	
12:	13: P-M	14: `P-M	
15:	16: DP-Net	17: DP-Plan	
↓	Press ENTER key to accept choice		
F1:	F2:	F3:	F4:
F5:	F6:	F7:	F8:

Grit/grain size parameters

OPTION: GRIT/GRAIN SIZE			
0:	1: #40	2: #50	
3: #60	4: #80	5: #100	
6: #120	7: #150	8: #180	
9: #220	10: #240	11: #280	
12: #320	13: #360	14: #400	
15: #500	16: #600	17: #800	
↓	Press ENTER key to accept choice		
F1:	F2:	F3:	F4:
F5:	F6:	F7:	F8:

Suspension type parameters

OPTION: SUSPENSION			
0: <input type="checkbox"/>	1: 15µm	2: 9µm	
3: 6µm	4: 3µm	5: 1µm	
6: .25µm	7:OP-U	8:OP-S	
9:OP-S	10:	11:	
12:	13:	14:	
15:			
Press ENTER key to accept choice			
F1:	F2:	F3:	F4:
F5:	F6:	F7:	F8:

Lubricant type parameters

OPTION: LUBRICANT			
0:	1:Water	2: Rec. water	
3:Blue	4:Red	5:	
6:	7:	8:	
9:	10:	11:	
12:	13:	14:	
15:			
Press ENTER key to accept choice			
F1:	F2:	F3:	F4:
F5:	F6:	F7:	F8:

F6 M-FUNC, Manual functions

F6=M-FUNC, Manual functions contains 3 control screens for timing, conveyors and cleaning of the dosing system.

Screen 1 of Manual functions

MANUAL FUNCTIONS			
CHOOSE ACTION:	ENTER PARAM.:		
Perform cleaning	Clean prg. : 2		
Perform dressing	Dressing [μ M]: 50		
Deliver specimen holder	Time : 0		
↓ <input type="text"/>			
Press ENTER key to accept choice			
F1:	F2:	F3:	F4:
F5:	F6:	F7:	F8:

Perform cleaning: activate the cleaning station. You may also change the cleaning program number.

Perform dressing: (only active on the grinding/polishing type of module). You may also change the dressing depth.

Deliver specimen holder: deliver the specimen holder on the conveyor belt and start the belt. You may change the running time of the conveyor belt (0-20 sec).

Screen 2 of Manual functions

MANUAL FUNCTIONS			
CHOOSE ACTION:	ENTER PARAM.:		
↑ <input type="text"/>			
Start delivery conveyor	Time: 0		
Start feed conveyor	Time: 0		
↓ <input type="text"/>			
Press ENTER key to accept choice			
F1:	F2:	F3:	F4:
F5:	F6:	F7:	F8:

Start delivery/feed conveyor: start feed or delivery conveyor. You may also change the running time of the conveyor belt (0-20 sec). The screen depends on the conveyor configuration of the preparation module.

Screen 3 of Manual functions

MANUAL FUNCTIONS			
CLEANING OF DOSING SYSTEM	P2	P1	
↑ <input style="width: 100%; height: 1.2em;" type="text"/>			
Lub. nozzles (0=no,1=yes)	0	0	
Susp. nozzles (0=no,1=yes):	0	0	
OP-S nozzles (0=no,1=yes):	0	0	
ACTION:		Start Start	
Type new parameter value			
F1:	F2:	F3:	F4:
F5:	F6:	F7:	F8:

Cleaning of dosing system: See the instructions on OPS/OPU suspension in Section 4, Maintenance. Cleaning of lubricants, suspension and OPS/OPU suspension should be done regularly. You only have to write the cleaning parameters once, as the module software saves the latest instructions. See details under daily and weekly cleaning.

Keys	Explanation
ARROW UP/DOWN	Select a function or a parameter.
PGDN/PGUP	Toggle between the 3 Manual functions screens.
ENTER	Activate a function or accept a parameter.

Method database

Method database screen

Through the Method database screen you can access the database of the MAPS preparation module directly.

Method: 20* (1, 1,20) METHOD DATABASE						
Description: Rails						
STEP	POS	DISC	#/μm	SUSP	LUB	T/μm
1	1	Stone	#60		Rec.	125
2	1	Stone	#60		Rec.	0:05
3	2	DP-Plan		15μm	Blue	5:05
4	3	DP-Plus		3μm	Blue	5:00
↓ 5	4	DP-Nap				0:35
Type new parameter value						
F1:EDIT		F2:		F3:COPY		F4:
F5:		F6:		F7:DELETE		F8:

Keys	Explanation
PGUP/PGDN	Scroll through the 30 methods in the MAPS preparation module database.
ARROW UP/DOWN	Select a method number or a method step.
ENTER	<i>When selecting a method number:</i> type a method number and press ENTER <i>When selecting a step number:</i> press ENTER and edit the desired step in the active method no.
F1 EDIT	Method Edit enables you to view and edit all the parameters of each of the 10 steps in the selected method.
F3 COPY	Copies the contents of the method to another method in the MAPS preparation module database.
F7 DELETE	Deletes a method in the MAPS preparation module database.

Method number format*)

The method number without paranthesis always refers to the method number of the internal database in the MAPS module. If the method has been transmitted from the MAPS Master PC program a number of the format [**X** (main group no.), **Y** (sub group no.), **Z** (method no.)] appears beside the internal method number. The PC database number will only appear if the following conditions apply:

- The method has been transmitted from the MAPS Master PC program.
- The data communication between the module and the PC (MAPS Master) is active.

F1 EDIT (Method database)

Editing of step parameters similar to Development mode (see also Section 10, Method parameter values).

Method: 20	Step: 1	METHOD EDIT	
POS DISC	#/µm	SUSP	LUB T/µm
1 Stone	#60		Rec. water 125
①	②	③	④
Time : --:--	--:--	--:--	Disc RPM :1500
S. rem: 25	100	0	Direction:→ →
Force : 200	500	0	Cleaning : 0
Susp. : 0	0	0	Dressing : 20
Lub. : ON	ON	ON	
Type new parameter value			
F1:	F2:	F3:INSERT	F4:
F5:	F6:	F7:REMOVE	F8:

Keys	Explanation
ARROW UP/DOWN/LEFT/RIGHT	Select a parameter.
ENTER	<ul style="list-style-type: none"> Pressing ENTER on ①, ②, ③, and ④ activates an option menu (see under Development mode). Otherwise ENTER accepts a parameter. If you know the parameter numbers by heart you may specify the number directly and press ENTER to accept without using the Option menu system.
PGUP/PGDN	Move between the method steps. The list of steps is circular, meaning that step 1 follows step 10.
F3 INSERT	Insert an empty step. WARNING! This will delete step number 10.
F7 REMOVE	Remove a step. This function moves <i>all the following steps 1 step back</i> and creates an empty step at the end of the list (step 10). WARNING! Lost (removed or deleted) data cannot be restored!

INSERT/REMOVE

F3=INSERT and F7=REMOVE change the order of the 10 steps in the method. When inserting a new step, you delete the last step (number 10). When removing a step, you move all steps, below the removed step, 1 number back.

F3 COPY

Copy the contents of the actual method to another method in the database of the MAPS preparation module.

METHOD COPY			
COPY FROM Method No.: 20			
COPY TO Method No.: 20			
Type new parameter value			
F1:	F2:	F3: COPY	F4:
F5:	F6:	F7:	F8:

Keys	Explanation
ENTER	Accept the method number.
F3 COPY	Copy the method to the destination

Please note that the additional method number in MAPS Master format [X,Y,Z] will not be copied to the new method in the module database (see also Section 2, Method number format).

F7 DELETE

Delete a method in the MAPS preparation module database.

METHOD DELETE
DELETE Method No.: 20
Use cursor keys to select parameter
F1: F2: F3: F4:
F5: F6: F7:DELETE F8:

Keys	Explanation
F7 DELETE	Delete the method. WARNING Deleted data cannot be restored!

Consumables

Refill/dispose liquids

Cleaning liquids

Jars for the soap solution and alcohol are placed on a trolley behind the recirculation tank.

It is important to connect the 2 different cleaning liquid supplies to the correct hoses:

Correct cleaning hose connections

- BLACK HOSE = Soap
- BLUE HOSE = Alcohol

A level indicator warns when the liquid level in the soap and alcohol jars subsides to 15 mm.

- Stop MAPS. The module will automatically stop if the trolley is removed while MAPS is running.
- Pull the trolley free of the unit.
- Unscrew the lid of the jar and pull the hose and sensor out of the jar.
- Fill the jar and remount the correct lid/hose on the jar (*Black hose for soap and blue hose for alcohol*).
- Push the trolley back into its place. You can now start MAPS again.

Alcohol separation unit

The jar of the MAPAL alcohol separation unit is placed on the trolley beside the soap and alcohol jars. A transparent hose connects the MAPAL jar to MAPS.

A level indicator warns when the jar is full.

Follow the above procedures for removing the trolley. Remember to put the right lid/hose on the jar after emptying the contents.

Recirculation cooling unit

The Recirculation cooling unit is placed on a trolley. A level indicator warns when the water is low.

- Stop MAPS. The module will automatically stop if the trolley is removed while MAPS is running.
- Pull the trolley free of the unit. Pull out the spiral cable from behind the electric box (see photo to the left).
- Click the red snap-lock on the lids of the soap and alcohol jars in the back of the trolley and disconnect the hoses (on the MAPAL waste alcohol jar, the hose is pulled out of the lid).

- Lift the motor with pump out of the cover and put it aside, carefully.

- Drive the trolley to the refill premises.
- Remove the tank cover and place a plastic insert in the tank (see Accessories, section 5). Always place the insert with the hole in front of the outlet hole in the tank.
- Fill the tank with 63 l of water and 2 l Struers Additive for cooling fluid. The water level should be 8-10 cm below the upper edge of the tank.

- Remount the tank cover.

Remember...

Check the correct direction of the cover

- Return the trolley to the MAPS cabinet and mount the motor and pump in the tank cover.
- Mount the hoses for soap and alcohol by inserting the hoses through the respective lids (*Black hose* for soap and *blue hose* for alcohol). Pull up the snap-lock. If you have a MAPAL waste alcohol jar, the hose (*transparent*) is screwed into the lid.
- Connect the spiral cable plug to the electric box on the trolley.
- Push the trolley back into its place. You can now start MAPS again.

Suspension and lubricant

The diamond suspension is contained in the small bottles and the lubricant is contained in the large bottles. The OPS/OPU suspension is contained in the large bottles.

- Manual check: the level is low if it sinks below the plate where the bottle is fixed.
- Automatic check: if you have the MAPDU level indicator, the display warns you if the level in the bottle is low.
- Stop MAPS.
- Suspension: unscrew the lid of the bottle and pull the lid with the hose out of the bottle.
- Lubricant: pull the lid of the bottle.
- Refill the bottle.

IMPORTANT

- Never fill the suspension bottle more than 2/3 full, as the suspension will be bubbled through with compressed air.
- Always clean the OPS/OPU-suspension bottle before refill in order to avoid clogging the system with sediments.

- Remount the bottle in the plate and put the lid on again. Make sure that it is the correct lid!
- Activate the predose program after every refill.

Changing preparation discs

- Stop the MAPS preparation module.
- Open the safety guards. Place the disc on the turntable. The 3 pins fit the corresponding holes in the turntable.

The disc must lie absolutely flat against the turntable to avoid deformations during method operation. Never use adhesive material or glue on the turntable. The contact surface between turntable and disc must be clean and free of loose particles.

We recommend manual predosing if the disc is new.

WARNING

The safety doors are closed while methods are being processed.

If you need access to MAPS, you must:

- Stop the method processing
or
- Wait until the specimen holder disc is in the cleaning station

Single unit without conveyors

- Once the process has been started by pressing **START** the single unit without conveyors runs until the disc has to be replaced by another type of disc. This is indicated by the message: **New consumables, press STOP or CONTINUE** which shows that the disc on the second work station needs to be replaced (see the instruction above).
- Open the safety guard and replace the disc with the correct type of disc. Close the safety guard and press **CONTINUE**.

An alternative way to change the discs, avoiding the "New consumables, press STOP or CONTINUE" pause:

- Wait until the specimens are in the cleaning station. Open the safety guard. The machine does not stop. Change the disc before the cleaning time expires (if the cleaning time should expire anyway the machine automatically stops).
- Close the safety guard again and press **CONTINUE**. The machine has now registered that the disc has been changed and will proceed without stop.

Changing grinding stone

When the grinding stone is worn down to a thickness of about 25 mm, MAPS issues a warning: "STONE AT MINIMUM" and prevents further grinding. To continue grinding, you must replace the grinding stone. This is done as follows:

- Stop MAPS momentarily by pressing STOP .

WARNING

Never change grinding disc/diamond disc when MAPS is turned off with the main switch as you risk damaging the disc and the module.

- Remove the cover over the grinding stone by turning the 4 finger screws and lifting the cooling pipe. Loosen the center screw in the flange on the grinding stone with an 8 mm hexagon spanner. You can now remove the stone.
- Immediately after having removed the cover, a new screen appears:

DISC CHANGED ?

Stone cover has been opened, change the module type according to the actual grinding disc/stone.
(Stone=0, Diamond disc=2)

Module type: **Grind/Polish**

Use numeric keys to select option

F1:	F2:	F3:	F4:
F5:	F6:	F7:	F8:

- Depending on the desired type of disc, you press 0=Stone or 2=Diamond.

WARNING

This procedure changes the module type!!!
(1=Polish is not available from this screen).

- Clean the rubber washer protecting the stone. Mount the grinding stone with the recessed side upwards. The center stud of the module fits into the grinding stone center hole. Place the upper flange with the cardboard washer in the recess and tighten the center screw firmly.
- Remount the cover. Remember to lock the safety switch key into its hole before remounting the 4 finger screws.
- Tilt the cooling pipe into position in the hole again.
- When the preparation module is started again the dresser automatically adjusts to the new stone.

3. Preparation of specimens

Preparation methods

PG Plane Grinding

Plane Grinding (PG) may be carried out on a grinding stone, PG-paper, SiC paper or, for hard materials, on Diamond Pad.

SiC paper is used on a wet grinding disc, according to the Knuth-Rotor principle. PG-paper or Diamond Pad is glued to a grinding/polishing disc.

FG Fine Grinding

Fine Grinding (FG) is carried out on a cloth with low resilience such as DP-Plan, DP-Pan or DP-Dur, or on Petrodisc-M. SiC-Paper may also be used.

DP Diamond Polishing

Diamond Polishing (DP) is carried out on a polishing cloth using diamond abrasive and a lubricant. The correct choice of cloth depends on the material.

OP Oxide Polishing

Oxide Polishing (OP) is especially suited for soft and ductile material as a final polishing step.

Struers Metalog Guide™

Struers Metalog Guide™ offers preparation methods for the most common materials, based on a simple analysis of two key properties: hardness and ductility. Finding the right method is easy, including choice of consumables. Always consult Struers Metalog Guide™ for the correct preparation method for the actual specimens.

Metalogram

Metalog Methods

Struers Metalog Guide™ contains 6 useful chapters:

A quick and safe guide to the right preparation method.

A complete catalogue of preparation methods, based on Struers' vast store of materialographic experience, and employing Struers' range of consumables.

Preparation Philosophy

The basics of modern specimen preparation, seen from a professional point of view.

Metalog Process

The materialographic preparation process from start to finish, logically explained.

Metalog Master

A combined trouble-shooting guide and supply of in-depth information about the processes of mechanical preparation, including an expert system for the solving of preparation problems.

Metalog Code

Quick access to the relevant consumables for the chosen preparation methods.

Struers Metalog Guide™

A complete guide to materialographic specimen preparation. Contact your local dealer for a free copy of Metalog Guide™.

4. Maintenance

Mechanical maintenance

Daily

If any of the polishing cloths become dirty then:

- Carefully brush with a clean, soft nail brush under lukewarm running water.
- Rinse with distilled water

When starting again do not forget to pre-dose with diamond suspension. See also Section 2, Operation. Worn-out cloths should be replaced.

Suspensions/lubricant

Refill suspension and lubricant bottles if the level sinks below the plate in which the bottle is placed. Do not fill the suspension bottle more than 2/3 full as they are bubbled through with compressed air between operations.

OPS/OPU-suspension

Clean the OPS/OPU-suspension bottle before refilling in order to avoid clogging of the system with sediments.

Grinding stone

Clean the drain every day.

*Every evening or before breaks
of more than 8 hours*

IMPORTANT

Grinding stone station

- Take off the protective cover above the grinding stone.
- Remove accumulated grinding grains and grinding dust from the stone station.
- Pull out the recirculation tank and place a bucket underneath the drain of the grinding station.
Clean tray and the entire drain system with a brush (using the "chimney sweeper method"). Flush with water.
- Remove the bucket and reposition the recirculation tank.
- Put the cover back. At "Start" MAPS automatically adjusts the diamond dresser. This is for safety reasons as it "thinks" that the stone has been exchanged.

OPS/OPU-suspension

Clean and empty the oxide suspension hose, pump and valves as follows:

- Choose the PROCESS MODE menu and press F2, Process control.
- Press F6, Manual functions and press PGDN twice. The Cleaning of dosing system screen appears.
- Replace the OPS/OPU suspension bottle with a bottle of water which is at room temperature.
- Activate the OP-S cleaning nozzles by typing "1" for one or both of the 2 module work stations. Move to the "Start" row with the arrows and press ENTER for the active work stations. The cleaning takes 2 min. You may stop the cleaning by pressing ⏏.

Every second day
The trolley jars

- Refill the alcohol jar. Use ethanol or possibly propanol.
- Refill the soap jar. Blend a soap solution according to the instruction on the label.
- Refill the recirculation tank until the level is 8 cm from the upper rim. Mix with Additive according to the instruction.

Every week or before breaks of more than 24 hours
Lubricant/Suspension

- Remove all discs.
- Clean and empty the lubricant/suspension tubes and nozzles as described for OPS solution under evening cleaning.
- Take out the water bottle and press ABRASIVE  until tube and nozzle are empty of water.
- Repeat this procedure for each diamond suspension bottle, starting with the smallest and ending with the largest grain size, so that the same bottle of water can be used and at the same time prevent the larger diamonds from mixing with the smaller ones.

Spill pans

- Take off polishing/grinding discs and splash rings.
- Remove accumulated grinding/polishing waste with a rounded scraper.
- Clean the spill pan with a small brush and a little water.
- Wash the splash rings.

Cleaning of cleaning station

- Clean the nozzles. If the nozzle holes are clogged, you can remove the waste by carefully using a thin needle.

Every month
Check and empty air filter

Behind MAPS are 2 air filters (coarse and fine). Through the holes in the screen you can view the condition of the filters.

- Check the condition of the filters.
- Activate the valve at the bottom of the filter to remove water and oil.

Recirculation cooling unit

Water and additive must be changed to avoid bacterial contamination which can emit an unpleasant smell.

Every time water is changed, the tank should be cleaned thoroughly to remove all grinding waste. By using Struers plastic inserts (see Accessories, section 5) the cleaning is kept to a minimum.

<p style="text-align: center;">Remember... to blend the cooling water with Additive for Cooling Fluid</p>
--

Cleaning station

- Dismount the splash ring and remove it from the cleaning station.
- Dismount the grating and clean it with a brush and a mild soapy solution.
- Clean the inner wall of the cleaning station with brush and soap. Rinse with a little water.
- Mount the grating and the splash ring again.

Plastic guards

Use a soft damp cloth and ordinary anti-static window cleaning agent.

Painted surfaces

Use a damp cloth and water with a dash of ordinary washing-up liquid.

WARNING!

Do not use alcohol, acetone or similar solvents.

Suspension nozzle

Dismantling and cleaning of suspension nozzle:

- Remove (and keep!) the 4 screws holding the nozzle plate.
- Dismount the suspension hose and remove the nozzle plate.
- Screw a M3 screw behind each nozzle point and pull the point out.
- Clean the nozzle plate and the points in warm (soapy) water.
- Clean the nozzle plate and the points with compressed air.
- Mount the nozzle points in the nozzle plate.

Remember...

to control the O-rings in the nozzle points
before they are pushed in place in the nozzle plate.


- Mount the hoses.
- Tighten the nozzle plate again with the 4 screws.

Software maintenance

Configuration menu

Change consumables, cleaning program or module ID numbers etc.

Keys	Explanation
ARROW UP/DOWN	Select a menu.
ENTER	Activate a screen/menu.

Discs/Cloths, configuration

Change disc/cloth parameters.

CONFIGURATION: DISC SIZES				
NO.	DESCRIPTION	TYPE	DISC	DIAMETER
0:			300	[mm]
1:	Stone	Stone		[mm]
2:	D. Disc	Diam.	356	[mm]
3:			300	[mm]
↓ 4:	'D.Disc	Diam.	230	[mm]
Type new parameter value				
F1:	F2:	F3:	F4:	
F5:	F6:	F7:	F8:	

Keys	Explanation
ARROW UP/DOWN/ LEFT/RIGHT	Select a parameter.
PGDN/PGUP	Scroll between the screens.
ENTER	Accept a new parameter value.

No. The "No." column shows the reference number assigned to a certain disc or cloth. This number is used under the "Edit" function to identify the type to be used.

Description/type The "Description" column shows all of the disc/cloth descriptions of the factory setting. You can change the descriptions using the optional MAPS Master PC program. The "Type" column shows the preset disc/cloth type.

Disc diameter Type the \emptyset mm of the disc/cloth and accept the choice by pressing ENTER. The position of the transportation arm depends on the diameter of the rotating disc (disc diameter).

WARNING

An incorrect diameter of the rotating disc (disc diameter) may result in damage to the specimen and/or the rotating disc. Be careful to check the Parameter value table in section 10, if in doubt.

Lubricants, configuration

Assign a lubricant specification (description, type and position) to a physical lubricant bottle.

No.	Descr.	Type	P3	P4	CONFIG. : LUBRICANTS
0			0	0	
1	Water	Water			
2	Rec. water	Rec.wa.	0	1	
3	Blue	Lubr.	11	11	
4	Red	Lubr.	0	21	
5			0	0	
6			0	0	
↓7			0	0	
Type new parameter value					
F1:	F2:	F3:	F4:		
F5:	F6:	F7:	F8:		

Keys	Explanation
ARROW UP/DOWN/ LEFT/RIGHT	Select a parameter.
PGDN/PGUP	Scroll between the screens.
ENTER	Type a new bottle number and press ENTER.

No. The "No." column shows the reference number assigned to a certain lubricant bottle. This number is used under the "Edit" function to identify the type to be used.

Description/type The "Description" column shows the lubricants specification of the factory setting. You can change the descriptions from the optional MAPS Master PC program. The "Type" column shows the preset lubricant type.

Dosing units (P1-P5) The Px numbers (P1 to P5) refer to the work station position according to the total number of modules (i.e. P3 = work station 1 on module no. 2). Each lubricant specification is linked to one or more bottle connections. The bottles and the work stations are coded according to Section 9, Technical data.

IMPORTANT
The bottles must be physically present by the start of the method.

Suspensions, configuration

Assign a suspension specification (description, type and position) to a physical suspension bottle.

No.	Descr .	Type	P3	P4	CONFIG. : SUSPENSIONS
0			0	0	
1	15µm	DP	0	0	
2	9µm	DP	0	0	
3	6µm	DP	12	12	
4	3µm	DP	13	13	
5	1µm	DP	0	22	
6	.25µm	DP	0	23	
↓7	OP-U	OP	0	31	
Type new parameter value					
F1:	F2:	F3:	F4:		
F5:	F6:	F7:	F8:		

Keys	Explanation
ARROW UP/DOWN/ LEFT/RIGHT	Select a parameter.
PGDN/PGUP	Scroll between the screens.
ENTER	Type a new bottle number and press ENTER.

No. The "No." column shows the reference number assigned to a certain suspension bottle. This number is used under the "Edit" function to identify the type to be used.

Description/type The "Description" column shows the suspension specification of the factory setting. You can change the descriptions from the optional MAPS Master PC program. The "Type" column shows the preset suspension type.

Dosing units (P1-P5) The Px numbers (P1 to P5) refer to the work station position according to the total number of modules (i.e. P3 = work station 1 on module no. 2). Each suspension specification is linked to one or more bottle connections. The bottles and the work stations are coded according to Section 9, Technical data.

IMPORTANT
The bottles must be physically present by the start of the method.

Cleaning programs, configuration In this screen you can change the cleaning data of the system.

```

Program no.: 4 CONFIG.: CLEANING PRG.
Total cleaning time [min:sec]: 1:30

STEP MEDIA      TIME  STEP MEDIA      TIME
1 Alcohol      10    6 Upper air      10
2 Lower air     10    7 Lower air      35
3 Alcohol       10    8 None           0
4 None          10    9 None           0
5 Alcohol       5     10 None           0
Press ENTER key to select option menu

F1:           F2:           F3:           F4:
F5:           F6:           F7:           F8:
    
```

Option menu of Cleaning programs

```

OPTION: CLEANING MEDIA

0:None
1:High pressure water
2:Low pressure water
3:Soap
4:Alcohol
5:Lower air
6:Upper air

Press ENTER key to accept choice

F1:           F2:           F3:           F4:
F5:           F6:           F7:           F8:
    
```

Keys	Explanation
ARROW UP/DOWN/LEFT/RIGHT	Select a parameter.
PGDN/PGUP	Scroll between 10 cleaning programs.
ENTER	Accept parameter or invoke the cleaning media Option menu.

See also [Cleaning station](#) in Section 9, Technical data.

Module configuration

Change the configuration of certain basic module functions. Should only be performed by service technicians approved by Struers.

Screen 1 of Module configuration

```

CONFIGURATION: MODULE

Module type: Polish/Polish
Module No.           (1 - 5):  2
No. of modules      (1 - 5):  2
No. of input belts  (0 - 3):  1
No. of output belts (0 - 3):  1
Language            (0=English,1=German):  0
Drying time         (0-255 sec): 60
Use numeric keys to select option

F1:      F2:      F3:      F4:
F5:      F6:      F7:      F8:
```

Screen 2 of Module configuration

```

CONFIGURATION: MODULE

LVDT Installed      (0=no, 1=yes):  1
Cloth detection sensitivity :4000
Stone cover instal. (0=no, 1=yes):  1
Cooling level on P1   (0 - 10):  0
Cooling level on P2   (0 - 10):  10
Susp./lub sensors    (0=no, 1=yes):  1
Comp. after dresser adj. [µm]: 250
Type new parameter value

F1:      F2:      F3:      F4:
F5:      F6:      F7:      F8:
```

Keys	Explanation
ARROW UP/DOWN/ LEFT/RIGHT	Select a parameter.
PGDN/PGUP	Scroll between the screens.
ENTER	Accept a parameter.

WARNING!

The Module configuration screen should generally **only** be operated by a technician authorized or trained by Struers

Screen 1 of Module configuration Module type etc.: default information on the physical MAPS preparation module set up. This information should not be changed by the user.

Drying time: if your MAPS preparation module is equipped with the optional *drying fan* unit (MAPFA) under the conveyor belt, this parameter will control the drying time for each specimen holder disc after finishing the process in the module.

Screen 2 of Module configuration The information in this screen should generally not be changed by the user, except for:
Cooling level: if your MAPS preparation module is equipped with a cooling arrangement (MAPOL), this parameter will control the water supply. See also Section 10, Method parameter values.

Comp. after dresser adj.: in the module configuration menu it is possible to state the number of μm to move the dresser stop as a compensation for the drilling movement made by the dresser needle during of dressing adjustment. Default value is 250 μm . The value should be increased if a soft stone is used and vice versa.

5. Accessories

Specification	Code word
Preparation modules	
<i>MAPS Preparation Module</i> with 1 plane grinding station and 1 fine grinding/polishing station (150/300 rpm), with Recirculation cooling unit (MAPRE)	MAPSE
<i>MAPS Preparation Module</i> with 1 plane grinding station and 1 fine grinding/polishing station (variable speed), with Recirculation cooling unit (MAPRE)	MAPON
<i>MAPS Preparation Module</i> with 2 fine grinding/polishing stations (150/300 rpm)	MAPFI
<i>MAPS preparation module</i> with 2 fine grinding/polishing stations (variable speed)	MAPVA
Conveyor belt	
<i>Conveyor</i> Feed module Connection module Delivery module	MAPEF MAPCO MAPLI
<i>Manual Feeding Station</i> for operating MAPS without conveyors	MAPEX (MAPST)
Water recirculation	
<i>Recirculation cooling unit</i> contents 65 l, with a rolling pallet	MAPRE
<i>Recirculating Connection Kit</i> for connection of 2 work stations to one Recirculation cooling unit	MAPIT
<i>Cooling Arrangement</i> for preparation disc	MAPOL
<i>Disposable Inserts</i> for MAPS and Exotom Recirculating Unit for collecting waste material. Package with 10 pcs.	EXOSP


Specification	Code word
Liquids supply	
<i>Dosing system</i> box with room for 3 bottle units. Bottle units (MAPNU or MAPOU) required	MAPTU (MAPTH)
<i>Bottle unit</i> for dosing system. With room for 2 DP-Suspensions and 1 Lubricant	MAPNU (MAPSU)
<i>Bottle unit</i> for dosing system. With room for 2 OP-Suspensions	MAPOU (MAPOP)
<i>Level indication unit</i> for monitoring the level of suspensions and lubricants in the bottles. Level indicators (MAPEV or MAPVE) required	MAPIN
<i>Level indicator</i> for bottle unit with DP-Suspension (MAPNU)	MAPEV (MAPDU)
<i>Level indicator</i> for bottle unit with OP-Suspension (MAPOU)	MAPVE (MAPPO)
<i>Stock removal unit</i> for measuring of the material removal during the plane grinding process and monitoring of abnormal vibrations	MAPOM
<i>Separation unit</i> for collection of waste alcohol from the cleaning station	MAPAL
Specimen handling	
<i>Drying fan</i> for drying of specimen holders	MAPFA
<i>Conveyor ring</i> for mounting on existing specimen holder: Ø160 mm Ø200 mm	MAPEL MAPIR
<i>Uniforce</i> Levelling device for levelling of specimens in specimen holders Ø 140, 160 and 200 mm. With pressure foot for positioning of specimens (MAXMO)	MAXFU

Specification	Code word
Specimen holder discs	
<p><i>Specimen holders for MAPS</i> Ø 160 mm, stainless, mounted with conveyor ring (MAPEL):</p> <p>6 specimens, Ø 12-40 mm (MAXCY) 12 specimens Ø 10-25.5 mm (MAXUD) 6 flat specimens max. 34x60 mm (MAXYK) 3 flat specimens max. 41x72 mm (MAXDI)</p>	<p>MAPCY MAPUD MAPYK MAPDI</p>
<p><i>Specimen holders for MAPS</i> Ø 200 mm, stainless, mounted with conveyor ring (MAPEL):</p> <p>3 flat specimens max. 54x94 mm (MAXOT) 5 specimens Ø 27-55 mm (MAXUK)</p>	<p>MAPOT MAPUK</p>
<p>Specimen holders for MAPS plane, solid specimen holder for adhesion of specimens or for custom-made shapes, stainless, mounted with conveyor ring (MAPEL)</p> <p>160 mm (MAXSO) 200 mm (MAXON)</p>	<p>MAPOS MAPNO</p>
Preparation stones and discs	
<p><i>Grinding stone</i> Ø 356 mm, grit 60, mainly for soft and ductile metals</p>	<p>ABGOF</p>
<p><i>Grinding Stone</i> Ø 356 mm, grit 150, suited for most kinds of metals</p>	<p>ABGAL</p>
<p><i>Petrodisc-M</i> grinding disc for plane one-step fine grinding, Ø 290 mm</p>	<p>PETIL</p>
<p><i>Polishing disc</i> Ø 300 mm, of aluminium. Extra plane, especially for MAPS</p>	<p>DEDMA</p>

Specification	Code word
Software	
<p><i>MAPS Master</i> PC software with method database/editor. Operates and controls the MAPS preparation modules. Converter for communication included.</p>	MAPMA
Tools	
<p><i>Compressor</i> Jun-air compressor for air supply to MAPS preparation modules. Type 24-50 with displacement 200 l/min. 1 compressor is recommended for each MAPS preparation module. 1 compressor can supply 2 MAPS preparation modules with intermittent use</p>	MAPJU
<p><i>Diamond Point for Dressing Tool</i></p>	PLAMT

Code words in parenthesis are valid only for serial numbers below 4590017.

Struers
consumables


Struers offers a wide range
of consumables for grinding and polishing.
Ask for separate leaflets from your Struers dealer.

6. Installation


MAPS must be installed by Struers or by an authorized service technician, trained by Struers for this specific task.

External hose connections

	Tube connections	Hoses delivered by Struers	Connections to supply source or outlet
Compressed air	Pipe branch Ø 8 mm	Ø 6.5 mm x 3 m	To be fitted with local standard connections
Water inlet	¾" British Standard pipe thread	Ø 10 mm x 2.5 m with standard connector	For ½" or ¾" British Standard pipe thread
Water outlet	Pipe branch Ø 32 mm	Ø 32 mm x 3 m	To be fitted with local standard connectors

Connection tube collar for optional exhaustion

1459_035.CDR

MAPS data communication and control

With a PC

With the *MAPS Master PC program* you can control a whole production line consisting of up to 5 MAPS preparation modules. The PC program can download methods to any of the modules (or to all the modules at once), start, stop or continue the process and edit a very large number of methods.

With the MAPS Master PC program you get the optimal production of preparations run with a maximum number of 30 batch lines each comprising up to 99 specimen holders.

Without a PC

Each MAPS preparation module can be operated direct from the front panel of each module. In the case of 2 or more modules, the modules do not communicate through cables but are activated by sensors on the conveyor belts.

In some cases this can be an advantage but more often the MAPS Master PC program will be a much easier way to operate the MAPS preparation modules as one automated production line.

7. Trouble-shooting

Error code listing

Error code listing for MAPS ver. 2.10.00.

Function Errors (Fatal errors)

- 8** " E8: Unable to mount specimen holder "
The specimen holder coupling has not come down completely while attaching to the specimen holder.
The "3/4 down" sensor has not been activated within the fixed period of time.
- 9** " E9: Holder not down, clean or dismount "
The specimen holder has not come down completely in connection with cleaning or decoupling.
The safety switch has not been activated within the fixed period of time.
- 10** " E10: Force on holder not achieved "
The specified force on the specimen holder has not been achieved.
The pressure-OK flag has not been raised within the fixed period of time.
- 11** " E11: Unable to move arm, holder not up "
The specimen holder head is not up, as initial position for turning the transportation arm.
The "holder up" sensor has not been activated.
- 12** " E12: Holder not up after UP-command "
The specimen holder has not come up after a up-command.
The "holder up" sensor has not been activated within the fixed period of time.
- 13** " E13: Holder not down after DOWN-command "
The specimen holder has not come down after a down-command.
The "1/2 down" sensor has not been activated within the fixed period of time.
- 14** " E14: Unable to move transportation arm "
The transportation arm has not moved.
The motor is defective, the power to the motor has been cut off or the motor card is defective.

- 17** " E17: Dresser stop not up before adj. "
The dresser stop (step motor + nut) has not been raised completely to upper position in connection with a dresser adjustment.
The dresser depth sensor "Stone worn down" has not been activated within the fixed period of time.
- 18** " E18: Dresser not up after UP-command "
The dresser has not come up in connection with an up-command or by the end of an adjustment procedure.
The dresser depth sensor "dresser up" has not been activated within the fixed period of time.
- 19** " E19: Dresser stop not reached "
The dresser has not come down after a down-command.
The dresser depth sensor "dresser down on stop switch" has not been activated within the fixed period of time.
- 20** " E20: Dresser not down during adjustment "
The dresser has not come down in connection with an adjustment procedure.
The dresser depth sensor "dresser up" has not been deactivated within the fixed period of time.
- 21** " E21: Dresser stop not reached during adj. "
The dresser stop has not been lowered completely to the switch elements during an adjustment procedure.
The dresser depth sensor "dresser down on stop switch" has not been activated within the fixed period of time.
- 22** " E22: Cleaning water pressure too low "
The water pressure is too low during cleaning.
The pressure switch for minimum pressure has not been activated within the fixed period of time.
- 23** " E23: Cleaning water pressure too high "
The water pressure is too high during cleaning.
The pressure switch for maximum pressure has been activated.
- 24** " E24: Abnormal vibrations during process "
The process has stopped due to a torn cloth.
The measured amplitude on the cloth is higher than the limiting value in the Module CONFIGURATION.

- 26** " E26: Grinding shield not over cleaning "
The movable lid over the stone not in position over the cleaning station.
The sensor "lid in cleaning position" has not been activated within the fixed period of time.
- 27** " E27: LVDT unable to detect holder down "
The removal sensor has not detected that the specimen holder has come down.
The signal from the removal sensor is interpreted as if the specimen holder were constantly moving. **Note** This error code has been omitted from version 1.06.00.
- 28** " E28: Grinding shield not over stone "
The movable lid over the stone has not been positioned perfectly over the stone.
The sensor "lid in stone position" has not been activated within the fixed period of time.
- 29** " E29: Detaching boom failure "
The boom for decoupling of specimen holders has not been raised (in decoupling position).
The sensor "boom up" has not been activated within the fixed period of time.
- 31** " E31: Unable to move arm, dresser not up "
The transportation arm has attempted a rotating movement with the dresser low at a time when the dresser was not supposed to be lowered.
- 32** " E32: Unable to find reference position "
The transportation arm cannot find its reference position.
Probably an error at the end-stop sensors, the reference sensor or the tachometer.
- 41** " E41: Reverse mounting of holder disc "
Reversed coupling of the specimen holder, the coupling has come all the way down before the cup has come up.
- 42** " E42: Disc type does not fit module type "
The preparation module type and the disc/cloth of the active method step do not match.
An ERROR 42 appears in the following cases:
- Preparation module type = pol/pol and support = "Stone"
 - Preparation Module type = dia/pol and support = "Stone"
 - Preparation Module type = plan/pol and support = "D. disc" or "Diamond Pad" and working position = 1.

WARNINGS (non-fatal errors)

- 2** " W2: One or both upper guards open "
One or both upper guards are open. Close the guards before restarting the preparation module.
This warning code has been omitted from ver. v.1.06.00. Instead a message without any reference number is displayed.
- 5** " W5: Permanent grinding cover not mounted"
The firm stone cover has been dismounted. Remount the stone cover before restarting the preparation module.
- 6** " W6: Waste alcohol jar full "
The waste alcohol jar is full and needs emptying.
- 7** " W7: Unable to deliver specimen holder "
Unable to deliver specimen holder
The transportation arm cannot deliver the specimen holder as a another specimen holder is already in the lifting cup thus blocking the deliverance (Warning 7 only applies to a single unit without conveyors).
- 15** " W15: Specimen holder screws too long "
The specimen holder does not fit the cup during coupling. The screws for fixing the specimens are too long. Replace with shorter screws.
- 25** " W25: Grinding stone worn down "
The process has stopped due to worn down grinding stone. Mount a new stone.
The dresser depth sensor "Stone worn down" has been activated as initial state for a down-command.
- 30** " W30: Step can not be startet, no LVDT "
A process step in "removal mode" has been attempted while the parameter "LVDT installed" was "0" in the configuration. Set the parameter to "1" if a LVDT is mounted in the preparation module.
- 33** " W33: Stopped by thermal relay F31/F32 "
The preparation module has stopped or will not start due to an overloaded specimen holder motor. Wait for a couple of minutes and then try to restart the preparation module.
If this method does not work or the WARNING 33 frequently appears, call a service technician.

- 34** " W34: Feed conveyor empty "
There are no more specimen holders on the feed conveyor. Place a new specimen holder on the conveyor and restart the preparation module.
The sensor at the end of the feed conveyor has not been activated within the fixed period of time.
- 35** " W35: Stopped by thermal relay F21/F22 "
The preparation module has stopped or cannot be started due to an overloaded disc motor on P1. Wait for a couple of minutes and then try to restart the preparation module.
If the above does not work, or the WARNING 35 frequently appears, call a service technician.
- 36** " W36: Stopped by thermal relay F23/F24 "
The preparation module has stopped or cannot be started due to an overloaded disc motor on P2. Wait for a couple of minutes and then try to restart the preparation module.
If the above does not work, or the WARNING 36 frequently appears, call a service technician.
- 37** " W37: Stopped by thermal relay F33 "
The preparation module has stopped or cannot be started due to an overloaded high pressure pump. Wait for a couple of minutes and then try to restart the preparation module.
If the above does not work, or the WARNING 37 frequently appears, call a service technician.
- 38** " W38: Stopped by thermal relay F16 "
The preparation module has stopped or cannot be started due to an overloaded conveyor (feed) motor. Wait for a couple of minutes and then try to restart the preparation module.
If the above does not work, or the WARNING 38 frequently appears, call a service technician.
- 39** " W39: Stopped by thermal relay F17 "
The preparation module has stopped or cannot be started due to an overloaded conveyor (delivery) motor. Wait for a couple of minutes and then try to restart the preparation module.
If the above does not work, or the WARNING 39 frequently appears, call a service technician.

40 " W40: Stopped by thermal relay F15 "

The preparation module has stopped or cannot be started due to an overloaded recirculation pump. The thermo relay does not have an automatic reset function, and it is therefore necessary to call a service technician.

43 " W43: No stone, or stone worn down "

The dresser has reached the "stone worn" in connection with an adjustment procedure. The stone has worn too thin to be accepted or no stone has been mounted.

The dresser depth sensor "dresser down on stop switch" is activated by the down movement and not by the turning movement of the "stop nut".

MESSAGES

The message numbers are internal codes not visible to the user.

- 1 **"Recirc. water low, press STOP or CONTINUE"**
The water level in the recirculation tank is too low. Refill with water. The preparation module does not stop when this message appears but goes into a waiting state. The process may be continued by pressing CONTINUE, or you may stop and refill/change the water.
- 2 **"Soap level low, press STOP or CONTINUE"**
The soap jar behind the recirculation tank is nearly empty. Refill with soap. The preparation module does not stop when this message appears but goes into a waiting state. The process may be continued by pressing CONTINUE, or you may stop and refill with soap.
- 3 **"Alcohol level low, press STOP or CONTINUE"**
The alcohol jar behind the recirculation tank is nearly empty. Refill with soap. The preparation module does not stop when this message appears but goes into a waiting state. The process may be continued by pressing CONTINUE, or you may stop and refill with alcohol.
- 4 **"Stopped by emergency stop"**
The preparation module has stopped after pressing the emergency stop or cannot be started as the emergency stop is still depressed.
- 5 **"Stopped by upper guard"**
The preparation module has stopped due to one of the guards being open during the process or it cannot be started due to one or more guards being open.
- 6 **"Stopped by emergency stop & upper guard"**
This message is analogous with MESSAGE 4 and 5.
- 7 **"Stopped by lower guard"**
This message is analogous with MESSAGE 4 and 5.

- 8 **"Stopped by emergency stop & lower guard"**
This message is analogous with MESSAGE 4 and 5.
- 9 **"Stopped by upper and lower guard"**
This message is analogous with MESSAGE 4 and 5.
- 10 **"Stopped by e. stop, upper & lower guard"**
This message is analogous with MESSAGE 4 and 5.
- 11 **"Guards closed and emergency stop released"**
This message is a receipt on releasing the emergency stop, pushing the recirculation tank back or closing the guards.
- 12 **"Refill Px: A,B,C, press STOP or CONT."**
The bottle in position Px (1-9), bottle no. A,B,C (11-3, 21-3, 31-3) is empty and should be refilled.
- 13 **"Not in process mode, method ignored"**
A method has been sent to the MAPS module while the module was set to Development mode. Switch to Process mode and try to transmit the method from the PC again.
- 31 **"Process not started: method empty"**
The preparation module cannot be started, as the method chosen does not contain any executable steps. (The message is only shown on a single unit without conveyors).
- 32 **"Command ignored: process in progress"**
The start command is ignored when the preparation module has already been started.
- 33 **"Process queue empty"**
The preparation module has stopped as the process queue is now empty.
- 37 **"Unable to clean: no sample holder"**
Manual cleaning cannot be started without a specimen holder in the transportation arm. Place a specimen holder on the feed conveyor and let the preparation module fetch the holder. Then start the cleaning process.

38 "New consumables, press STOP or CONTINUE"

The consumables for the actual method do not match the cloth lubricant and/or suspension of the previous method on the working station.

The preparation module does not stop when the message appears but goes into a waiting state.

- Provided you have changed the consumables (most often the cloth) in advance, press CONTINUE to start the method.
- If you have not changed to the correct cloth, lubricant and/or suspension according to the method, you must do so before pressing CONTINUE.

Note: Besides starting the preparation module, CONTINUE registers the consumables of the new method. The information is used to check the next methods for possible mismatch.

A feature mainly of interest to the users of the single unit without conveyors: it is possible to change a cloth while the preparation module performs a cleaning process. You can thus avoid the mismatch MESSAGE 38 by pressing CONTINUE while the cleaning station is activated, thereby acknowledging the change of cloth in advance and ensuring that the preparation module does not stop.

40 "Command ignored: no stone in step no. 1"

Manual dressing of the stone cannot be performed due to one or more of the following causes:

- The previous process on POS 1 was not with a grinding stone.
- The actual method does not indicate a grinding stone in step 1.
- The dressing depth of the manual dressing function is zero.

42 "Process started"

If the user presses START or CONTINUE this message appears.

43 "Process stopped"

If the user presses STOP this message appears

44 "Process already stopped"

If the user presses STOP when the preparation module has already been stopped, this message appears.

45 "Predosing started"

If the user presses PREDOSE, this message appears.

46 "Manual process started"

This message appears when one of the manual functions of the screen picture MANUAL FUNCTIONS is activated.

47 "Stone worn out, change it"

The dressing depth sensor has detected "stone worn" depth after dressing the stone. For safety reasons you must not use the stone anymore. The method will not be stopped immediately but next time a grinding or dressing process is attempted, the preparation module will stop with a WARNING 25.

8. Spare Parts

Ask your Struers dealer
for information about MAPS spare parts.

9. Technical data

Preparation module Electrical Data

Module type	Voltage V / Hz	Current consumption in A
<i>Plane grinding</i>	3 x 200 V / 50 Hz	37.5
	3 x 200 V / 60 Hz	36.0
	3 x 220 V / 50 Hz	35.0
	3 x 220 V / 60 Hz	33.5
	3 x 380 V / 50 Hz	20.0
	3 x 415 V / 50 Hz	19.5
	3 x 440 V / 60 Hz	18.5
	3 x 480 V / 60 Hz	16.5
<i>Grinding/ polishing</i>	3 x 200 V / 50 Hz	21.0
	3 x 200 V / 60 Hz	21.0
	3 x 220 V / 50 Hz	20.5
	3 x 220 V / 60 Hz	20.0
	3 x 380 V / 50 Hz	12.0
	3 x 415 V / 50 Hz	11.0
	3 x 440 V / 60 Hz	11.0
	3 x 480 V / 60 Hz	9.5

Mains Cable Recommendation	<i>Voltage/frequency:</i>	<i>Min. Fuse size</i>	<i>Minimum cable size @ Min. fuse</i>	<i>Max. Fuse size</i>	<i>Minimum cable size @ Max. fuse</i>
<i>Plane grinding</i>	3 x 200 V / 50 Hz	32	3x2,5mm ² + PE	32	3x2,5mm ² + PE
	3 x 200 V / 60 Hz	30	3xAWG12 + PE	30	3xAWG12 + PE
	3 x 220 V / 50 Hz	32	3x2,5mm ² + PE	32	3x2,5mm ² + PE
	3 x 220 V / 60 Hz	30	3xAWG12 + PE	30	3xAWG12 + PE
	3 x 380 V / 50 Hz	16	3x1,0mm ² + PE	25	3x1,0mm ² + PE
	3 x 415 V / 50 Hz	16	3x1,0mm ² + PE	25	3x1,0mm ² + PE
	3 x 440 V / 60 Hz	15	3xAWG14 + PE	25	3xAWG14 + PE
	3 x 480 V / 60 Hz	15	3xAWG14 + PE	25	3xAWG14 + PE
<i>Grinding/polishing</i>	3 x 200 V / 50 Hz	16	3x2,5mm ² + PE	32	3x2,5mm ² + PE
	3 x 200 V / 60 Hz	15	3xAWG12 + PE	30	3xAWG12 + PE
	3 x 220 V / 50 Hz	16	3x2,5mm ² + PE	32	3x2,5mm ² + PE
	3 x 220 V / 60 Hz	15	3xAWG12 + PE	30	3xAWG12 + PE
	3 x 380 V / 50 Hz	10	3x1,0mm ² + PE	25	3x1,0mm ² + PE
	3 x 415 V / 50 Hz	10	3x1,0mm ² + PE	25	3x1,0mm ² + PE
	3 x 440 V / 60 Hz	10	3xAWG14 + PE	25	3xAWG14 + PE
	3 x 480 V / 60 Hz	10	3xAWG14 + PE	25	3xAWG14 + PE
	Important: Local standards may overrule the recommendations for the main supply cable. If necessary, please contact a qualified electrician to verify which option is suitable for the local installation setup.				
Residual Current Circuit Breaker	type B, 30 mA (or higher) is REQUIRED.				

Connections

Power

3-phase, consumption: see table above

Compressed air

6 bar (90 psi) - consumption approximately 200 l/min

Water

Main water input max. 10 l/min, 1-9.9 bar, outlet to normal drain

Computer (PC)

RS485 multidrop communication line

Main data

Height	Approximately 1905 mm
Length	Approximately 1535 mm
Depth	Approximately 1115 mm
Weight	Approximately 850 kg
Noise level	68 dBA
Safety standard	IEC 204 / EN 60204-1 (VDE 0113)

Special demands

- The module must be placed on an even floor
- The room temperature must be 15°C to 30°C
- Exhaustion recommended capacity max. 1000-2000/min

Specimen holder discs

160 mm (plus conveyor ring)
200 mm (plus conveyor ring)

Transportation arm

Motor for transportation arm	150 rpm, 0.37 kW (0.5 HP)
Operating pressure	50-700 N in 10 N steps, automatic adjustment in the phases

Plane grinding work station

Main motor	5.5 kW (7.4 HP)
Rotational speed	1450 rpm
Grinding stone	Ø 356 mm Grinding width 115 mm
Dressing	Automatic dressing of grinding stone
Recirculation cooling unit	60 l/min Container capacity 65 l

Fine grinding and polishing work station

Main motor (2 speeds)	0.55 kW (0.7 HP) / 1.1 kW (1.5 HP) Rotational speed (2 speeds): 150/300 rpm
Main motor (variable speed)	1 kW Rotational speed (variable speeds): 60-600 rpm
Grinding/polishing disc	Ø 300-350 mm

Cleaning station

Cleaning time User-defined
Cleaning method Water, detergent, alcohol, air
Cleaning programs 10 programs each comprising 10 steps.

Built-in computer facilities

Display LC display (13 lines of 42 characters)
Data storage 30 methods of 10 steps each

*External computer facilities
(accessory)*

Environment DOS
Hardware 80286 AT or higher
Database editor Load and save method from MAPS Master on a PC
Control functions Total control of all functions from a PC

Method programming

Steps 10 per method group
Process time 00:00:01 to 00:10:59 (H:M:S)
Force 50-700 N in 10 N steps
Rotation Direction: co-rotation or counter rotation
Speed: 150/300 rpm (2 speeds, fine grinding/polishing)
Speed: 60-600 rpm (variable speed, fine grinding/polishing)

Liquids Water: ON/OFF
Diamond suspension dosing: 20 levels or OFF
Lubricant dosing: 20 levels or OFF

Cleaning 10 possible cleaning programs per preparation module
Specimen holder disc The transportation arm position depends on the specimen holder disc diameter and the preparation disc diameter

Text Comments and definition names

Conveyor belt

Belt connections

Power 3-phase
Cabling To preparation module

Belt main dimensions

Length Approximately 1650 mm
Height/Width For specimen holder discs of 160/200 mm with conveyor ring
Weight Approximately 30 kg

General data, Conveyor belt

Motor 0.18 kW
Speed 0.13 m/s
Control Automatic drive controlled by the preparation module

Belt capacity 6 specimen holder discs

Recirculation cooling unit

Cooling unit capacity

Motor 0.27 kW
Water tank capacity 65 l
Weight 20 kg (empty)
85 kg (full)

Dosing system

Bottle unit capacity

DP-Suspension: 250 ml
DP-Lubricant: 1 l
OPS/OPU Suspension: 1 l

Alcohol separation unit capacity 5 l

MAPS Master, PC kit

Hardware

Communications adapter with serial RS232 to RS485 converter box and cable

Software

Method database/editor and communications program

Database capacity

8,000 preparation methods

Code numbering of MAPS modules

Work station codes

Each work station in the MAPS preparation module has an identification code (position number). The left/right station concept makes it easier for the user to understand the work station numbering system.

Module No.	Left station Pos. No.	Right station Pos. No.
1	P2	P1
2	P4	P3
3	P6	P5
4	P8	P7
5	P10	P9

Dosing system/Bottle unit codes

Each work station (polishing discs) can be equipped with an optional dosing system, including 1-3 optional bottle units. Each Dosing system has room for 3 bottle units. 2 types of bottle unit are available:

DP-Unit 3 bottles (1 lubricant and 2 suspensions)

OP-Unit 2 bottles (2 oxide suspensions)²

The individual bottles in the bottle units are identified by the pos. no. of the work station followed by a 2-digit code number:

1st digit The bottle unit number


2nd digit The bottle number in the Bottle unit

Example: **P4 12**: bottle unit 1, Bottle No. 2 of the left Dosing system in MAPS preparation module no. 2

² The OP-Unit is solely designed for use in Bottle unit 3 in the Dosing system placed at the left work station!

Bottle unit numbering

The drawing shows the positions and the numbering system of the 1-3 bottle units within the Dosing system.


Construction descriptions

	Table of contents	Page
General description	Specimen holder disc.....	78
	Transportation arm.....	79
	Plane grinding station.....	79
	Dresser for stone.....	80
	Fine grinding and polishing station.....	80
	Cleaning station	81
Optional units	Recirculation cooling unit	82
	Alcohol separation unit.....	83
	Dosing system.....	83
	Dosing system level indicator.....	83
	Conveyor belt	84
	Drying fan.....	85

Specimen holder disc

Struers specimen holder discs are designed for simultaneous preparation of several specimens fixed at the same level. The coupling stud for raising and lowering the disc has a ball joint which allows the stud to be moved up to 10° in vertical position. The position of the specimens in the specimen holder disc itself compensates for minor irregularities in specimens and/or preparation disc.

Collection


The specimen holder disc is attached to the transportation arm with a snap-coupling. The stud in the middle of the specimen holder disc has a slot which fits the snap-coupling. 3 balls in the snap-coupling slide into the slot during mounting to snap and lock the specimen holder disc.

After locking the stud, 3 driving pins slide into the 3 corresponding holes in the specimen holder disc.

If the snap-coupling fails to snap the stud of the specimen holder disc, a sensor activates an error code in the microprocessor and the module stops.

Fastening screws


The specimen holder disc is placed on the conveyor belt. When MAPS starts, the specimen holder disc is transported to the collection position. The transportation arm moves above the specimen holder disc and a ring raises the specimen holder discs to control the position of the fastening screws. The MAPS preparation module automatically stops if the screws protrude .

WARNING!

The fastening screws of the specimen holder disc may damage the cleaning station if they protrude!

Always remember to use screws with an appropriate length!



Delivery

When the specimens have been processed by MAPS, the transportation arm moves the specimen holder disc to the delivery position. The specimen holder disc is lowered onto the conveyor belt. A pneumatic device locks the conveyor ring of the specimen holder disc to the conveyor belt. The snap-coupling then unlocks the stud of the disc while moving away from the conveyor belt.


The single unit without conveyors still feeds and delivers the specimen holders similar to the standard MAPS preparation modules.

Transportation arm

The transportation arm of MAPS moves the specimen holder discs. It is also equipped for grinding and polishing purposes. The outer part of the transportation arm with the snap-coupling moves up and down, driven by a powerful pneumatic cylinder.

The transportation arm turns around its vertical axis, from the feeding position, over the grinding/polishing and cleaning stations and onto the delivery position. A sensor constantly monitors the position of the transportation arm, in order to place the specimen holder discs correctly.

The force of the specimen holder disc against the preparation disc is adjusted by changing the pressure of the pneumatic cylinder.

The specimen holder disc is rotated by an asynchronous motor with gear. The motor is placed behind the transportation arm and transmits the force to the specimen holder disc through a tooth belt.

A disc brake mounted below the transportation arm locks the transportation arm to keep it steady during the preparation process.

See also Section 10, Method Parameter values.

The plane grinding work station is equipped with a grinding stone with automatic diamond dresser and a Recirculation cooling unit.

The bearing unit consists of 2 heavy ball bearings, a roller bearing and an angle bearing, ensuring a firm and quiet movement.

The motor is a powerful asynchronous 5.5 kW motor. Power is conveyed to the grinding stone by 3 V-belts. The rotational speed of the stone is 1450 rpm.

The stone is surrounded by heavy protection cover against loose specimens. 4 finger screws are all that are needed to remove the cover when the stone is ready to be replaced. A safety switch prevents the stone from rotating without its cover. A spring lifts the cover above the safety switches provided that the finger screws are not fastened properly.

Dresser for stone


To ensure perfect working conditions of the grinding stone, a diamond dresser has been mounted on the transportation arm. The diamond dresser is fastened to a compressed air cylinder with a mechanical stop. After each grinding operation, the arm turns to place the diamond dresser against the edge of the stone. The stone starts rotating and the piston rod with the diamond tool is lowered to a predefined depth. The arm (and the dresser) swings over the grinding stone and removes a thin layer of the surface.

The multi-point dressing diamond of the dresser dispenses the diamond grains from the cylinder. During the dressing several diamonds attack simultaneously, providing an even surface. The diamond cylinder is used to its full extent, with worn-out diamonds falling out of their own accord and being replaced by new diamonds. The adjustment of dressing depth is done by a step motor.

The diamond dresser automatically checks the actual height of the stone, when booting the MAPS module or when the stone cover has been dismantled. A sensor detects whether the stone has completely worn down and must be replaced.

Fine grinding/polishing station

See also Section 2, Operation, Change of grinding stone

The fine grinding/polishing work station features built-in water intake and a dosing system for automatic dosing of lubricant and diamond suspension.

The station is built on the same type of ball bearings as the grinding station. The motor may be a 2-speed motor or a variable speed motor. The power from the motor is conveyed to the turntable by 2 V-belts.

<i>Cleaning station</i>	The cleaning station is placed between the 2 work stations. The station is made of stainless steel. The fully programmable cleaning station is provided with the following nozzle types:
Soap	The specimen holder disc is cleaned underneath with a soap solution from 3 soap nozzles of the ejector type. Compressed air sucks the soap from the soap jar. The compressed air is filtered to avoid dirt and oil in the soap spray.
Alcohol	The cleaning station has 3 alcohol ejector nozzles to clean and facilitate the drying of specimens. As in the case of soap, the alcohol spray is filtered for dirt and oil, in order avoid damaging the specimens.
High pressure water	A vane pump, driven by an asynchronous motor, increases the water pressure to 10 bars, distributing the water through 6 nozzles. 5 of the nozzles clean the underneath of the specimen holder disc, while the 6th nozzle cleans the side.
Low pressure water	A low pressure nozzle washes the upper side of the specimen holder disc to remove dirt between conveyor ring and specimen holder disc. 2 low pressure nozzles simultaneously wash the side and the upper part of the conveyor ring.
Air drying	In the cleaning station there are 5 nozzles for air drying using cold compressed air. As in the case of soap and alcohol, the air is filtered for dirt and oil. 3 of the nozzles blow on the underneath of the specimen holder disc and the 4th nozzle blows on the upper side. The 5th nozzle dries the conveyor ring.
Diamond dresser	<p>A MAPS preparation module for plane grinding (PLAN/POL) is provided with a diamond dresser (the models MAPSE or MAPON for plane grinding). A nozzle for this special purpose cleans the dresser with air after each dressing process.</p> <p>Level indicators in the soap and alcohol jars, placed on the trolley behind the recirculation tank, warn when there are only 15 mm soap or alcohol left in the jar.</p>

Recirculation cooling Unit
The Trolley

The trolley delivered with the plane grinding and fine grinding/polishing preparation module (Code word: *MAPSE*) consists of a water tank, a motor with a centrifugal pump and a sedimentation system for purification of water.

A level indicator checks the amount of cooling water in the tank.

The Recirculation cooling Unit

The Recirculation cooling unit cools the plane grinding disc during method processing. One or more fine grinding/polishing stations may be equipped with optional Recirculation cooling units as well.


Recirculating Connection Kit

This optional kit connects both preparation modules to a shared Recirculation cooling unit. The kit consists of a cross-over valve for inlet from the Recirculation cooling unit and a cross-over valve from the spill pan in the work stations to a Recirculation cooling unit or the drain.

Remember...

Protect both the module and your specimens by using Struers Anti-corrosion additive in the cooling water!

Trolley with Recirculation cooling unit (MAPRE)


Alcohol separation unit

To protect the local environment, we recommend the use of an optional Alcohol separation unit. The unit consists of a third jar to be placed in the trolley beside the soap and alcohol jars. With an Alcohol separation unit the surplus alcohol from the cleaning station is led to a separate alcohol jar to avoid contamination of the public sewer systems.

The Alcohol separation unit is equipped with a max. level indicator.

Dosing system

The MAPS preparation modules are delivered with custom designed dosing systems. The dosing systems consist of a metal box with plastic cover and 1-3 bottle units which in turn come in 2 versions, for lubricant/diamond or OPS/OPU suspension.


Bottle unit types	
Lubricant/ diamond suspension	OPS/OPU suspension
<p>OP-Unit</p> <p>1459_028.CDR</p>	<p>OP-Unit</p> <p>1459_028.CDR</p>

See also the previous section of Bottle unit code numbers.

Dosing system level indicator

Level indicators (accessories) may be mounted beneath the bottle units to detect whether the contents of the bottles drops below 15 mm. The indicators work according to a capacitive principle, and are not in contact with the liquid.

Conveyor belt

The optional MAPS conveyor belts automatically move the specimen holder discs between the MAPS preparation modules. Each conveyor module holds up to 6 specimen holder discs, allowing the belt to function as a buffer.

Up to 5 preparation modules may be interconnected with conveyor modules in each preparation module, combining the preparation units to form a single fully automated system. The conveyor belts may be mounted both in parallel and at a right angle to the preparation module to connect the individual modules, suiting the dimensions and arrangement of the laboratory room.

The specimen holder discs are transported by a conveyor ring on 2 rubber belts. The conveyor belts are controlled by the preparation module. If more preparation modules are interconnected, the first preparation module controls the feed module and the first connection module. The next preparation modules only control the attached delivery module.

Feed module


The conveyor belt feeding the first preparation module is equipped with a sensor verifying that the specimen holder disc actually reaches the module.

If no specimen holder disc arrives when START \diamond is activated, the conveyor belt will run until a specimen holder arrives or until the belt has completed its full circuit.

Connection module

The connection module between the preparation modules is further equipped with a sensor, which is able to detect whether there is sufficient space to deliver another specimen holder disc. Only the delivering preparation module can activate the conveyor belts with the specimen holder discs. A sensor placed at the next preparation module checks that a specimen holder disc is ready for processing. If no specimen holder discs arrive, the preparation module is on standby until the previous preparation module delivers a specimen holder disc.

Delivery module

The delivery conveyor belt has a sensor to confirm that there is free space for a new specimen holder disc. If there is enough space, the transportation arm delivers the specimen holder disc, to be transported out of the preparation module on the delivery conveyor belt.

An optional drying fan may be installed beneath one or more of the delivery Conveyor modules. A grate replaces the first part of the conveyor belt. The drying fan is mounted below the grate. When a specimen holder disc arrives the drying fan starts drying for a preprogrammed amount of time (See Section 4, Software maintenance). While drying, the conveyor belt pauses until the specimen holder disc on the grate has been processed.

The drying fan may also be mounted on the special single unit without conveyors version.

Safety


MAPS is equipped with approved safety guards for protection of the operator. The grinding and polishing stations only work when these guards are closed. The safety guards cannot be opened unless the specimen holder disc is in the cleaning station or the preparation module has been stopped; in all other cases the emergency stop will be activated.

Exhaustion

Alcohol, which is hazardous to breathe, is used during polishing and cleaning. MAPS is prepared for connection to an exhaustion system. The connecting plastic collar is placed on the top of MAPS.

Emergency stop

The emergency button stops the module and the attached conveyor belts, both input and output. Other preparation modules connected to the modules you have stopped are not affected directly, but will stop automatically when they run out of specimen holder discs.

Other safety features

- Safety cover over the grinding stone
- Safety guards around the work stations
- Recirculation cooling unit on trolley
- Safety switches on the doors

10. Method parameter values

The parameters of a preparation method can be divided into 3 main groups: method parameters, step parameters and phase parameters.

- A method in the MAPS system consists of 3 method parameters and 10 process steps.
- Each process step consists of 9 step parameters and 3 process phases.
- The process phases are named pre-phase, main phase and final phase, each in turn consisting of 5 phase parameters.

A method has a total of 243 parameters ($3 + 10 \times (9 + 3 \times 5) = 243$).

Method parameters

Method

Method number in the module database

A value between 1 and 30, i.e. the database of the module contains a total of 30 methods. If you type a method number in any screen (except for the Process Control screen), the whole parameter set of the chosen method is loaded into active database memory. The method number is the parameter which identifies the whole parameter set of the method.

Description

Method description

A text string of max. 28 characters, e.g. the name of the method. This text can only be edited from the MAPS Master PC database.

X,Y,Z

Method number in the MAPS Master PC database

A special 3 digit method number in parentheses indicates the method number of the MAPS Master PC database. The format is [X,Y,Z], similar to the MAPS Master format of main group (X), sub group (Y) and method (Z).

This parameter in the MAPS module screen is only visible if the MAPS Master PC database is connected to the module. It is not possible to edit the MAPS Master method number from the MAPS preparation module.

Method parameters, method level

```

Method: 20* ( 1, 1,20)  METHOD DATABASE
Description: Rails

STEP POS DISC      #/µm SUSP    LUB      T/µm
  1   1  Stone      #60
  2   1  Stone      #60
  3   2  DP-Plan          15µm  Blue     5:05
  4   3  DP-Plus          3µm   Blue     5:00
↓  5   4  DP-Nap

```

Type new parameter value

```

F1:EDIT    F2:          F3:COPY    F4:
F5:         F6:          F7:DELETE  F8:

```

Method parameters, step level

```

Method: 20    Step: 1  METHOD EDIT
POS DISC      #/µm  SUSP    LUB      T/µm
  1  Stone      #60

```

```

Time   :  --:--  --:--  --:--  Disc RPM :1500
S. rem:    25   100    0  Direction:→→
Force  :    200  500    0  Cleaning  :  0
Susp.  :     0    0    0  Dressing  :  20
Lub.   :     ON   ON   ON

```

Type new parameter value

```

F1:         F2:          F3:INSERT  F4:
F5:         F6:          F7:REMOVE  F8:

```

Step parameters

POS

Work position

A value between 0 and 10 refers to the work station attached to the process step. If the value is 0, the step will be inactive.

Grinding/polishing disc

A value between 0 and 59 refers to the disc description of the step. The value is used as an index key to a text table with all the disc descriptions. The descriptions may be redefined from the MAPS Master PC database.

Each disc description refers to 2 other parameters; a disc type and a disc diameter. The disc type may only be redefined from the MAPS Master PC database.

Disc type: is used by the machine to identify the disc type. This helps prevent accidental dressing of diamond discs and the like.

Disc diameter: is used to determine the position of the transportation arm during the process. The diameter of the disc may be changed in the configuration menu, Disc Sizes.

The main purpose of choosing a disc description is to decide cloth for the step and the desired position of the transportation arm.

/ μm

Grain/grit size of disc

A value between 0 and 49 refers to a grit description for the step. The value is used as an index key to a text table with all the grit descriptions. The descriptions may be redefined from the MAPS Master PC database. The parameter is only descriptive and has no physical functions in the module.

SUSP

Suspension

A value between 0 and 15 refers to a suspension description for the step. The value is used as an index key to a text table with all the suspension descriptions. The descriptions may be redefined from the MAPS Master PC database.

Each suspension description refers to 3 other parameters: 1 suspension type and 2 bottle number. The suspension type may only be redefined from the MAPS Master PC database.

Suspension type: defines the type of suspension (diamond or oxide). When you choose oxide suspension, water will be used after finishing the preparation method.

Suspension bottle number: defines the physical position of the 2 suspension bottles of each dosing system. The bottle numbers may be changed in the configuration menu, Suspensions.

The suspension description simply tells the module where to find the suspension bottle in the dosing system.

LUB

Lubricant

A value between 0 and 15 refers to a lubricant description for the step. The value is used as an index key to a text table with all the lubricant descriptions. The descriptions may be redefined from the MAPS Master PC database.

Each lubricant description refers to 3 other parameters: 1 lubricant type and 2 bottle numbers. The lubricant type may only be redefined from the MAPS Master PC database.

Lubricant type: defines the type of lubricant (mains water, recirculating water or lubricant from a bottle).

Lubricant bottle number: defines the physical position of the 2 lubricant bottles of each dosing system. The bottle numbers may be changed in the configuration menu, Lubricants.

The lubricant description simply tells the module where to find the lubricant bottle in the dosing system.

RPM

Rotational speed of the disc

A value of 150, 300 or 1500 defines the rotational speed of the disc. If the value you type differs from the 3 predefined values, the value will be adjusted to the nearest predefined value. The rotation speed of a plane grinding work station is always 1500 rpm.

Direction

Rotational direction of the specimen holder

A value of 0 or 1 defines the rotational direction of the specimen holder in relation to the disc.

Rotation direction	
Value	Symbol/explanation
0	←→ (counter rotation)
1	→→ (co-rotation)

Cleaning

Cleaning program number

A value between 0 and 10 defines the cleaning program to be performed after a method step.

Cleaning program (step parameter)	
Value	Explanation
0	No cleaning after the method step
1	Cleaning program 1 <u>including</u> cleaning of the dresser tip
2-10	Cleaning program 2 to 10

All 10 cleaning programs may be edited in the configuration menu, Cleaning programs. New cleaning programs may also be defined from the MAPS Master PC database.

Dressing

Dressing depth

A value between 0 and 400 defines the dressing depth of the stone in μm after each method step. The dressing depth is subdivided in steps of 75 μm . The parameter is only available if the disc type of the step is "Stone".

Dressing depth	
Value	Repeat dressing
0-74	1 time
75-149	2 times
150-224	3 times
225-299	4 times
300-374	5 times
375-400	6 times

Phase parameters

Time

Process time

2 values defining the timing of the process phase in minutes and seconds.

Minutes: between 0 and 10 min.

Seconds: between 0 and 59.

If a stock removal value (see below) other than 0 has been chosen, the process timer will be inactive ("---:---").

S. rem

Stock removal

A value between 0 and 2000 defines the amount in μm to be removed from the specimen (0 to 2.0 mm). The process is said to run in *stock removal mode* if the total stock removal is higher than 0 (pre-phase + main phase + final phase). If stock removal is 0, the process runs in *time mode*. You cannot use both Process time and Stock removal in the same method step.

Phase modes		
Mode	Process Time	Stock removal
Time mode	XX:YY	0
Stock removal mode	---:---	(1-2000)

Force

Force

A value between 0 and 700 Newton defines the vertical force on the specimen holder during the process. The force is divided into steps of 10N each. Value between the fixed steps is rounded down to the nearest step (e.g. 257N down to 250N).

Susp / Lub.

Dosing level for suspension (Susp) and lubricant (Lub.)

A value between 0 and 20 defines the suspension/lubricant dose per time unit. The scale has been calculated to increase/decrease the dosing level differently depending on the suspension/lubricant type:

Type	Approximate change of dosing* in %	Active dosing	Pause between dosings
Diamond suspension	18%	0.3 sec	2.7-74.7 sec
OPS/OPU suspension	17%	0.15 sec	0.85-19.85 sec
Lubricant**	21%	0.15 sec	0.35-19.85

* The percentage is the increase/decrease each time you change from one level to another.

** If the lubricant type is mains water or recirculating water, the numerical dosing parameter changes to the text "ON" for continual dosing.

The following table shows the percentage of dosing time for each of the 20 levels:

Dosing level in percentage			
Level	Dosing time in % of process time		
	Diamond sus- sension	OPS/OPU sus- sension	Lubricant in bottle
0	0.00%	0.00%	0.00%
1	0.40%	0.75%	0.75%
2	0.47%	0.88%	0.91%
3	0.56%	1.0%	1.1%
4	0.67%	1.2%	1.4%
5	0.79%	1.4%	1.6%
6	0.93%	1.7%	2.0%
7	1.1%	1.9%	2.4%
8	1.3%	2.3%	2.9%
9	1.6%	2.7%	3.5%
10	1.8%	3.1%	4.3%
11	2.2%	3.6%	5.3%
12	2.6%	4.2%	6.4%
13	3.1%	5.0%	7.7%
14	3.6%	5.8%	9.4%
15	4.3%	6.8%	12.0%
16	5.1%	7.9%	14.0%
17	6.0%	9.4%	17%
18	7.1%	11%	20%
19	8.5%	13%	25%
20	10%	15%	30%

Positioning the transportation arm during the process

The stop position of the transportation arm moving towards the work station is defined by 2 parameters: the specimen holder diameter and the disc diameter defined in the method step. The disc diameter may be changed in the configuration menu. Disc sizes.

Automatic transportation arm position: the specimen holder periphery is automatically positioned to be flush with the periphery of the disc (see sketches). This position guarantees that the specimens remain inside the rim of the disc.

Forced transportation arm position: in the case of Diamond Pad or Petrodisc-M, however, the specimens should sometimes exceed the rim of the disc, in order to avoid uneven wear of the disc. You have 2 possible procedures:

- ① Define a disc diameter larger than the actual diameter of the disc.
- ② Define a specimen holder diameter less than the actual diameter of the specimen holder.

Change transportation arm position permanently: both of the above procedures forces the transportation arm to stop at the desired position. If you want the transportation arm position changed permanently, in connection with certain disc types, change the disc diameter of the disc type according to procedure ①. The specimen holder will always be positioned towards the far edge of the disc to ensure maximum room for the dosing jets.

Cleaning programs

A cleaning program is the process of cleaning the specimens and the specimen holder after a preparation step. The program itself consist of 10 sequential cleaning steps. Each cleaning step has 2 parameters, a cleaning medium and a process time:

Cleaning program parameters		
Parameter type	Value	Explanation
Cleaning media	0	No cleaning (all valves inactive)
	1	Flush with high pressure water
	2	Wash with mains water
	3	Soap
	4	Alcohol
	5	Compressed air underneath
	6	Compressed air on the upper side
Process time	0-255 sec	

You can provoke a pause in the cleaning process by setting the cleaning medium to 0 (no cleaning) and the process time to a duration larger than 0.

Similarly you can skip a cleaning step by setting the process time to 0.

The MAPS module database can handle 10 different cleaning programs, numbered 1 to 10. The cleaning programs are defined in the configuration menu, Cleaning programs. If you want to experiment with the cleaning programs we recommend that you use the empty programs, instead of overwriting the predefined (and tested!) Struers cleaning programs.

Note Cleaning program #1 also cleans the dresser tip, unlike the other cleaning program numbers.

Disc cooling

A value between 0 and 10 defines how much cooling water to use per time unit. During disc cooling, the cooling water is on for 10 sec. Depending on the cooling level there will then be a pause of between 0 and 35 sec. The table below shows the percentage of the process time during which the disc cooling (water on) is active:

Disc cooling level in percentage	
Level	Cooling water on in % of process time
0	0,0% (OFF)
1	22.2%
2	25.0%
3	28.6%
4	33.3%
5	40.0%
6	50.0%
7	57.1%
8	66.7%
9	80.0%
10	100.0%

11. Index

Abrasive	13
Additive cooling fluid	44
Air filters	44
Alcohol	80
Alcohol separation unit.....	82
ARROW	12
Auto stop	23
Batch	17
Bottle unit	75
Bottle unit codes.....	74
Bottle unit types.....	82
Bottle units	74
Changing grinding stone	39
Cleaning program.....	50
Cleaning station	8, 80
Collection	77
Communication	58
Configuration menu.....	46
Continue.....	13
Conveyor belt	8, 83
COPY	33
Database.....	9
DEL	12
DELETE	34
Delivery	77
Development mode	6, 22
Dosing system.....	74, 82
DP-Unit.....	74
Dresser.....	79
Dresser for stone.....	79
Drying fan.....	84
Edit	25
Emergency stop	12
ENTER	12
Error code	59
ESC.....	12
Exhaustion	84
Fastening screws	77
Fine grinding/polishing station.....	79
Function keys	12
Key	12
Level indicator	80, 81, 82
Lock system	13
Lubricant	13
Main menu	16
Manual functions	28
Manual predosing.....	14
MAPS Master PC program.....	58
Mechanical maintenance	42
	97

Menu structure	10
MESSAGES	65
Method	9
Method database	30
Method editor	9
Method number format.....	20, 31
Method parameter values	85
Method steps.....	8
Method number format.....	18
Module	5, 6
Module configuration.....	51
Nozzle plate	45
Nozzle types.....	80
OPS/OPU suspension.....	82
Option menu.....	26
OP-Unit	74
P-CTRL	20
PGDN/PGUP	12
Plane grinding station.....	78
Position	9
Position number	9, 74
Predose.....	13
Preparation methods	8, 41
Process	99
Process Control.....	22
Process Edit	22, 25
Process mode	6, 17
Process queue	17
Recirculating cooling unit	81
Safety	84
Safety switch	78
Security system.....	13
Soap.....	80
Software maintenance	46
Specimen holder disc.....	77
Start.....	13
Step parameters.....	9
Stop.....	13
Struers Metalog Guide_	41
Technical data	70
Transportation arm.....	78
Version	15
WARNINGS	62
Work station codes.....	74
Work stations	6

Quick Reference

Processing

Process finished methods

- Choose Process mode in the Main menu.
- Type a method number and press ENTER.
- Type the actual specimen holder disc diameter.
- Type the number of actual specimen holder discs to be prepared.
- Press START ◊.

Process a single or test method

- Choose Development mode in the Main menu.
- Type a method number and press ENTER.
- Type the actual specimen holder disc diameter.
- Press START ◊.

Editing

Edit a method

- Choose Method database in the Main menu.
- Type a method number and press ENTER.
- Select the desired step with arrow and press ENTER or F1=EDIT.
- Fill out all necessary *Method parameter values* for all the necessary *steps* in the method.

Edit a running method

- While running a method in Development mode, press F1=EDIT.
- Change the parameters.

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Herewith declares that

<i>Product Name:</i>	MAPS
<i>Type No:</i>	459
<i>Machine Type:</i>	Automatic Grinding and polishing machine

is in conformity with the provisions of the following directives:

Safety of Machinery 98/37/EEC according to the following standard(s):
EN ISO 12100-1:2005, EN ISO 12100-2:2005, EN ISO 13850:2006, EN 60204-1:2006.

EMC-Directive 2004/108/EEC according to the following standard(s):
EN 61000-6-1:2002, EN 61000-6-3:2002.

Low Voltage Directive 2006/95/EEC according to the following standard(s):
EN 60204-1:2006.

Supplementary Information The equipment complies with the American standards:
FCC part 15, subpart B and UL508.

The above has been declared according to the global method, module A

Date: 24.10.2008

Christian Skjold Heyde,
Vice President, R&D and Production, Struers A/S

Hersteller Struers A/S
Pederstrupvej 84
DK-2750 Ballerup, Danmark
Telefon +45 44 600 800

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Christian Skjold Heyde,
Vice- President, R & D et Production, Struers A/S

MAPS

**Fully automatic,
computer-controlled
preparation system**

- **Tailor-made solution**
- **High volume preparation**
- **Excellent preparation quality**
- **Outstanding reproducibility**
- **Considerable savings in manpower and consumables**
- **Very easy operation**

MAPS is the ultimate solution for automatic preparation of metallographic specimens. A modular system which is tailored exactly to the user's specific requirements. MAPS takes care of the entire preparation process - from plane grinding to final cleaning and drying. Preparation quality is excellent, and due to exact control of preparation processes and consumption of abrasives and lubricants, the reproducibility is second to none. Another consequence of controlled preparation

is considerable savings in preparation costs.

Modularity means flexibility. MAPS fits into a medium-sized lab, and it is an obvious choice for high capacity preparation. Whatever the requirements are, a MAPS system can be designed to meet the needs.

Modular Automatic Preparation Systems = MAPS

A Preparation Sequence - Step by step

MAPS consists of modules, each with two work stations, one cleaning station and a specimen holder arm. Accessories such as dosers, recirculation, hot-air dryers, level indicators for consumables, etc. differ according to the configuration of each MAPS system. Conveyers transport the specimen holders from unit to unit. Suspensions and lubricants are stored in easily accessible modules. Up to three modules can be connected to a preparation station, each module containing two suspension bottles and one lubricant bottle. Alternatively, an oxide polishing module for two oxide suspensions may be added.

More than 30 combinations

The MAPS modules can be configured in at least 32 different ways, depending on the requirements. This is the consequence of modularity taken to its logical extreme: each MAPS system is configured on the basis of a thorough analysis of the user's specific needs. However, a MAPS system can at any time be expanded by addition of one or more modules, if the need for larger capacity arises.

Very high efficiency

A steady flow of specimen holders can pass through MAPS. When one module takes care of polishing, another is busy plane grinding the next batch of specimens. A special queue function makes it possible to prepare different materials according to different preparation methods – fully automatically. The queue function ensures an extremely fast throughput and eliminates all bottleneck problems, even at high volume preparation.


When MAPS controlled modules, 30 preparations are available, 10 predefined and 20 user-methods

1. Six specimen holders queue up in the conveyer leading into MAPS

5. When the grinding is completed, the specimen holder is carried to the cleaning chamber

2. Specimen holders are fed into the first position - plane grinding - by conveyer belt

6. After fine grinding the specimen holder moves on to polishing

3. The specimen holders are collected from the conveyer belt, using a special lifting device

7. The finished specimen holder is moved to the conveyer belt, where it can be collected by the metallographer for microscopic examination

4. Grinding is now in operation. The next specimen holders are ready for collection

is
n the
on methods
d methods
defined

Reduce preparation costs


More time for demanding work

MAPS takes care of the entire preparation process, even when it comes to the most complex specimens, leaving the operator time for other chores. In many cases, materialographic preparation has been a bottle-neck. MAPS eliminates this problem and leaves plenty of time for the one task which has not yet been satisfactorily automated: microscopic examination and evaluation of a materialographic structure.

Reduce preparation costs with MAPS

Automation not only means unique reproducibility and high quality – it also implies large cost savings. Within the world of materialography it is gradually realized that the initial cost of automatic equipment is more than counterbalanced by savings in time, labour and consumables. Materialographic preparation employs costly diamond products. Due to precise, automatic control of consumables consumption, MAPS takes care that exactly the right amount is used – no more, no less.

Cost saving analysis: investment in a MAPS system, depreciated over a period of three years, compared to an existing Abra system.

PC-controlled operation

While MAPS is working, the operator may monitor the process from outside the production area via a PC

MAPS can be programmed directly on the separate modules. A maximum of 30 preparation programs may be stored in each module's memory. However, MAPS becomes an even more versatile tool when controlled by a PC. For this purpose, a PC-program, MAPS Master, has been developed. MAPS Master allows for a controlled batch processing and supervision of the preparation process. Also included in MAPS Master is a user's database with facilities for creating, printing and storing up to 8,000 individual preparation methods. MAPS Master controls procedures, dosing, process times and all movements of specimen holder arm and conveyers. Switching to a different preparation procedure is just a matter of selecting from the PC database and pressing start.

High capacity specimen production


Today, quality control is gradually moving from isolated laboratories into the production area. This makes new demands to the preparation equipment within the fields of reliability and safety. MAPS is prepared for a rough life – working around the clock in an industrial environment. Sensors constantly monitor the process, ensuring a smooth and uninterrupted preparation process. In case of shutdown the operator is specifically advised which consumables to refill, etc. While MAPS is working in the production area, the operator may monitor the process from an office via the PC. With MAPS as the foundation of current quality control, adjustments in the production can be carried out with minimal delay.

Safety


MAPS sets new standards for safety. With a totally automatic system every aspect of personal and machine safety has to be included. MAPS complies with the strictest international safety standards.

Two typical MAPS configurations

A two-module configuration for fully automated, high-volume preparation. The two modules are controlled by a PC and linked by transfer conveyers. The two work stations in Module 1 perform plane grinding, fine grinding and specimen cleaning. Module 2 carries out the two polishing steps, plus cleaning and drying. A queue function makes it possible to line up different specimen holders. It also eliminates stoppages: the transport system is non-synchronous, so the conveyers serve as accumulators of specimen holders.

A basic configuration with one preparation unit, supplemented by an AbraPlan-10 plane grinding machine. Once the specimen holders are loaded, the preparation process takes place automatically. Like any MAPS configuration, this system may later be expanded if the need arises.

Struers A/S
 Pederstrupvej 84
 DK-2750 Ballerup, Denmark
 Phone +45 44 600 800
 Fax +45 44 600 801
 struers@struers.dk

Technical data

Connections	Power: Compressed-air supply: Water supply: Standard PC:	Three-phase 6 bar (90 psi) – consumption approx. 50 l/min Mains water input, output to normal drain RS 485 multidrop communication line
Main dimensions	MAPS Preparation Module Height: Width: Depth: Weight:	1900 mm 1470 mm 1080 mm 850 kg
Specimen holders	Size:	160 mm dia. and 200 mm dia. Abra specimen holders with conveyer ring
Specimen holder arm	Specimen holder motor: Operating force:	150 rpm, 0.37 kW (0.5 HP) 50-700 N in steps of 50 N, automatic adjustment in the phases
Plane grinding work station	Main motor: Rotation speed: Grinding stone: Dressing: Recirculation cooling unit:	5.5 kW (7.4 HP) 1450 rpm Diameter 356 mm, grinding width 115 mm Automatic dressing of grinding stone 60 l/min. Container capacity 65 l
Fine grinding and polishing work station	Main motor: Rotation speed: Grinding/polishing disc:	550 W (0.7 HP)/1110 W (1.5 HP) 150/300 rpm Diameter 230-350 mm
Cleaning station	Programmable automatic cleaning:	- Cleaning time - High-pressure water rinsing - Cleaning with detergent - Cleaning with alcohol - Air-drying
Operating system	Memory capacity: Display:	30 preparation programs: 10 fixed; 20 open. 13-line LCD with 42 characters
Conveyor	Width: Length: Capacity:	For 160 and 200 mm dia. Abra specimen holders with conveyor ring 1650 mm 6 specimen holders
Dosing system	Bottle for DP-Suspension: Bottle for DP-Lubricant: Bottle for OP-Suspension:	250 ml 1000 ml 1000 ml
MAPS Master 2.0	Communications unit: PC-software Data base capacity: Requirements:	RS485 8,000 preparation programs PC running Windows XP

USA
Struers Inc.
 24766 Detroit Road
 Westlake, OH 44145-1598
 Phone +1 440 871 0071
 Fax +1 440 871 8188
 info@struers.com

CANADA
Struers Ltd.
 7275 West Credit Avenue
 Mississauga, Ontario L5N 5M9
 Phone +1 905-814-8855
 Fax +1 905-814-1440
 info@struers.com

SWEDEN
Struers A/S
 Smältvägen 1
 P.O. Box 11085
 SE-161 11 Bromma
 Telefon +46 (0)8 447 53 90
 Telefax +46 (0)8 447 53 99
 info@struers.dk

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Specifications

	Code
MAPS Preparation Module With 1 plane grinding station and one fine grinding/polishing station (150/300 rpm), with recirculating cooling unit (MAPRE) With 2 fine grinding/polishing stations (150/300 rpm),	MAPSE MAPFI
Conveyor feed module connection module delivery module	MAPEF MAPCO MAPLI
Recirculation Cooling Unit contents 65 l, with a rolling pallet	MAPRE
Recirculation Kit for connection of 2 work stations to one recirculation cooling unit	MAPIT
Cooling Arrangement for preparation disc	MAPOL
Dosing System box with room for 3 bottle units. Bottle units (MAPSU or MAPOP) required	MAPTH
Bottle Unit for dosing system, with room for 2 DP-Suspensions and 1 lubricant for dosing system, with room for 2 OP-Suspensions	MAPSU MAPOP
Level Indication Unit for monitoring the level of suspensions and lubricants in the bottles. Level indicators (MAPDU or MAPPO) required	MAPIN
Level Indicator for bottle unit with DP-Suspension (MAPSU) for bottle unit with OP-Suspension (MAPPO)	MAPDU MAPPO
Stock Removal Unit for measuring of the material removal during the plane grinding process and monitoring of abnormal vibrations	MAPOM
Separation Unit for collection of the waste alcohol from the cleaning station	MAPAL
Manual Feeding Station For MAPS without conveyor	MAPEX
Drying Fan for drying of specimen holders	MAPFA
Conveyor Ring for mounting on specimen holder, 160 mm dia. for mounting on specimen holder, 200 mm dia.	MAPEL MAPIR
MAPS Master 2.0 Software for operation and control of the MAPS preparation modules. Complete with modem for communication and PC-software with method database/editor. Requires Windows XP SP2 or higher.	MAPMO

Struers' equipment is in conformity with the provisions of the applicable International Directives and their appurtenant Standards. (Please contact your local supplier for details)

Struers' products are subject to constant product development. Therefore, we reserve the right to introduce changes in our products without notice.

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Struers GmbH