# Specimen holders

## User Guide

# Table of Contents

<b>1</b>	<b>English</b> .....	<b>3</b>
<b>2</b>	<b>German</b> .....	<b>11</b>
<b>3</b>	<b>French</b> .....	<b>21</b>
<b>4</b>	<b>Spanish</b> .....	<b>31</b>
<b>5</b>	<b>Japanese</b> .....	<b>41</b>
<b>6</b>	<b>Chinese</b> .....	<b>51</b>

# Specimen holders

## User Guide

**Original Instructions**

Doc. no.: 12607025-01\_A\_en  
Date of release: 2021.02.15

---

**Copyright**

The contents of this manual are the property of Struers ApS. Reproduction of any part of this manual without the written permission of Struers ApS is not allowed.

All rights reserved. © Struers ApS 2021.04.12.

---

# Table of Contents

<b>1</b>	<b>Intended use</b> .....	<b>4</b>
<b>2</b>	<b>Device description</b> .....	<b>4</b>
<b>3</b>	<b>Specimen holders safety precautions</b> .....	<b>4</b>
3.1	Read carefully before use .....	4
<b>4</b>	<b>Overview</b> .....	<b>5</b>
4.1	Specimen holders for round or rectangular specimens .....	5
4.2	Specimen holders with clamps .....	5
4.3	Specimen holders without holes/plane .....	6
<b>5</b>	<b>Instructions for use</b> .....	<b>6</b>
5.1	Leveling the specimens .....	6
5.2	Using specimen holders with clamps .....	6
5.3	Using specimen holders without holes .....	7
<b>6</b>	<b>Maintenance</b> .....	<b>7</b>
6.1	Weekly .....	8
6.2	Monthly .....	8
<b>7</b>	<b>Disposal</b> .....	<b>8</b>
<b>8</b>	<b>Manufacturer</b> .....	<b>8</b>

## 1 Intended use

The device is intended to be used in combination with: Struers equipment for grinding/polishing:

- AbraPol: 160 and 200 mm diameter specimen holders
- AbraPlan: 160 and 200 mm diameter specimen holders
- Tegramin-30: 160 mm diameter specimen holders
- Tegramin-25: 140 mm diameter specimen holders
- TegraForce-100: 140 mm diameter specimen holders

**Do not use the device for the following**

Any other equipment than stated in the section Intended use.

## 2 Device description

The device is a holder which secures the specimen during the preparation process.

The operator secures the specimen in the device and places it in the equipment.

The operator takes the specimens out of the device when the process is finalized.

## 3 Specimen holders safety precautions

### 3.1 Read carefully before use


1. Read this User Guide, as well as the Instruction Manual for the machine with which you are going to use the specimen holder.

## 4 Overview

### 4.1 Specimen holders for round or rectangular specimens

#### Teardrop type

Use this type of specimen holder if you are working with irregular, cold and hot mounted or unmounted specimens.

#### Rectangular type


These specimen holders are designed for round or rectangular specimens, as long as the specimen has parallel sides in the longitudinal plane and a height of minimum 15 mm.

#### Dual-clamping type

Some specimen holders make it possible for you to fix two specimens at a time by tightening a single screw.

Dual-clamping allows a certain tolerance for the diameter of the specimen - thus making these specimen holders well-suited for all specimens.


### 4.2 Specimen holders with clamps

Use this type of specimen holder if you need to clamp wire ends or plates.


This type of specimen holder has three individual detachable clamps that make it easier to fix several specimens.

You can also prepare both faces and even cut the specimens before preparation.

### 4.3 Specimen holders without holes/plane

Use this type of specimen holder if you are working with large specimens. The specimens must be plane parallel, and you must tape the specimens to the specimen holder using double adhesive tape.

You can also use this type of specimen holder if you need tailor-made solutions.

## 5 Instructions for use

Use Uniforce or a leveling plate for accurate leveling of specimens in a specimen holder. For more information on how to use Uniforce, see the Instruction Manual for this unit.

Uniforce is a device for fast and accurate leveling of specimens in a specimen holder.

### 5.1 Leveling the specimens

1. Place the specimen holder on the Uniforce, or on a leveling disc.
2. Place at least three specimens symmetrically around the center of the specimen holder.
3. Clamp the specimens by tightening the screws.
4. Always use screws with a length that can fit inside the holder.
5. If you are working with small specimens, use inserts to reduce the hole size.
6. After clamping, make sure that the specimens are firmly fixated.

#### Placing the specimen holder in the machine

See the Instruction Manual for the machine you are using.

**Note**

Make sure that the screws do not stick out of the specimen holder, as this will damage the machine.

### 5.2 Using specimen holders with clamps

1. Remove the clamps from the specimen holder.
2. Insert the specimens in the clamps.

3. Fasten the specimens using the screws. It can be necessary to insert a piece of cardboard between the specimens and the clamps and tighten the screws again.
4. Fasten the clamp in a cut-off machine, and do the following:
  - Cut the specimens in the clamp to the same length, approximately 5 mm from the edge.
  - Place the three clamps in the specimen holder and fasten them in in Uniforce or any other leveling device.All three clamps are prepared simultaneously.

#### Placing the specimen holder in the machine

See the Instruction Manual for the machine you are using.

**Note**

Make sure that the screws do not stick out of the specimen holder, as this will damage the machine.

### 5.3 Using specimen holders without holes

1. Make sure that the specimen holder surface is:
  - Clean and dry
  - Plane
  - Degreased
2. Fasten the specimens to the surface of the specimen holder using Struers double adhesive tape. Make sure that the specimens are:
  - Placed symmetrically
  - Of the same height (max specimen height: 14 mm)

#### Placing the specimen holder in the machine

See the Instruction Manual for the machine you are using.

**Note**

Make sure that the screws do not stick out of the specimen holder, as this will damage the machine.

## 6 Maintenance

Proper maintenance is required to achieve the maximum up-time and operating lifetime of the device.

## 6.1 Weekly

1. Clean the specimen holder thoroughly.
2. Check the coupling of the specimen holder. The coupling must be able to move a bit. If this is not the case, contact Struers Service.

## 6.2 Monthly

1. Remove the screws used for clamping the specimens.
2. Clean and lubricate the screws and replace them.

# 7 Disposal

**Note**  
Always follow local regulations.

Specimen holders are made of stainless steel/aluminium and do not contain any component classified as hazardous.

Specimen holders can be disposed of as metal waste.

# 8 Manufacturer

Struers ApS  
Pederstrupvej 84  
DK-2750 Ballerup, Denmark  
Telephone: +45 44 600 800  
Fax: +45 44 600 801  
www.struers.com

### Responsibility of the manufacturer

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations.

The manufacturer assumes no responsibility for errors in the text and/or illustrations in this manual. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the supplied version of the equipment.

The manufacturer is to be considered responsible for effects on safety, reliability, and performance of the equipment only if the equipment is used, serviced, and maintained in accordance with the instructions for use.

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