Struers Cooling Units Instruction Manual



Manual No.: 15767001

Date of Release 2017.06.23



Struers Cooling Units Instruction Manual

Always state Serial No and Voltage/frequency if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the Date and Article No of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

Original instructions. The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2017.

Struers

Pederstrupvej 84 DK-2750 Ballerup Denmark Telephone +45 44 600 800 Fax +45 44 600 801



Struers Cooling Units Safety Precaution Sheet

To be read carefully before use

- 1. The unit must be installed in compliance with local safety regulations.
- To achieve maximum safety and lifetime of the machine, use only original Struers consumables.
- Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
 Do not use flammable cooling fluid.

Take care; the cooling fluid may be hot!

- **4.** Use of gloves and safety goggles is recommended when filling and cleaning the tank.
- 5. The machine must be disconnected from the mains prior to any service.
- **6.** The recirculation pump must be disconnected from the power supply prior to removal from the cooling unit.
- All safety functions must be intact and in working order. Damaged or malfunctioning safety functions must be replaced and/ or repaired
- Keep the area around the tank clean and free from spills. Risk of falling
- **9.** Take care when closing the cover use the handle. Do not insert hands through the openings into the water inlet ducts on the cover.
- **10.** If any malfunction or unusual noise is detected, the apparatus should be stopped immediately and technical service called.

The equipment should only be used for its intended purpose and as detailed in the Instruction Manual.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).



Disposal

Equipment marked with a WEEE symbol a contain electrical and electronic components and must not be disposed of as general waste.

Please contact your local authorities for information on the correct method of disposal in accordance with national legislation.

Table of Contents	Page
1. Getting Started	
Struers Cooling Units Components	2
Getting Acquainted with Struers Cooling Unit	
Assembling the Struers Cooling Unit	
Mounting the Cooling Unit's Components	9
Cooli-1 Control Unit: Rear Panel	10
Connecting Components to the Cooli Control Unit	
Preparing the Cooling Unit for Use	
Filling the Tank	11
Connecting to the Power Supply	
Single-phase Supply2-phase Supply	
Filter Bag	
Filter Tube	
Connecting the Cooling Unit to the Machine	
2. Operation	
Cooli-1 Control Panel	16
Control Panel Indicators/ Functions	
Control Farier maleators/ Farietions	10
3. System Components and Accessories	17
4. Troubleshooting	
Problem Solving	10
Location of Fuses	
Location of Fuses	10
5. Maintenance	
Daily	
Monthly	
Replacing the Cooling Fluid	
Emptying the Cooling Unit Tank	
Refilling the Cooling Unit Tank	23
6. Technical Data	24

1. Getting Started

Struers Cooling Units Components

The cooling unit will be delivered as a number of separate components, which will need to be assembled before use.

Main components


- 1 Wheel-mounted cooling unit tank
- 1 Recirculation pump
- 1 Filter Bag (with adaptor plate for 100 I and 150 I tanks)
- 1 Control unit
- 1 24 V / CAN control cable
- 2 Mains cables
- 1 Cable connection box

Optional Components (to be mounted on the tank)

Magnetic filter Water level indicator


Getting Acquainted with Struers Cooling Unit

Take a moment to familiarize yourself with the location and names of some of the cooling unit's components.



- ① Wheel mounted tank
- ② Control unit
- 3 Recirculation pump
- 4 Filter Bag
- S Magnetic filter


Struers Cooling Units Instruction Manual



- ① Wheel mounted tank

- Wheel modified tall
 Magnetic filter
 Gekka connection
 Recirculation pump
 Control unit
 Filter Tube tray

Struers Cooling Units Instruction Manual




- ② Control unit
- 3 Gekka connection
- Recirculation pump Magnetic filter
- Water level gauge
- Filter tray
- ® Drain valve

Assembling the Struers Cooling Unit


50 I Tank

■ Mount the components on the tank.




Mounting positions for:

- Water level indicator
- ② Magnetic filter
- 3 Not used
- 4 Filter unit
- S Cable holders
- 6 Adapter ring for small pump7 Recirculation pump
- ® Cooli control unit



- Mounting positions for:① Adapter ring for small pump② Recirculation pump
- 3 Not used
- 4 Not used
- S Water level indicator
- 6 Magnetic filter
- ⑦ Cable holders
- ® Cooli control unit
- 9 Filter tray
- XL Filter unit

150 I Tank



Mounting positions for:

① Magnetic filter

- ② Cooli control unit
- ③ Recirculation pump④ XL Filter unit
- ⑤ Filter tray

Mounting the Cooling Unit's	
Components	

Note:

A number of these components are optional.

Control Unit

Mount the Cooli control unit by slotting the unit's bolts into the holes provided and then securing (from underneath) with the supplied nuts.

Note:

The control unit's control panel must be mounted facing the cooling unit tank's handle

Recirculation Pump

Mount the recirculation pump in the hole provided.

Adapter ring:

Small pump: use the adapter ring to reduce the hole size. Large pump: remove the adapter ring to increase the hole size.

Magnetic Filter

Mount the Magnetic filter in the hole provided.

Water level indicator

Mount the Water level indicator in the hole provided.

Filter Bag 50 I tank

■ Mount the Filter Bag in the filter unit.

100 I tank


Place the adapter plate in the filter unit then mount the Filter Bag.

150 I tank

■ Insert the positioning rods in the Filter Bag and mount in the filter unit. Place the plate on top of the bag.

Cooli-1 Control Unit: Rear Panel

The diagram below shows the sockets to which the cables have to be connected.



Note:

The 24 V / CAN (A) is connected to the machine with which the cooling unit is used. The 24 V cable used is supplied with the Cooli control unit.

Connecting Components to the Cooli Control Unit

After mounting the components on the tank, connect the cables to the Cooli control unit:

■ Locate the correct socket and plug in the cable from the component.

The pin and hole locations of the sockets and cables are all different from each other, so that it's only possible to connect the right cable to the right socket.

Preparing the Cooling Unit for Use

Filling the Tank

Insert a clean plastic liner in the tank.

Important!

Ensure that the liner lies flat on the base of the tank so that it will not block the pump.

Important!

The cooling unit will be very heavy when full.

Before filling the tank, check that the cooling unit is in position.

If this is not possible, ensure that the tank is directly in front of its final position (e.g. the cut-off machine's cooling unit compartment), with the tank's wheels in-line with the sides of the compartment and the control unit at the front so that it is ready to be pushed into position without being moved to the left or right.

Note

To avoid corrosion, Struers recommends the use of Struers Additive, Corrozip in the cooling fluid (percentage stated on the Additive container). Remember to top up with Struers Additive each time you refill with water. For machines which mainly cut copper and copper alloys, use Corrozip-Cu.

Fill with cooling fluid comprised of:

50 I tank: 42.7 I water and 1.3 I Corrozip

(A total volume of 44 I is recommended)

100 I tank: 97 I water and 3.0 I Corrozip150 I tank: 145.5 I water and 4.5 I Corrozip

Note

Do not overfill the tank.

Avoid spills whilst moving the cooling unit.

Connecting to the Power Supply

Always remember to switch the power off when installing electrical equipment.

Important

Check that the mains voltage corresponds to the voltage stated on the type plate on the machine

The Cooling Unit is shipped with 2 types of Mains cables:

Single-phase Supply



The 2-pin (European Schuko) plug is for use on single-phase connections.

If the plug supplied on this cable is not approved in your country, then the plug must be replaced with an approved plug. The leads must be connected as follows:

Yellow/green: earth Brown: line (live) Blue: neutral

2-phase Supply



The 3-pin (North American NEMA 6-15P) plug is for use on 2-phase power connections.

If the plug supplied on this cable is not approved in your country, then the plug must be replaced with an approved plug. The leads must be connected as follows:

Green: earth
Black: line (live)
White: line (live)



Connect the mains cable to the power cable on the control unit.

Important!


Place the cable box around the connection of the two cables. This protects the cables from water (Protection Rating IP44).


WARNING!


The output voltage from this cable is 200 – 240V and not 110V. DO NOT use this cable to connect equipment that uses an 110V power supply. Failure to adhere to this may result in material damage.

Hang the cables in the cable hooks mounted on the cooling unit tank (this is to prevent the cable trailing on the ground and being damaged).

Filter Bag








■ Insert the water outlet hose from the connected machine into the hole with the correct size.

Filter Tube



- Insert 5-10 cm of the filter tube through the quick-release clamp.
 - Tip: the end of the rigid connector pipe can be removed whilst mounting the filter tube. Lubricate the sealing ring with grease or soap to facilitate re-insertion.
- Mount onto the outlet hose/pipe, approx. 5 cm from the end of the pipe.
- Place the filter tube in a U-shape in the Cooli filter unit.
- Check that there are no folds in the tube.









During the first cutting operation:

■ Check that the filter tube expands along its full length as it fills with water.

Important!

If the tube is twisted or has folds, **stop cutting** and reposition the tube. Do NOT use the Filter Tube when dry-cutting.

Important!

Do NOT re-use the Filter Tube.

Struers Cooling Units Instruction Manual


Connecting the Cooling Unit to the Machine

- Connect the cooling unit's water outlet hose to the quick coupling on the machine.
- Connect the 24 V / CAN control cable to the Cooli control unit by plugging one end into the machine's control socket (see machine manual for details) and the other end into the socket on the rear panel of the Cooli control unit.
- Place the cooling unit underneath the outlet on the connected machine, or lead the outlet hose from the machine into the filter unit on the cooling unit tank.

2. Operation

Cooli-1 Control Panel

The diagram below illustrates the Cooli-1 control panel.



Control Panel Indicators/ Functions

The pump LED's will only light up after the cooling unit has been started.

Name	Key	Function
1		Constant Green: Power on. Normal operation.
POWER		Constant Red: Cooli is on Standby.
(2)		Constant Green: Normal operation.
PUMP	Constant Red: Error. Pump failure (Overload).	

3. System Components and Accessories

Please refer to the *Struers Cooling System brochure* for details of the range available.

The use of Struers consumables is recommended. Other products (e.g. coolants) may contain aggressive solvents, which dissolve e.g. rubber seals. The warranty may not cover damaged machine parts (e.g. seals and tubes), where the damage can be directly related to the use of non-Struers consumables.


4. Troubleshooting

Problem Solving

The following table describes the problems that users may potentially experience when using the cooling units and contains advice on how to solve them.

Error	Explanation	Action
Water leaking.	Leak in recirculation water hose.	Check the hose and tighten the hose clamp.
	Water overflow in the recirculation water tank.	Remove the excess water in the tank.
Cool unit stops and cannot be restarted.	Blown fuse(s)	Replace using appropriate 15 A, slow-blow fuse(s). See illustration, below.
Samples, cooling unit or equipment corroded.	Insufficient additive for cooling fluid.	Add Struers Additive for cooling fluid to the cooling fluid, using the correct concentration. Check with a refractometer. Follow the instructions in the Maintenance section.

Location of Fuses



5. Maintenance

To ensure a longer lifetime for your equipment Struers strongly recommends regular cleaning.

WARNING!

Accumulated dirt and swarf (cutting debris) can restrict or cause damage to the water inlet valves.

Important

Disconnect the power before performing any maintenance work.

NB

Use of gloves and safety goggles is recommended when filling and cleaning the tank.

Daily

Checking the Cooling Fluid

Check and if necessary replace the cooling fluid. (see the section on Replacing the Cooling Fluid for instructions) Top up with cooling fluid if required; the cooling fluid level should be 25 mm below the upper edge of the tank.

IMPORTANT!

Replace the cooling fluid immediately if infected by algae or bacteria.

Remember to add Struers Additive Corrozip: One part of Additive for 33 parts of water.

To check the concentration of additive, use a refractometer.

Concentration = $1.9 \times Brix value$.

The concentration of additive should always be between 2.7 and 3.3 %. Add Corrozip if the concentration is too low.

Filter Bag

- Check and if necessary empty and clean the Filter Bag.
 - Disconnect the power.
 - Remove the cooling unit from the machine's compartment.
 - Remove the Filter Bag from the tank and turn it upside down over a waste container of suitable size.
 - Let the debris fall down into the waste container and clean the Filter Bag.
 - Replace the Filter Bag on the tank.

Important!

Swarf must be disposed of according to the current safety regulations for handling and disposal of swarf/ additive for cooling fluid.

Filter Tube

Check and if necessary replace the Filter Tube.
 (Change the Filter Tube when it is full of swarf)

Important!

Used Filter Tubes must be disposed of according to the current safety regulations for handling and disposal of swarf/ additive for cooling fluid.

Please note:

Depending on which metals are being cut, it is possible that the combination of the metallic swarf (cutting debris) from metals with a large difference in electro positivity (a large distance apart in the electrochemical series), could result in exothermic reactions when 'favourable' conditions are present.

Therefore it is always good practice to bear in mind which metals are being cut and the amount of swarf produced.

Examples:

The following are examples of combinations which could result in exothermic reactions if a large amount of swarf is produced during cutting/ grinding on the same machine, and when favourable conditions are present:


Aluminum and Copper Zinc and Copper

Magnetic filter

- Check and if necessary clean the Magnetic Filter.
 - Remove the filter from its mounting.
 - Slide the plastic outer tube off the magnet.
 - Use a stiff brush to clean the plastic tube.
 - Refit the plastic tube.
 - Place the filter back in its mounting.

Static Filter

- Check and if necessary clean the Static filter.
 - Remove the pump.
 - Loosen the wing nut and remove the Static filter.
 - Use a stiff brush to clean the mesh.
 - Rinse the static filter and re-mount.



Monthly

 Clean painted surfaces, and the control unit with a soft damp cloth and common household detergents.
 For heavy duty cleaning, use Struers Cleaner (Cat. No. 49900027).



Note:

Ensure that no detergent or cleaning agent residue is flushed into the cooling unit tank; excess foaming will occur.

Replacing the Cooling Fluid

Replace the cooling fluid in the Recirculation Cooling Unit at least once a month.

REMEMBER...

Avoid skin contact with the additive for cooling fluid.

Emptying the Cooling Unit Tank

IMPORTANT!

The cooling fluid will contain additive and cutting swarf and may **NOT** be disposed of into a main drain.

Cooling fluid must be disposed of in compliance with local safety regulations.

- Empty the cooling unit using the recirculation pump:
 - Disconnect the cooling unit's water outlet hose from the quick coupling on the equipment.
 - Place the open end of the hose in a suitable container.
 - Start the equipment.
 This will start the recirculation pump, to pump the cooling fluid out of the tank.
 - Stop the equipment when the tank is empty.
 - Empty the cooling unit completely, remove the plastic liner and clean out all water and debris from the tank.
 - Clean the tank and connected tubes thoroughly.
 If the cooling fluid has been infected with bacteria or algae, flush the tank and tubes with a suitable antibacterial disinfectant e.g Struers Unitclean.

150 I tank

Alternatively, the tank can be emptied using the drain valve (the disposable liner must be moved to allow the cooling fluid to flow freely) or by using an industrial dredge pump.

Refilling the Cooling Unit Tank

IMPORTANT!

The cooling unit will be very heavy when full.

Before filling the tank, check that the cooling unit is in position.

If this is not possible, ensure that the tank is directly in front of its final position (e.g. the cut-off machine's cooling unit compartment), with the tank's wheels in-line with the sides of the compartment and the control unit at the front so that it is ready to be pushed into position without being moved to the left or right.

- Refill the cooling unit with water using an external water hose or via direct connection to the water mains. Stop refilling when the water level is 25 mm from the upper edge.
- Fill with cooling fluid comprised of :

50 I tank: 42.7 I water and 1.3 I Corrozip

(A total volume of 44 l is recommended)

100 I tank: 97 I water and 3.0 I Corrozip 150 I tank: 145.5 I water and 4.5 I Corrozip

NB!

Never refill by pouring the water into the cutting chamber/ turntable compartment, as it would be impossible to check the increasing level in the tank.

IMPORTANT!

Always maintain the correct concentration of Struers Additive, Corrozip, in the cooling fluid (percentage stated on the container of the Additive). Remember to add Corrozip each time you refill with water. For machines which mainly cut copper and copper alloys, use Corrozip-Cu.

Use only Struers Additive in the Recirculation Cooling Unit.
Struers Additive is specially selected for use with Struers equipment.
Other additives may not be compatible with some of the components of the equipment.

6. Technical Data

Subject	Specification		
	Metric/International	US	
50 I tank		·	
Height	26 cm	10.2 "	
Width	52 cm	20.7 "	
(with extended handle)	(53 cm)	(21.0 ")	
Depth	46 cm	18.1 "	
Volume	50	13.2 gallons	
100 I tank			
Height	26 cm	10.2 "	
Width	73 cm	28.7 "	
(with extended handle)	(74 cm)	(29.0 ")	
Depth	67 cm	26.4 "	
Volume	100 I	26.4 gallons	
150 I tank			
Height	74 cm	29.1 "	
Width	83 cm	32.7 "	
Depth	76 cm	29.9 "	
Volume	150 l	39.6 gallons	
Small pump			
Flow	60 l/min at 1 bar	16 gallons/min at 1 bar	
Power		0.16 kW	
Large pump			
Flow	125 l/min at 1 bar	33 gallons/min at 1 bar	
Power		0.7 kW	
Operating environment			
Surrounding temperature	5-40°C/41-104°F		
Humidity	0-95% RH non conden	sing	



English

Declaration of Conformity

Manufacturer Struers ApS

Pederstrupvej 84

DK-2750 Ballerup, Denmark Telephone +45 44 600 800

Herewith declares that

Name: Struers Cooling Units

Cat. No.: 057665xx, 057666xx, 057667xx, 057668xx,

057670xx, 057671xx, 057672xx

Function: Recirculating unit for cooling water for a cut-off

machine or a similar type of equipment

Type No.: 576

fulfils all the relevant provisions of the:

Machinery Directive according to the following standard(s): 2006/42/EC EN ISO 12100:2010, EN ISO 13849-1:

EN ISO 12100:2010, EN ISO 13849-1:2015,

EN ISO 13850:2008, EN 60204-1:2006/AC:2010.

and is in conformity with the:

EMC Directive according to the following standard(s):

2014/30/EU EN 61000-3-2:2014, EN 61000-3-3:2013, EN 61326-1:2013.

RoHS Directive according to the following standard(s):

2011/65/EU EN 50581:2012.

Supplementary Information The equipment complies with the following standards:

NFPA70:2014, NFPA79:2012, FCC 47 CFR Part 15, AS/NZS 4251.1

The above has been declared according to the global approach, module A.

Authorized to compile the Technical File:

Klavs Tvenge

Director of Business Development

Struers ApS

Pederstrupvej 84

DK-2750 Ballerup, Denmark Date of Issue: 2017.10.05

Doc. No.: 15767901

Rev.: A



Pederstrupvej 84 DK-2750 Ballerup Denmark

Struers Cooling Units



Manual No.: 15767001

Date of Release @17.06.23

Spare Parts and Diagrams





Pederstrupvej 84 DK-2750 Ballerup Denmark

Struers Cooling Units Spare Parts and Diagrams

Always state Serial No and Voltage/frequency if you have technical questions or when ordering spare parts.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manuals may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manuals may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2017.

Struers

Pederstrupvej 84 DK-2750 Ballerup Denmark Telephone +45 44 600 800 Fax +45 44 600 801

Spare Parts and Diagrams

Table of contents	Drawing
Reservoir 100L, complete	15760010G
Reservoir 50L, complete	
Tank 150 I, assembly	
Magnetic filter, complete	
Filter console welded	
Filter Bag Adapt. w. XL Static Filter	
24V Control box, complete	
Door for small Control box, assembly	
Big Pump	
1x100V/50Hz	15760090E
1x120V/60Hz	
1x220V/50Hz	15760092F
1x240V/60Hz	
Small Pump	
1x100V/50Hz	15760094C
1x120V/60Hz	
1x220V/50Hz	15760096C
1x240V/60Hz	
1x100-120V/50-60Hz	15760098C

Some of the drawings may contain position numbers not used in connection with this manual.


Diagrams	Diagram Nr.
Wiring diagrams	
W1 cable	15763501C
Cooli-1	15763593E
Contactor	15763581B
Other Diagrams	
Cooli-1 PCB, 2 pages	15763004


Spare Part list for Struers Cooling Units


Drawing	Pos.		Cat no.
15760010		Reservoir 100L, complete	
	10	Body of reservoir-100L-welded	05766905
	20	Cover of reservoir-100L-welded	15760555
	30	Wheel 2478PJH075P60, 2 pcs	2GB00110
	100	Wheel 2470PJH075P40, 2 pcs	2GB00111
15760011		Reservoir 50L, complete	
	10	Body of reservoir-50L-welded	05766906
	30	Cover of reservoir-50L-welded	15760702
	100	Wheel 22198UOI045P60-67X24, 2 pcs	2GB00112
	110	Wheel 1430UOI038P60-55x39, 2 pcs	2GB00113
15760086		Tank 150 I, assembly	
	10	Tank 150 I, welded	15041081
	30	Cover of reservoir-150L, welded	15041084
	20	Wheel Fixed Ø65mm, 4 pcs	2GB00130
15760014		Magnetic filter, complete	
	10	Magnetic separator, glued	15760610
	60	O-ring 15.60-2.40 72 NBR 872	21024034
15760660		Filter console welded	
		Filter bag	05766928
		Adapter plate for filter bag (05766928) on 100 I tank	05766927
15760089		Filter Bag Adapt. w. XL Static Filter	
	10	Filter Box, welded	15760180
	20	Filter Bag Rod	15760197
	30	Top Closure	15760114
	40	Filter Bag Assembly	15760006
15760028		24V Control box, complete	
	20	Door for small control box,asm.	15760029
	80	CAP for chassis connector IP68	2XM90000
	90	Cap, fem. sock. 08-0202-17-00	2XM90202
	110	Sealing strip Black 3x15 4030	2IP00302
	150	15AT FUSE MDA 6.3 x 32 250V, 2 pcs	2FU17550
	160	Fuse holder bayonet No.800	2XS19800
	170	Fuse cap bayonet No.816	2XS19816


Spare Part list for Struers Cooling Units


Drawing	Pos.		Cat no.
15760029	20	Door for small Control box, assembly PCB Cooli-1, Tested	15763003
15760092/ 15760091/ 15760092/ 15760093	30	Big Pump GEKA Coupling ¾" int. thread	2NF60002
15760094/ 15760095/ 15760096/ 15760097/ 15760098	13	Small Pump GEKA Coupling ¾" int. thread	2NF60002
		Loose part for Pumps	
	10	Hose Clamp ø12-20 / 9.0-C6 W2, 2 pcs	2NS21220
	20	Hose clamp Norma S24/9Zy, 2 pcs	2NS12409
	30	Reinforced Pressure Hose 5/8 ø16x3, 2.5 m	2NU29322
	40	Reinforced Hose, PVC ½"-Ø12.5, 2.5 m	2NU29316
	50	GEKA Pipe branch ½", 2 pcs	2NF60000
	60	GEKA Pipe branch ¾", 2 pcs	2NF60001
		Static Filter Insert	05766934
		Struers Cooling Units, "loose parts"	
		Mains Cable, 0.75mm², Schuko	2WC04668
		W16 Cooli Connection Cable	15763516
		Mains Cable, AWG16, Nema 5-15P	2WC02520





























The Pump is delivered with Terminal block from the supplier. The Terminal block in the Wire set shall not be used. (Scrapped)


\neg									
	Е	2016-03-30	Text about Terminal block added.				BRY	2016-03-30	JJO
	Α								
	Rev	Crea. date		Revision description		Draw. Init	Appr. date	Appr. Init	
	1/2 OL		Material:			Scale: Format:		Tolerance: DS/I	SO 2768 - mK
Struers				1:5		A4	Weight:	g	
	Pederstrupvej 84 DK-2750 Ballerup		ID: Description:					Rev:	
	Copen Denma Phone	hagen	15760090 Big pump 1x100 50 Hz CSA, complete						E











Crea_date pyly-min-dd | Revision description | Draw. Init | Appr. date pyly-min-dd pyly-mi

























ASSEMBLING of the Simple PCB (12050426-3) detail







Pederstrupvej 84 DK-2750 Ballerup Denmark