# **Magnutom** Instruction Manual



Manual No.: 16147001 Revision A

Date of Release 2018.01.16



#### Important

READ the instruction manual carefully before use. Keep a copy of the manual in an easy-to-access place for future reference.

Intended use:

For automatic cutting of metallic or other solid materials for further material inspection and only to be operated by skilled/trained personnel. The machine is only designed to be used with abrasive cut-off wheels specially designed for this purpose and this type of machine.

#### Do NOT use the machine for:

Cutting of materials other than solid materials suitable for materialographic studies. In particular, the machine must not be used for cutting any type of explosive and/or flammable material, or materials which are not stable during machining, heating or pressure. The machine may not be used with cutting wheels which are not compatible with the machine requirements (e.g. saw-blade or toothed cutting wheels).

The machine is for use in a professional working environment (e.g. a materialography laboratory).

#### **Table of Contents**

#### Page

User's Guide	9
Reference Guide	57
Quick Reference	102
Declaration of Conformity	103

obligations: Instruction Manuals: Struers Instruction Manual may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manual may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to changes without notice. The manual may mention accessories or parts not included in the present version of the equipment.

Original instructions. The contents of this manual is the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2018.

Struers Pederstrupvej 84 DK-2750 Ballerup Denmark Telephone +45 44 600 800 Fax +45 44 600 801

Always state Serial No and Voltage/frequency if you have technical questions or when ordering spare parts. You will find the Serial No. and Voltage on the type plate of the machine itself. We may also need the Date and Article No of the manual. This information is found on the front cover.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal



## Magnutom Safety Precaution Sheet

### To be read carefully before use

- 1. The machine must be installed in compliance with local safety regulations.
- 2. The machine must be placed on a stable floor which is able to carry at least 3000 kgs / 6600 lbs (Magnutom, Coolimat and specimens). The machine must be levelled by means of the adjustable legs provided.
- **3.** When lifting the machine by the built-in truck lifting point, ensure that the boom is properly secured with the locking pins provided.
- **4.** Before transport, secure the drive shaft with the locking system provided.
- 5. The operator(s) must read the Safety and User's Guide sections of this manual and the relevant sections of the manuals for any connected equipment and accessories. The operator(s) must read the Instructions for Use and, where applicable, Safety Data Sheets for the applied consumables.
- Use only intact cut-off wheels. The cut-off wheels must be approved for a spindle speed of min. 1950 rpm / 42 m/s (Magnutom-400), 1850 rpm / 50 m/s (Magnutom-500).
- 7. The machine is not for use with saw-blade type cut-off wheels.
- B. Do not use the machine for cutting materials that are flammable or unstable during the cutting process (e.g. combustible or explosive materials).
   Do not use the machine for cutting materials that are not suitable for materialographic cutting.
- **9.** Observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
- **10.** The workpiece must be securely fixed in the a quick-clamping device or similar. Large or sharp workpieces must be handled in a safe way.
- **11.** Use of working gloves is recommended as workpieces may be both very hot and produce sharp edges. Wearing of gloves is also recommended when flushing and cleaning the machine.
- **12.** Use of safety shoes is recommended when handling large or heavy workpieces.
- **13.** Use of safety goggles and gloves are recommended when using the flushing hose.

- **14.** Do not work on or around cutting table when the table is repositioned using the positioning joystick
- **15.** For maximum safety and lifetime of the machine, use only original Struers consumables.
- 16. Never look directly into the laser beam. (Laser straight edge option).
- **17.** Struers recommend the use of an exhaust system as the cutting materials may emit harmful gasses or dust
- **18.** The machine emits only moderate noise. However, the cutting process itself may emit noise, depending on the nature of the workpiece Use hearing protection if exposure to noise exceeds levels set by local regulations.
- 19. When a recirculation cooling unit is used, observe the current safety regulations for handling, mixing, filling, emptying and disposal of the additive for cooling fluid.
  Do not use flammable cooling fluid.
  Use of gloves and goggles are recommended.
  Do not use cooling liquid other than water and Struers additives for cooling fluid.
- 20. The machine must be disconnected from the mains prior to any service.
- **21.** All safety functions of the machine must be intact and in working order. The machine must be installed in compliance with local safety, regulations.
- **22.** Ensure that the cut-off wheel is secured before working on or around the cutting table.
- 23. Only use the flushing gun for cleaning *inside* the cutting chamber.
- 24. If any unusual noise is heard during operation, refrain from further use of the machine, and contact a Struers Service Technician.

The equipment is designed for use with consumables supplied by Struers. If subjected to misuse, improper installation, alteration, neglect, accident or improper repair, Struers will accept no responsibility for damage(s) to the user or the equipment.

Dismantling of any part of the equipment, during maintenance, service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.).

The equipment should only be used for its intended use and as detailed in the Instruction Manual.

# **Icons and typography**

The following icons and typographic conventions are used in this instruction manual:

Icons and Safety Messages



### **ELECTRICAL HAZARD**

indicates an electrical hazard which, if not avoided, will result in death or serious injury.



### WARNING

indicates a hazard with a medium level of risk which, if not avoided, could result in death or serious injury.



### CAUTION

indicates a hazard with a low level of risk which, if not avoided, could result in minor or moderate injury.

**General Messages** 



#### Important or Note indicates a risk of damage to property, or the need to proceed with special care.



*Information or Tip* indicates additional information and tips.

#### Colour Inside Logo



The 'colour inside' logo on the cover page of this Instruction Manual indicates that it contains colours which are considered to be useful for the correct understanding of its contents. Users should therefore print this document using a colour printer.

Typographic conventions

Bold typeindicates button labels or menu options in<br/>software programsItalic typeindicates product names, items in software<br/>programs or figure titlesBulletsindicates a necessary work step



### Disposal

Equipment marked with a WEEE symbol  $\stackrel{\boxtimes}{=}$  contains electrical and electronic components and must not be disposed of as general waste.

Please contact your local authorities for information on the correct method of disposal in accordance with national legislation.

# **User's Guide**

### **Table of Contents**

Page

# **1.Getting Started**

-	
Checking the Contents of Packing	11
Unpacking and Placing Magnutom	11
Getting Acquainted with Magnutom	12
Hold-to-run operation	13
Power Supply	16
Compressed Air Connections	17
Emptying the Water / Oil Filter	17
Mounting the Cut-off Wheel	18
Connection to an External Exhaust System	19
Noise level	19
Connecting a Recirculation Cooling Unit	20
Connecting a Coolimat-200 Recirculation Cooling Unit	20
Connecting other Struers Cooling Systems	20
Connecting Other External Filter Units	21
Moveable X-Table (option)	21

# 2. Operation

Using the Controls	22
Control Panel of Magnutom	22
Sleep Mode	23
Using the Joystick	23
Signals	24
Cutting chamber lights	24
Signalling light (option)	24
Navigating the Software	26
Setting the Language	26
Main Menu	27
Changing the Language	28
Editing Numeric Values	30
Editing Alphanumeric Values	31
Editing Cutting Method Names	33
Changing Cutting Mode and Cutting Parameters	35
Feed speed	37
Cut length	37
Cut-off wheel	38
Wheel Speed	38
Cutting Mode	39
Sampling Method	39
Changing the Cut-off Wheel	40
Clamping the Workpiece	41
Cutting on Magnutom	42
Starting the Cutting	42

The Cutting process screen	43
Manual Stop	44
Re-starting cutting	44
Additional cooling	44

# **3. Maintenance**

General Cleaning	45
Daily	45
AxioWash	46
Cleaning the Cutting Chamber	47
Maintenance of Cut-off Wheels	49
Storage of Bakelite Bonded Al <sub>2</sub> O <sub>3</sub> Cut-off Wheels	49
Maintenance of Diamond and CBN Cut-off Wheels	49
Weekly	50
Cleaning the Cutting Chamber	50
Monthly	50
Replacing the Cooling Fluid	50
Lubricating the Movement Mechanisms	50
Maintenance of Cutting Table	52
Lubricating the Guide Shafts	52
Lubricating the Drive Shaft Spindle	53
Adjusting the Laser	54
Yearly	55
Inspection of Cover	55
Testing Safety Devices	56

# 1. Getting Started

#### Checking the Contents of Packing

# . .. .\_. .

Unpacking and Placing Magnutom

- The packing box contains the following items:
- 1 Magnutom
  - (high capacity cut-off machine)
- 4 Drain pipes for water outlet
- 3 Elbow pipes for water outlet
- 1 Fork spanner 36 mm
- 1 Key for hardware compartment
- 1 Key for door lock
- 1 Coolimat connection kit
- 1 Instruction Manual Set
- Move the crate as close as possible to the desired location, using a fork lift.
- Remove the sides of the crate.
- Remove the bolts securing Magnutom to the pallet.
- Using the forklift, lift Magnutom approx. 5 cm
- Remove the pallet and the plastic covering
- Lower Magnutom slowly to the ground.
- Using the special rollers, move the machine to its final position.
- Adjust the 6 legs to level the Magnutom.
- Allow Magnutom to reach room temperature before connecting electrical power.
- Remove the lifting slots from under Magnutom and store in a safe place.



Getting Acquainted with Magnutom

Take a moment to familiarise yourself with the location and names of the Magnutom components.

- ① Signalling light (optional)
- ② Control Panel
- ③ Emergency stop
- ④ Water outlet
- ⑤ Demister (optional)
- © Laser straight edge (optional)
- ⑦ Safety lock release
- ⑧ Side door
- In the second second
- ① Access door
- 1 Hold-to-run button



MAIN SWITCH

The main switch is located on the right hand side of the machine.

Turn clockwise to switch on the power.



The EMERGENCY STOP is located on the front of the machine. Emergency Stop

- Push the red button to Activate.
- Turn the red button clockwise to Release.



#### Important

Do not use the Emergency stop for operational stop of the machine during normal operation. BEFORE releasing (disengaging) the Emergency stop, investigate the reason for activating the Emergency stop and take any necessary corrective action.

Hold-to-run operation

When moving the cut-off wheel whilst the protection hood or side door are open, press the Hold-to-run button on the front of Magnutom before operating the joystick.



### Inside the Cutting Chamber



- ② Flexible water jets
- ③ Guard for cut-off wheel
  ④ Removable grate (for catching specimens/ large pieces of debris)
- S Cutting table
- © Support block (only for machine transport then must be removed)
- ⑦ Flushing gun



- ④ Water outlet
- ⑤ 24 V connection for Recirculation Cooling Unit
- 6 Compressed air inlet
- Quick coupling for water inlet

15

**Power Supply** 



### **ELECTRICAL HAZARD**

- Switch the power off when installing electrical equipment.
- The machine must be earthed.
- Check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine. Incorrect voltage may result in damage to the electrical circuit.
- Open the electric connection box and connect a 4-lead or 5-lead cable<sup>1</sup> in the following way:
  - Yellow/green: Earth (ground) Brown: Line (live) L1, L2,
  - Brown: Line (live) L1, L2, L3 Blue: Neutral (this terminal is c
    - Neutral (this terminal is only used for a mechanical connection if there is a Neutral wire in the cable)
- The other end of the cable can be fitted with an approved plug or hard-wired into the mains, according to the electrical specifications and local regulations.



#### Important

Check that the mains voltage corresponds to the voltage stated on the type plate on the side of the machine.



#### Note!

The doors can only be opened when the machine is connected to a power supply and the main power switch is on. To open the doors when the power is not connected, use the triangle key to release (de-activate) the safety lock.



*Tip* Magnutom cannot operate if the safety lock is de-activated. Remember to re-activate the safety lock before operating Magnutom.

<sup>&</sup>lt;sup>1</sup> Please see the section on *Technical Data* at the rear of the Instruction Manual for recommended cable specifications.

**Compressed Air Connections** 

To connect compressed air:

Connect an 8 mm compressed air hose to the compressed air inlet.

*Important* The air pressure must be between 6 bar (87 psi) and 9.9 bar (143 psi).

Emptying the Water / Oil Filter

Magnutom is fitted with a water / oil filter that removes excessive amounts of these substances from the compressed air supply. As a result of this, it is necessary to empty the filter periodically:

Hold a container under the release outlet valve and press the release valve.



**Release Outlet Valve** 



# Connection to an External Exhaust System

Struers recommends the use of an exhaust system as workpieces may emit harmful gases when cut. When a central exhaust system is not available, the optional Demister is recommended.

Magnutom is prepared for connection to an exhaust system via a 160 mm (approx.  $6\frac{1}{4}$ ") fitting on the top of the cabinet. Recommended minimum capacity for exhaust system: 700m<sup>3</sup>/h / 25,000 ft<sup>3</sup>/h at 0 mm /0" water gauge.

 Mount an exhaust hose from your local exhaust system onto the flange.



**Noise level** 

Approx.  $67^2$  dB (A) measured at idle running, at a distance of 1.0 m/39.4" from the machine Use hearing protection if exposure to noise exceeds levels set by local regulations.

<sup>&</sup>lt;sup>2</sup> "The figures quoted are emission levels and are not necessarily safe working levels. Whilst there is a correlation between the emission and exposure levels, this cannot be used reliably to determine whether or not further precautions are required. Factors that influence the actual level of exposure of the workforce include the characteristics of the workroom and the other sources of noise, i.e. the number of machines and other adjacent processes and the length of time for which an operator is exposed to the noise. Also, the permissible exposure level can vary from country to country. This information, however, will enable the user of the machine to make a better evaluation of the hazard and risk." (ref. European standard EN ISO 16089:2015)

Connecting a Recirculation Cooling Unit To ensure optimal cooling, Magnutom should be fitted with a Recirculation Cooling Unit.

Struers Coolimat-200 is designed for use with large cut-off machines such as Magnutom. Coolimat-200 is available as a Band Filter or a Static Filter unit.

Note



### Before connecting the cooling unit to the Magnutom, follow the instructions in the Cooling Units Instruction Manual to prepare it for use.



### CAUTION

During cutting, the cooling fluid coming out of the water outlet may be very hot.

Connecting a Coolimat-200 Recirculation Cooling Unit

- Lead the long drain pipe through the opening in the wall of the cabinet, and then connect it to the drain outlet under the cutting table.
- Connect to Coolimat-200 using the pipes and connectors supplied.
- Magnutom, connect the other end to the pump on the cooling unit.
- Connect the 24 V / CAN control cable (supplied with Coolimat-200) to the 24 V socket on the right hand side of Magnutom, and the other end to the Cooli control unit.

Connecting other Struers Cooling Systems

- Slide the cooling unit into the compartment under the Magnutom.
   Connect the 24 V / CAN control cable (supplied with the Cooling System) to the 24 V socket on the right hand side of Magnutom, and the other end to the Cooli control unit.
- Connect the water inlet tube to the quick coupling on the side of Magnutom, connect the other end to the pump on the cooling unit.
- Close the compartment door.

Connecting Other External Filter Units



#### Important

Please contact a qualified electrician to verify that the external filter unit can be used with Magnutom. The electrical diagrams in the Spare Parts section of the manual can be used for identification of the different wires. The pressure of the cooling fluid supplied to Magnutom must be max. 2 bar.

- Mount an elbow pipe on the drain outlet under the cutting table.
- Lead the drain pipe through the opening in the wall of the cabinet, and then connect it to the elbow pipe.
- Connect to the external unit using the pipes and connectors supplied.
- Connect the water inlet tube to the quick coupling on the side of Magnutom; connect the other end to the pump on the cooling unit.
- Connect a 24 V electric (available as an accessory, part no 15483549) to the 24 V socket on the right hand side of Magnutom, and the other end to the external unit.
   The electric 24 V signal (max 200 mA) can be used to start a pump or open a magnetic valve for the water supply. It might be necessary to use a relay to control connected equipment.

### Moveable X-Table (option)



The X-table option is a motor driven, movable table. The table can be moved left to right by moving the joystick left or right. The X-table option is required to use the Magnutom's MultiCut functions.

# 2. Operation





Signals Cutting chamber lights	The Cutting chamber lights will flash if Magnutom has been idle for 30 seconds after a cutting process is completed. (The lights will not flash if the operator has interacted with Magnutom)		
Signalling light (option)			
Signalling light	The optiona current conc	l signalli dition of	ng light mounted on top of Magnutom signals the the machine.
	<mark>Green</mark> I	Magnuto	om is busy cutting.
	Yellow /	Attentior	n required
		Cutting p	process is complete and Magnutom is idle
	Red /	An unex immedia	pected stop has occurred and requires te attention
Beeps	The "Signall attention to "	ling light the som	sound"is a configurable signal for directing e of the state changes to "Emergency" and
	Emergency	state:	signaled by a repeated series of a beep and a short pause
	Stopped sta	te:	signaled by a repeated series of a beep and a long pause
	The beeps o	can be s	topped by acknowledging the pop-up message.

# **Overview of Signals\***

State	" <b>Operating</b> " Green light	" <b>Stopped</b> " Yellow light	" <b>Attention</b> " Yellow light	" <b>Emergency</b> " Red light
	Operation is in progress.	Operation is stopped. A caution or marginal condition is detected, and attention is required.	Stopping is impending. A caution or marginal condition is detected, and attention is required.	Operation is stopped. A protective device is engaged, or immediate action is required to deal with a hazardous condition.
Event				
Start or Continue is pressed	State change ⇒ Operating	State change ⇒ Operating		
Process is complete			State change ⇔ <mark>Stopped</mark>	State change ⇔ Stopped
Manual stop is pressed			State change ⇔ <mark>Stopped</mark>	State change ⇔ Stopped
A door is opened			State change ⇔ <mark>Emergency</mark>	State change ⇔ <mark>Emergency</mark>
Emergency stop is pressed			State change ⇔ <mark>Emergency</mark>	State change ⇔ <mark>Emergency</mark>
Thermal overload	Beep ON	State change ⇔ Emergency	State change ⇔ Emergency	State change ⇔ Emergency
		Beep ON	Beep ON	Beep ON

\* when optional Signalling light is fitted.

Navigating the Software	Use the TURN/PUSH knob to navigate the software:		
	Push knob to select a menu or item.		
	Turn knob to move the cursor or to adjust settings.		
	Press knob to store modified settings.		
	ESC Press Esc to return to the previous menu (without saving any modified settings).		
	The <i>Main menu</i> is the highest level in the menu structure. Cutting methods can be selected and edited and/ or started.		
	Other functions such as defining user cut-off wheels or changing operation mode can be carried out in the <i>Maintenance</i> and <i>Configuration</i> menus. Details of these functions can be found in the <b>Reference Guide</b> section of the Instruction Manual.		
Setting the Language	When the Magnutom is switched on for the first time, a screen on the control panel prompts you to select your preferred language. Then to set the date and time.		
	<ul> <li>From the pop-up menu, select your preferred language by turning the knob then enter the setting by pushing the knob.</li> </ul>		
	Press Esc to move to the Main Menu.		

Whenever the Magnutom is switched on the screen displays the version number of the software.

Magnutom 500	Version 0.71 X
<u>SERVICE INFO:</u> Total operating time: Time since last service: Time until next service: Memory module installed	0 h 0 h 220 h I

Additionally some service information is supplied; the total operation time and the time since the last service of the machine.

The display will then change to the last screen shown before Magnutom was switched off.

Main Menu

During normal operation, immediately after start up, where the splash screen is displayed, the software goes to the screen that was used before the machine was switched off. Thus you can continue exactly where you left last time the machine was used.

The Main Menu is the highest level in the menu structure. From this menu, you can enter all the other menus. To enter the *Main menu*:

- Press Esc until the Main menu is displayed.





### **User options**

Display brightness	40
Language	English
Keypad sound	On
Units	Metric
Number of methods	10
Time [hh:mm:ss]	12:59:32
Date [yyyy-mm-dd]	2012-05-31
Default value	

### ţ

Push knob to activate the Select language pop-up menu.

Turn knob to select the language you prefer.

### Ļ

User op	Select language	1
Display Languag Keypad Units Number Time [hh Date [yy	English Deutsch Français Español 日本語 中文 Italiano Русский 한국어 Polski	40 English On Metric 10 :59:32 ·05-31
Default v	·	J



...

Push knob to accept the language.

The Configuration menu now appears in the language you have chosen.

Check if there are any other settings that need changing in the Options menu. If not, Push ESC to return to the Configuration menu.

Otherwise use the Turn/Push knob to select and change the required parameters.

Editing Numeric Values

Turn knob to select the value to be changed, e.g. Display brightness:

### ↓

User options	
Display brightness	40
Language	English
Keypad sound	On
Units	Metric
Number of methods	10
Time [hh:mm:ss]	12:59:32
Date [yyyy-mm-dd]	2012-05-31
Default value	

### ↓

Ť

Push knob to edit the value.

A scroll box appears around the value.

		User options	
		Display brightness Language Keypad sound Units Number of methods Time [hh:mm:ss]	40 English On Metric 10 12:59:32 0010_05_21
		Default value	2012-03-31
	Ļ	Note: If there are only two options, the po displayed. Pressing the knob (Ente the 2 options.	pup box is not r) will toggle between
	$\bigcirc$	Turn knob to increase or decrease to toggle between the two options).	he numeric value (or to
	↓ •	Push knob to accept the new value. the changes, preserving the original	(Pressing Esc, aborts value.)
Editing Alphanumeric Values	$\mathbf{\hat{O}}$	Turn knob to to select the text value Keyboard sound:	e to be changed, e.g.
	Ļ		

### **User options**

Display brightness	40
Language	English
Keypad sound	On
Units	Metric
Number of methods	10
Time [hh:mm:ss]	12:59:32
Date [yyyy-mm-dd]	2012-05-31

### ţ

1

Push knob to toggle between the 2 options.

- 1	
л.	
¥	·

User options	
Display brightness	40
Language	English
Keypad sound	Off
Units	Metric
Number of methods	10
Time [hh:mm:ss]	12:59:32
Date [yyyy-mm-dd]	2012-05-31
Default value	

If there are more than two options, a popup box is displayed. Turn knob to select the correct option.

Note

# ↓

ţ

Press Esc to accept the option and return to the previous menu

Or turn knob to select and edit other options in the menu.

(...) Push knob to open the *Cutting methods* menu. ↓ Cutting methods 1. Method 2.5mm/s 📀 10851 2. Method 0.5mm/s 😣 10851 ⊕ 3. Method 0.5mm/s 🙆 10851 0 •))) Rename Turn knob to select a method. Ť

Up to 10 cutting methods can be saved in the database.

#### Editing Cutting Method Names

F1 Press F1 to open the *Text editor* menu.

Text editor
Current text: Method Edited text: My Method
abcdefghijk Imnopqrstuv wxyzSpace-/#µ, 1234567890.
DEL + - Save & Exit
UPPER CASE



Save the changes in the database by selecting *Save* & *Exit*.

#### 1. Method Α 100mm В С Omm/s н 1 1500rpm D **0S51** Ε G F **.**|....| Set feed speed

Cutting methods can be edited from the editing menu.

# Changing Cutting Mode and Cutting Parameters

- A Cutting method number
- **B** Feed speed
- **C** Cut Length
- **D** Cut-off wheel
- E Cut-off wheel speed
- **F** Cutting mode
- **G** MultiCut setup (only available with x-table option)

Editing the Parameters

The cutting parameters can be set both prior to and during cutting. From the *Cutting methods* menu, select a cutting method:



Push knob to open the editing menu.



Turn knob to highlight the parameter to be edited e.g. feed speed.





Push knob to select the parameter. A pop-up menu appears.



Turn knob to adjust the value of the parameter.





Push knob to save the setting.


Can be set to values between 0.1 - 10 mm/sec, in steps of 0.1 (0.002-0.2 "/s). Default value: 0.5mm/s

There are two ways to set the stop position: Setting a *Cut length* or *AutoStop*.

- Select **Cut-off wheel** and Push the knob.
- The Cutting mode setup menu will appear.

Can be set to values between 1- 385 mm, in steps of 1 mm. Default value: 100 mm.



The available Cut length will differ depending on the cutting mode used. If the value selected is greater than this, a message will appear suggesting the Cut length which is available.

Cut length can also be set when editing the Cutting mode.

Adjust the stop position to '0'mm (Auto). When the *AutoStop* function is selected, the machine automatically stops when the workpiece has been cut through.

For a detailed description of the *Stop Settings* and how they can be used, see the **Reference Guide** section of the Instruction Manual.

Feed speed



Cut length



Setting a stop position



AutoStop

Cut-off wheel



■ Select **Cut-off wheel** and Push the knob.

Select the cut-off wheel from the pop-up menu. When Wheel Selection Mode is set to 'Intelligent', the Material guide menu will appear.

Select the hardness of your Material in the pop-up menu and Magnutom will suggest the suitable Struers cut-off wheel, and automatically retrieve the recommended rotational speed for that specific wheel.



 Select the diameter of the cut-off wheel: Magnutom-500 432 or 508mm, Magnutom-400 432 mm only

Can be set to values between

in steps of 50 rpm. Default value: 1500 rpm. 1000- 1950 mm (Magnutom-400), 1000- 1850 mm (Magnutom-500),



Cutting Mode



4 cutting modes are available:

Cutting mode can also be selected in the *Cutting mode setup* menu when editing the Cut length.

lcon	Cutting Mode	Typical use:
⊕↓	Z Vertical cut	Vertical cut-length can be set. Cutting of larger samples
$\bigcirc \rightarrow$	Y Horizontal cut.	Horizontal cut-length can be set. Cutting of smaller samples
<b>⊕</b> ↓_	<b>Z+Y</b> Combined cut	Vertical and Horizontal cut-length can be set. Cutting of extra-large samples
() TT	AxioCut Combined cut	Combined cut. Cutting of extra-large and extra hard samples:- the cut-off wheel moves downwards in predefined steps

(Maximum cutting capacity is achieved using Y+Z or AxioCut).

For a detailed description of *AxioCut* and how it can be used, see the **Reference Guide** section of the Instruction Manual.

Sampling Method



5 sampling methods are available:

Mode

lcon

	Off	Single cut.
a a a	MultiCut 1	Cut up to 9 slices of equal thickness
¢ b a	MultiCut 2	Cut up to 9 slices of varying thickness
t i c b a 0	MultiCut 3	Cut up to 9 slices of varying thickness counted from a common 0-position
	MultiCut 4	The thickness of the cuts are set by moving the joystick then pressing Enter to set the position of the cut

For a detailed description of the *MultiCut* option and how it can be used, see the **Reference Guide** section of the Instruction Manual.

**Changing the Cut-off Wheel** 

- Remove the nut with a fork spanner (36 mm).
- Remove the flange and the cut-off wheel.
  - Mount the new cut-off wheel.

#### Important

Conventional cut-off wheels based on Al<sub>2</sub>O<sub>3</sub>/SiC abrasives should be placed between two cardboard washers, to protect the cut-off wheel and flanges. For maximum precision with diamond or CBN cut-off wheels, do not use cardboard washers.

Mount the flange and nut. Tighten carefully and close the guard. 

### Clamping the Workpiece

X-Table Models (option)

Position the x-table to allow room for the required number of slices or the size of the workpiece to be cut. If this is not necessary, the tables should be kept close together to support the sample as much as possible during cutting.

Note
The Max load on the X-table is 300 kg / 660 lbs.

All Models

- Move the cut-off wheel backwards to the rear of the cutting chamber for optimum access to the cutting table.
- Clamp the workpiece with the clamping device of your choice e.g. a quick clamping device.
  - Place the workpiece between the clamp and the backstop.
  - Push the clamp towards the workpiece, and lock the quick clamping device with the locking handle.
- Ensure that only one of the quick clamping devices is tight, the other device should only press lightly. Use support tools if the geometry of the workpiece makes support necessary.

Before starting to cut, check that the cut-off wheel or cut-off wheel guard will not come into contact with the clamping device during the cutting process.

When using MultiCut programs, check that the cut-off wheel will not be obstructed by the workpiece or the clamping device during its Xaxis movement.

Cutting on Magnutom	For a detailed description of the <i>AutoCut</i> option and how it can be used, see the <b>Reference Guide</b> section of the Instruction Manual.
<i>Starting the Cutting</i> AutoCut: Off	Position the cut-off wheel by slowly moving it to within 1 - 2 mm of the sample.
	Important Take care when moving the cut-off wheel. If moved too quickly and it comes into contact with the workpiece, it is possible that the cut-off wheel will break.
	<ul> <li>Close the protection hood.</li> <li>Press START <sup>(1)</sup>.</li> <li>The cut-off wheel starts rotating, the cooling fluid starts running and the cut-off wheel will slowly move down into the workpiece at the pre-set Feed speed.</li> </ul>
AutoCut: Simple	<ul> <li>Press START <sup>(1)</sup>.</li> <li>The cut-off wheel will advance towards the workpiece at a maximum speed of 5mm/s.</li> <li>After contact with the workpiece, the cut-off wheel is automatically retracted 2 mm, to be ready for cutting.</li> <li>The cut-off wheel will then continue to move into the workpiece at the pre-set Feed speed.</li> </ul>
AutoCut: Program	<ul> <li>When AutoCut has been pre-programmed:</li> <li>Press START Φ.</li> <li>The cut-off wheel will advance at high speed to the pre-programmed start postion.</li> <li>After reaching this position the speed is reduced until contact is made with the workpiece.</li> <li>After contact with the workpiece, the cut-off wheel is automatically retracted 2 mm, to be ready for cutting.</li> <li>The cut-off wheel will then continue to move down into the workpiece at the pre-set Feed Speed.</li> </ul>

### The Cutting process screen

The Cutting process screen shows information about the cutting process including:

- Cutting Parameters
- Motor Information
- Cutting Mode and a Countdown timer





Feed speed



Motor load

Cutting Mode	Cutting mode MultiCut option off/selected
Timer Countdown	O2:44 An estimate of time remaining for the cutting process. (When AutoStop is selected, a countdown timer will not be displayed.)
Manual Stop	<ul> <li>Magnutom automatically stops cutting when the workpiece is cut through but can be stopped at any time during operation by pressing STOP <sup>∞</sup>.</li> <li>Press STOP <sup>∞</sup> to manually interrupt the cutting process.</li> </ul>
Re-starting cutting	<ul> <li>Press START I to resume cutting.</li> </ul>
Additional cooling To activate jets, turn valve	<ul> <li>When cutting hollow and/or thin workpieces, two flexible water jets are provided for workpiece cooling.</li> <li>Position the cooling jets to the left and right of the cutting area.</li> <li>Turn the valve on the jet to a position parallel to the hose to activate the cooling jets. The cooling liquid will flow as soon as the cutting starts.</li> <li>When cutting is finished, return the valve on the jet to horizontal.</li> </ul>
	Note When using the flexible water jets, cooling liquid is diverted from the integrated water jets positioned over the cut-off wheel. The flexible water jets should not be used when cutting workpieces of large diameter as cooling will be less efficient. They are designed to supply a more localized cooling e.g. the internal surface of hollow workpieces.





### AxioWash

Time

Cut-off wheel movement

The AxioWash cleaning program is an efficient way to automatically clean the cutting chamber. The following parameters can be adjusted: time, cut-off wheel movement or X-table movement (with X-table option).

Can be set to values between 30 sec - 30 min, in steps of 30sec. Default value: 15 minutes

Main me	enu			
Í	A	kioWash		
	Select	AxioWash	mode	5:00
		Off		ment
	Full	Z moveme	ent	
	Reduce	d Z move	ement	
Please remove workpiece and clamping tools before starting the cleaning programme if necessary.				
	Esc Close	Â	Start	
·				

Off:

cut-off wheel arm stays in the current position.

Full Z movement:cut-off wheel arm cycles backwards and<br/>forwards and between the lowest and the<br/>highest Z position.

Reduced Z movement: cut-off wheel arm cycles backwards and forwards and between two Z positions – the z-position of the cut-off wheel arm when AxioWash is started and the highest Z position.

Reduced Z movement is useful as it does not require the clamping tool or workpiece to be removed before AxioWash. Before using this AxioWash function, move the cut-off wheel to the lowest Z position where there will be no collision with clamping tools/ workpiece.

When an X-table is fitted, AxioWash can be programmed to run with X-table movement set to **On** or **Off**.



Clean the cutting chamber.

X-table movement

Cleaning the Cutting Chamber

 Clean the cutting chamber , especially the cutting table the Tslots and the bellows.
 Clean automatically (using AxioWash) and then manually (using the flushing gun).

Automatic cleaning: AxioWash

To start the AxioWash function:

- Remove the workpiece and tools from the cutting chamber.
- Close the adjustable cleaning nozzles.
- Close the protection hood and the side door.
- Press the AxioWash key on the Control Panel to activate the AxioWash PopUp. If required, adjust the AxioWash parameters - time, cut-off wheel movement or X-table movement – then press the AxioWash key again to start cleaning.

The AxioWash program will then run for the preset time.

Manual cleaning

When AxioWash is finished:

- Press the Flush key on the control panel to activate the recirculation pump.
- Lift the flushing gun from its holder.
- Point the flushing gun towards the bottom of the cutting chamber.
- Turn on the water by pressing the rear of the nozzle.
- Clean the cutting chamber thoroughly.
  - Clean the spaces between the T-slots on the cutting table.
  - Clean the space around the cutting table.
  - Clean the bellows thoroughly, in particular the folds where swarf can accumulate.
- Turn off the water by pressing the Flush key again.
- Return the flushing gun to its holder.



### CAUTION

Risk of wet floor/ slipping hazard.

 Do not clean the roof of the cutting chamber or directly on the doors with the flushing gun.

### Maintenance of Cut-off Wheels

Storage of Bakelite Bonded Al<sub>2</sub>O<sub>3</sub> Cut-off Wheels

Maintenance of Diamond and CBN Cut-off Wheels

This type of cut-off wheel is sensitive to humidity. Therefore, do not mix new, dry cut-off wheels with used, humid ones. Store the cut-off wheels in a dry place, horizontally on a plane support.

The precision of diamond and CBN cut-off wheels (and thus the cut) depends on how carefully the following instructions are observed:

- Never expose the cut-off wheel to overload, such as heavy mechanical load, or heat.
- Store the cut-off wheel in a dry place, horizontally on a plane support, preferably under light pressure.
- A clean and dry cut-off wheel does not corrode. Therefore, clean and dry the used cut-off wheel before storing. If possible, use ordinary detergents for the cleaning.

Weekly	The machine should be cleaned regularly, in order to avoid damaging effects to the machine and the specimens from abrasive grains and matal particles.
	<ul> <li>Clean painted surfaces, and the control panel with a soft damp cloth and common household detergents. For heavy duty cleaning, use Struers Cleaner (Cat. No. 49900027).</li> <li>Clean the cover with a soft damp cloth and a common household anti-static window cleaning agent. Do not use harsh or abrasive cleaning agents.</li> </ul>
	<i>Note:</i> Ensure that no detergent or cleaning agent residue is flushed into the cooling unit tank; excess foaming will occur.
Cleaning the Cutting Chamber	<ul> <li>Remove the clamping device(s).</li> <li>Thoroughly clean and lubricate the clamping device(s).</li> <li>Store the clamping device(s) in a dry place or replace on the cutting table after cleaning.</li> </ul>
	<ul> <li>Clean the cutting chamber thoroughly:         <ul> <li>Remove the grates at the rear of the cutting chamber and clean.</li> <li>Clean along the length of the guide shafts with the flushing gun and a brush to remove accumulated swarf.</li> <li>Clean under the cutting table with the flushing gun and a T-slot cleaner (Cat. No. 05486910) or brush to remove accumulated swarf.</li> </ul> </li> </ul>
Cooling Unit	Check the level of the cooling liquid after 8 hours use or at least every week.
Monthly	Deplose the coefficient is the Desire vertice. Coefficient list of
Replacing the Cooling Fluid	Replace the cooling fluid in the Recirculation Cooling Unit at least once a month.
Lubricating the Movement Mechanisms	Magnutom is equipped with an automatic lubrication system for the movement mechanisms.
X-Table Models (option)	<ul> <li>To perform lubrication of the X-table:</li> <li>■ Move the X-table fully to the left, then fully to the right.</li> </ul>

Maintenance of **Clamping Devices** 



*Important.* It is recommended to thoroughly clean and lubricate the Quick Clamping Device and Vertical Quick Clamping Device at regular intervals.

Maintenance of the Clamping Devices is part of the regular yearly Struers Service.

Maintenance of Cutting Table

The stainless steel bands forming the cutting table should be replaced if they have become worn or damaged. The bands are available as spare parts. To allow humidity to escape from the cutting table and chamber, it is recommended to leave the hood open when the machine is not in use.

Lubricating the Cutting Table Part of Struers ServiceGuard

To maintain the optimum performance of Magnutom, lubricate the cutting table at regular intervals (approx every 100 hours). Check the service info displayed on the screen at start-up to monitor the actual number of usage hours. After lubrication of the x-table, note the date and number of service

After lubrication of the x-table, note the date and number of service hours on the Maintenance Log table.

A bottle with oil for lubricating the guide shafts and slide bearings is supplied with Magnutom. When all the oil has been used, refill with eg. Shell TELLUS oil S100.

- Remove the stainless steel bands located over the guide shafts (see illustration).
- Remove the M6 screws.



- Fill the holes with approx. 10ml Oil or until the reservoir is full.
- Re fit the screws.

Lubricating the Guide Shafts (models with optional X-table)

### Lubricating the Drive Shaft Spindle

Part of Struers ServiceGuard

The oil in the spindle of the cut-off wheel drive shaft should be replaced after 1 year or 1,500 hours.

To replace the oil in the spindle:

- Move the drive shaft as far forwards and as low down as possible.
- Place a bowl under the M8 drain plug (on the bottom of the drive shaft).
- Remove the plug from the oil reservoir, then unscrew the drain plug on the bottom and allow the oil to drain.
- When the spindle is empty, replace the drain plug.
- Fill the spindle by pouring 140ml of oil into the reservoir. (Shell Morlina S2 B32 is recommended (requires 2 bottles of Cat. No. R6080852)).
- Replace the plug.



*Important* If the oil in the reservoir appears milky; this indicates that water has entered the oil reservoir. Call a Struers Service Technician to clean the reservoir.

\_

Adjusting the Laser Option Rotating the laser: Adjust the position of the laser so that the laser beam is aligned with the cut-off wheel.



Loosen screw and rotate laser until the laser beam is parallel to the cut-off wheel.



Tighten screw to fix laser in the laser arm.

- The laser arm can be raised to access the cutting chamber from above.



Loosen screw to lift the laser arm. Tighten screw to fix the position of the laser arm.



Rotating the laser arm:

Yearly Inspection of Cover Part of Struers ServiceGuard

The protective cover consists of a metal frame and a composite material (PETG) screen that protects the operator. In the event of damage, the screen will be weakened and offer less protection.

Visually inspect the cover and the screen for signs of wear or damage (e.g. dents, cracks, damage to edge sealing).



*Important* Carry out inspection at more regular intervals if Magnutom is used for more than one 7 hour shift a day.

## Replacing the Screen in the Cover



To ensure its intended safety, the PETG screen must be replaced every 5 years<sup>3</sup>. A label on the screen indicates when it is due to be replaced.

The cover screen should be **replaced immediately** if it has been weakened by collision with projectile objects or if there are visible signs of deterioration or damage.

A label on the cover indicates when the screen is due to be replaced.

Inspection of the cover and replacement of the screen are part of ServiceGuard, the Struers range of service plans.

<sup>&</sup>lt;sup>3</sup> Replacement of the screen is required to remain compliant with the safety requirements in the European standard EN 16089.

Testing Safety Devices

The doors have a safety switch system to prevent the cut-off wheel from starting while the cover is open. Furthermore, a locking mechanism prevents the operator from opening the doors until the cut-off wheel stops spinning.

- Start a cutting process.
- Open the door



### WARNING

Do NOT use the machine with defective Safety Devices. Contact Struers Service.

# **Reference Guide**

### **Table of Contents**

Page

### **1.Advanced Operations**

Configuration Menu	59
User options menu	59
Changing Operation Mode	60
New Pass Code	61
Process options menu	62
Return position:	62
ExciCut	63
AxioWash time	63
Additional Cutting Distance (AutoStop)	63
AutoCut	64
Setting an AutoCut position:	65
Water flow checking	66
Max. Z-axis position	66
Cutting start pos. Y	67
Cutting start pos. Z	67
Chamber flash warning	67
Signalling light	67
Signalling light sound	67
User defined cut-off wheels menu	68
Maintenance menu	69
Service functions	69
Reset configuration	69
Cutting Actions	70
ExciCut	70
Sampling methods	71
MultiCut 1	71
MultiCut 2	72
MultiCut 3	73
MultiCut 4	74
Changing the Cut-off Wheels During MultiCut	76
AxioCut cutting mode	77
Stop Settings	78
AutoStop	78
Additional Cutting Distance (AutoStop)	.78
Cut length	78
Positioning Screen	79
Setting a Relative Zero Position	79
Removing a Relative Zero Position	81
Reference Positions	81
OptiFeed function	82
Clamping Irregular Workpieces	83
Safety Features	84
Optimising Cutting Results	85

2.Struers Knowledge	86
3. Accessories	87
4. Consumables	88
<b>5.Trouble-Shooting</b> Error Messages Messages Errors	89 92 92 92
<b>6.Service</b> Service Information Lubricating Moving Parts	98 99
7.Technical Data	

Mains Cable Specification	100
External Short Circuit Protection	100
Residual Current Circuit Breaker	100
Cutting Capacity	101

### **1. Advanced Operations**

### **Configuration Menu**

A number of different settings and parameters can be set or adjusted from the *Configuration* menu.

There are 3 under menus:

- User options
- Process options
- User defined cut-off wheels

User options	
Display brightness	40
Language	English
Keypad sound	On
Units	Metric
Number of methods	10
Time [hh:mm:ss]	12:59:32
Date [yyyy-mm-dd]	2012-05-31
Default value	

Default value: To re-set the setting to default value, press the F1 key on the Control Panel.

Display brightness:	The brightness of the display can be adjusted to suit individual preferences Can be set between 20-100 Default value: 40
Language:	The language can be set to English (default), German, French, Spanish, Japanese, Chinese, Italian, Russian or Korean.
Keypad sound	Can be set to On or Off
Units:	The units displayed can be set to Metric (mm/s, mm) (default) or Imperial (mil/s, inch).
Number of methods	Up to 300 cutting methods can be stored. The default number: 10.

### User options menu

Wheel selection mode	Can be set to Int Default value: Int Manual: Intelligent:	elligent or Manual telligent the cut-off wheel is selected manually from the selection tree in the <i>Edit method</i> screen a cut-off wheel is automatically suggested based on the material (hardness) selected by user. The recommended RPM and thickness is automatically set.
Operation mode:	It is possible to select three different operation modes:	
	Configuration:	Full functionality
	Development:	No access to parameters in CONFIGURATION menu, except Display contrast
	Production:	Access to START, STOP, Stop position and movement of cut-off wheel, and to Display contrast in the CONFIGURATION menu

Changing Operation Mode

To change the operation mode, go to the *Configuration* menu and then the *Options* menu. Select **Operation mode** to access the *Operation mode* menu.

)	Push knob to select Pass code.
	I SIIIIII
	Enter pass code
	000 <b>ộ</b>
	1 Next position

ţ

-

Use the F1 key and the knob to enter the current pass code (The default pass code is '2750'.):

- Use the F1 key to select digits.
- Turn knob to change the digits and press knob to enter the pass code.

↓





### **New Pass Code**

A New pass code can also be selected from the *Operation mode* menu.





*Note* When a pass code is set the operator has 5 attempts to enter the correct pass code after which the Magnutom will be locked. Re-start Magnutom using the Main Switch then enter the correct Pass Code.



### Important

Remember to make a note of the new Pass code as settings can no longer be changed without the Pass code. **Process options** 

**Return position** 

Excicut amplitude Excicut frequency

Add. cut dist. (Autostop)

AxioWash time

Default value

ExciCut

AutoCut

### Process options menu

### Return position:

After cutting or after pressing STOP  $\odot$ , the return movement of the

Start

2.0 Hz

15:00

3 mm

0ff

0n 1.5 mm

cut-off wheel can be set to three different functions:

- Start: Magnutom automatically retracts the cut-off wheel to the original position of the cut-off wheel, at the time you pressed START  $\oplus$  (default).
- Zero: Magnutom automatically retracts the cut-off wheel to the rear of the chamber (cutting mode Y) or move up in the chamber (cutting modes Z, ZY, AxioCut).
- *Stay:* The cut-off wheel stays in position after the cut is completed.

Top & Rear: Magnutom automatically retracts the cut-off wheel to the rear of the chamber and to its top position.

Select return position	
Start	
Zero	
Stay	
Top & Rear	



#### Important

Use the *Stay* function for bakelite bonded diamond or CBN cut-off wheels, as retraction might destroy the rim of the cut-off wheel.



### Important

When using MultiCut the *Stay* function cannot be used. With MultiCut 4, *Zero* position will automatically be used.

ExciCut	ExciCut cutting action can be set to On or Off. Default value - On Both the <b>amplitude</b> and the <b>frequency</b> of the ExciCut action can be tailored to provide optimal cutting of specific workpieces:	
	ExciCut amplitude	Can be set in the range 0.5 – 4mm, in steps of 0.5 Default value1 5 mm
	ExciCut frequency	Can be set in the range 0.2 – 2.5Hz, in steps of 0.1 Hz. Default value 0.5 Hz
	<i>Tip</i> Change the amplitude and frequency values gradually and in small steps; the ExciCut action is very powerful and making larger adjustments in the values may result in damage to the workpiece or the cut-off wheel.	
	For further details on the <i>ExciCut cutting action</i> , see the section <i>Cutting actions</i> .	
AxioWash time	Can be set in the range 0.30 – 30:00 in steps of 0.30mins. Default value 15 mins	
Additional Cutting Distance (AutoStop)	An additional distance can be set to ensure that the workpiece is completely cut through when the <i>AutoStop</i> is selected. Can be set in the range $0 - 20$ mm, in steps of 1 mm. Default value 2 mm	

AutoCut

Can be set to Off, Simple or Program. Default value – Off



- Off: The operator moves the cut-off wheel towards the workpiece before START is pressed.
  Simple: The cut-off wheel will automatically move towards the workpiece at a maximum speed of 5mm/s when START is pressed.
  After contact with the workpiece, the cut-off wheel is automatically retracted 2 mm, to be ready for cutting. The cut-off wheel will then continue to move down into the workpiece at the pre-set Feed speed.
- *Program:* When AutoCut is set to Program, the *Edit method* menu will change to indicate that a programmed **AutoCut position** is set.

When an AutoCut position is set:

The cut-off wheel will advance at high speed to the preprogrammed start position when START is pressed. **After reaching this position the speed is reduced** until contact is made with the workpiece. After contact with the workpiece, the cut-off wheel is automatically retracted 2 mm, to be ready for cutting. The cut-off wheel will then continue to move into the workpiece at the pre-set Feed speed.

The main purpose of AutoCut is for cutting of large/long workpieces where the cutting arm preferably should be at the far back while loading. The programmable Return Position is set to Top & Rear for maximum access. In order to make the cut as quickly as possible, the start position (just ahead of the workpiece) can be set using AutoCut Program.



*Information* AutoCut is automatically set to Off whilst using MultiCut, to avoid damage to the cut-off wheel and workpiece. Setting an AutoCut position:

From the Edit method menu



### F1 Press the F1 key.

The *AutoCut* pop-up screen will appear (note: AutoCut must be set to Program from the Process options menu)





Use the Joystick to move the cut-off wheel to the required AutoCut start position.

ESC Press **Esc** to leave the menu without changing the AutoCut start position.



Push knob to set the AutoCut start position to the present position of the cut-off wheel.



**Note** AutoCut position is set for the current cutting mode selected. If cutting mode is then changed, the AutoCut position is automatically reset to zero (for safety reasons). When AutoCut position is set to 0.0mm, AutoCut Program will function in

the same way as AutoCut Simple.

	Frocess options		
	Water flow checking Max. Z-axis position Cutting start pos. Y Cutting start pos. Z Chamber flash warning Signalling light Signalling light sound	On 230 mm Off Off Disabled Disabled Off	
Water flow checking	Can be set to On or Off Default value - On When Water flow checking is set checked at the start of and during is below a set value, then the cutt Water flow checking Off can be u Magnutom if the water flow is ade measuring correctly.	to On, the flow of cooling liquid is a cutting process. If the water flow ing process will be stopped. sed to continue operating quate but the flow sensor is not	
Max. Z-axis position	<ul> <li>Can be set to Off or between 100 Default value – 240 mm</li> <li>This can be used to limit the Z mo function is useful when using spe with irregular geometries.</li> <li>To set a <i>Max. Z-axis position</i>: <ul> <li>Mount the clamping tool.</li> <li>Use the joystick to position the position.</li> </ul> </li> </ul>	<ul> <li>Can be set to Off or between 100 to 240 mm</li> <li>Default value – 240 mm</li> <li>This can be used to limit the Z movement of the cutting arm. This function is useful when using special clamping devices or workpieces with irregular geometries.</li> <li>To set a <i>Max. Z-axis position</i>: <ul> <li>Mount the clamping tool.</li> <li>Use the joystick to position the cutting arm in the maximum Z position.</li> </ul> </li> </ul>	

- Read the **Abs. Z**-value from Positioning screen and change the *Max. Z-axis position* in the Process options menu to this value.

Cutting start pos. Y	Can be set to Off or between 0 to 385 mm Default value - Off	
Cutting start pos. Z	Can be set between 0 to 230 cm or Off Default value - Off	
	Setting the Cutting start position is useful when cutting workpieces of the same size. When Start is pressed, the cut-off wheel will quickly move to the set position.	
Chamber flash warning	Can be set to On or Off Default value - Off	
Signalling light	Can be set to On or Off Default value - Off Please refer to the section on <i>Signalling light</i> (option) on page 24 for details.	
Signalling light sound	Can be set to On or Off Default value - Off Please refer to the section on <i>Signalling light</i> (option) on page 24 for details.	

User defined cut-off wheels menu

Up to 10 user cut-off wheels can be defined for use on Magnutom. The parameters for the wheel size, rotational speed and width can be set using the *Edit wheel* menu.



- Size Can be set in the range 300 508mm Default value 508mm
- RPM Can be set in the range 1000 1850 rpm (Magnutom-500) 1000 – 1950 rpm (Magnutom-400) Default value 1500rpm
- Width Can be set in the range 1.0 6.0 mm Default value 3.5 mm

Maintenance menu

The Maintenance menu has 2 under menus

- Service functions
- Reset configuration



Service functions



functions menu. Basic service functions can also be carried out.

Service information and statistics can be viewed in the Service

*Note* The Service functions menus are in English only. Using the same names/ terms is useful when communicating with your local Service Technician or Struers Customer Service department.

Reset configuration

Magnutoms user and process options can be set to default values by selecting **Reset configuration**. (The cutting methods and cut-off wheel database will be unaffected).

Cutting Actions ExciCut ExciCut cutting action is ideal for fast cutting of very hard materials (HV >400). The oscillating movement of the cut-off wheel has two main advantages: there is less risk of damage to the workpiece and less risk of the motor overheating. ExciCut is also recommended when cutting softer, non-hardened materials

ExciCut cutting action is selected as default for each Cutting mode.

The function can be disabled in the Configuration menu.

### Sampling methods

MultiCut 1



(MultiCut can only be selected when an optional X-table is fitted.)

Magnutom's five sampling methods are: Single Cut (MultiCut Off),

MultiCut 1 permits the cutting of several samples of equal width. To set up cutting using MultiCut 1 mode:



MultiCut 1, MultiCut 2, MultiCut 3 and MultiCut 4.

This parameter sets the width of the samples that will be cut.

This parameter sets the number of samples that will be cut.

Select this parameter if you need to make an initial cut, before you start cutting the samples that you need. This cuts a scrap sample, which you will not use. For example, if the workpiece has an uneven edge that would make it unsuitable as a first sample.

This parameter is calculated by the Magnutom to tell you the length of the X-displacement needed to cut your samples, based on the parameter settings.

)	<b>Note</b> Sample width + Thickness of cut-off wheel x No. of samples.	
	Initial cut: ☑ Required X-displacement =	(Sample width + Thickness of cut-off wheel) x (No. of samples +1)
	Initial cut: □ Required X-displacement =	(Sample width + Thickness of cut-off wheel) x (No. of samples)

No. of samples

Initial cut

**Required X-displacement** 

MultiCut 2

Sample width

No. of samples

Initial cut



MultiCut 2 permits the cutting of several samples of different widths.

To set up cutting using MultiCut 2 mode:



This parameter sets the width of the samples that will be cut.

Select this parameter if you need to make an initial cut, before you start cutting the samples that you need. This cuts a scrap sample, which you will not use. For example, if the workpiece has an uneven edge that would make it unsuitable as a first sample.

This parameter indicates the number of samples that will be cut.



Required X-displacement

*Note* Pressing the F1 key in this menu will clear all of the samples and their values.

This parameter is calculated by the Magnutom to tell you the length of the X-displacement needed to cut your samples, based on the parameter settings.

For samples 1 to n: Required X-displacement = (Sample width 1 + Thickness of cut-off wheel) + (Sample width 2 + Thickness of cut-off wheel) + (Sample width n + Thickness o
MultiCut 3



MultiCut 3 permits the cutting of several samples of different widths at different relative distances from the 'zero', or starting position. The distances are manually entered into the Magnutom.

To set up cutting using MultiCut 3 mode:

Cutti	na nasition:			
	ing pooreioni			
Cut :	t zoro nos '			0
No of	f samples:		0	
Req.	X-displaceme	nt:	<b>0.0</b> mm	
	lear all			cotup

Select this parameter to make an initial cut at zero position. Otherwise Magnutom will immediately move to the position for sample 1 and start cutting at that position.

This parameter indicates the number of samples that will be cut.



This parameter is calculated by the Magnutom to tell you the length of the X-displacement needed to cut your samples, based on the parameter settings.

Required =	The last relative cutting position entered
------------	--

Cutting position

Cut at zero position

No. of samples



Required X-displacement

MultiCut 4

— • <b>-</b> •
----------------

MultiCut 4 permits the cutting of several samples of different widths at different relative distances from the 'zero', or starting position. The distances are entered into the Magnutom by using the X-table to position the workpiece under the cut-off wheel where the sample is to be cut and then recording this position.

0

0.0 mm

0.0 mm

Edit

Cut-off wheel position is also recorded so different starting-height positions are possible.

MULTICUT 4

To set up cutting using MultiCut 4 mode:

Setting the cutting parameters	<ul> <li>Use the joystick to move the X-table to the position where the first cut is to be made.</li> <li>Position the cut-off wheel about 2 mm above the workpiece.</li> <li>Press the knob to insert the current position as the cutting position.</li> <li>Repeat the steps to insert the cutting positions for all samples.</li> </ul>
X-pos.: Y-pos. or X-pos.: Z-pos.	The different cutting positions of both the X-table and the cut-off wheel . Y or Z position (depending on which cutting mode is selected). X-pos.: Y-pos.: cutting mode Y X-pos.: Z-pos.: cutting modes Z, ZY, AxioCut
Go to start	This key moves the X-table until the cut-off wheel is in the start position.
No. of samples	This parameter indicates the number of cuts that will be made.

No of samples:

Absolute X-pos.:

Absolute Z-Pos.:

E Remove pos.

Absolute Y-pos. / Absolute Z-pos. The current Y or Z position of the cut-off wheel (depending on which cutting mode is selected). **Absolute Y-pos.:** cutting mode Y **Absolute Z-pos.:** cutting modes Z, ZY, AxioCut

Note

Pressing the F1 key in this menu will clear all of the samples and their values.

#### Important

When MultiCut 4 is selected and F1 is pressed, Stop position is also set to *Auto*.

Stop position can be changed to a specific position. However, this position is defined from the position of the cut-off wheel at the first cutting position. If the cut-off wheel is in a different starting position for the following cuts, the cut length will remain the same as for the first cut and will be more or less than the cut length required. . *Auto* setting is therefore recommended.

Example for Z cutting mode:

MULTICUT 4	
X-Pos.: Rel. Z-Pos.:	Height:
1. 0.0 mm 0.0 mm	100 mm 📗
2. 22.4 mm 16.0 mm	90 mm
3. 47.8 mm 28.0 mm	105 mm
4. 75.0 mm 10.0 mm	80 mm 📋
Go to start	
No of samples:	4
Absolute X-pos.:	75.0 mm
Absolute Z-Pos.:	10.0 mm
E Remove pos.	Edit 🕒
	mouricon set

AutoCut: Off

After each cut the cutting arm will:

- Move to top position (to ensure there is no collision with the workpiece during movement in the X direction).
- Then move to the set X position.
- Move quickly to the predefined cut-off wheel position.
- Starts cutting with the predefined Feed speed.

AutoCut: Simple or Program

#### After each cut the cutting arm will:

- Move to top position (to ensure there is no collision with the workpiece during movement in the X direction).
- Then move to the set X position.
- Move quickly to the predefined cut-off wheel position then search for workpiece with relatively high speed.
- Starts cutting with the predefined Feed speed.

#### Changing the Cut-off Wheels During MultiCut

When cutting several specimens it may be necessary to change the cut-off wheel during a MultiCut program.

- Press STOP to manually interrupt the cutting process. The cut-off wheel arm goes to the start position.
- Change the cut-off wheel.
- Position the cut-off wheel approximately 1- 2mm above the workpiece.

## !

#### Important

Do NOT move the cut-off wheel back to its position immediately before STOP was pressed.

When START is pressed, Magnutom registers this as a new start position and will return to this position after the cut is finished. The cut-off wheel will be damaged if it is still inside the workpiece when the X-table moves to the next cut position.

- Press START  $\diamondsuit$  to resume cutting.
  - A Pop-up will appear: Press Enter to continue with the programmed MultiCut batch.



#### Note

The cut-off wheel will start rotation and move towards the workpiece. If desired, the process can be sped up by using the joystick (feed speed will be 3 times the preset feed speed).

#### AxioCut cutting mode



AxioCut cutting mode offers cutting of large and very hard materials. The vertical cut length and height of the workpiece can be set as well as the depth of each downwards step (in the example below, set to 5 mm).



When in AxioCut cutting mode, the vertical feed speed is the same as the horizontal feed speed.



*Note* AxioCut does not function when AutoCut is selected. In practice, AxioCut cutting mode should not be used with MultiCut.

AutoStop.

Stop Settings



AutoStop

Cut length

Additional Cutting Distance (AutoStop)

*Note* When AutoStop is selected, a countdown timer will not be displayed.

There are two ways to set the stop position: Setting a Cut length or

When the *AutoStop* function is selected, the machine automatically stops when the workpiece has been cut through. *AutoStop* is based on changes in the cutting motor load where a decrease in load

However, Magnutom has a very powerful motor so small changes in the load may not be detected. As a result the AutoStop function may not detect when the workpiece has been cut through. This may particularly be the case when cutting soft materials, smaller workpieces, pipes or workpieces with changing cross-section or

If *AutoStop* does not work properly use the *Stop position* function instead or set an Additional Cutting Distance to compensate.

An additional distance can be set to ensure that the workpiece is completely cut through when the *AutoStop* function is used. This is

An Additional cutting distance is set in the Process options menu.

Cut length is used when a specific stop position is desired. When cutting pipes or other workpieces with changing cross sections, the cut-off wheel may retract before the workpiece has been cut through. To overcome this, a specific Cut length can be set.

 Clamp the workpiece and position the cut-off wheel just above the workpiece.

This position is automatically set to 0 (zero). Accordingly, as soon as Start has been pressed, the actual position of the cut-off wheel, becomes a relative starting point (zero), from where the cutting depth is calculated.

 Select the STOP parameter and set the desired stop position using the knob.

Magnutom will now stop when it reaches the pre-set stop position.

Remember to take the wear of the cut-off wheel into consideration.

Note

important when using the MultiCut option.

indicates that the workpiece has been cut through.

when cutting at low Feed Speeds.

#### **Positioning Screen**

The *Positioning* screen shows the Y and Z positions of the cutting arm and the X position of the X-table (with x-table option). The screen is displayed for 3 seconds when the joystick has been

- activated.
   To keep the screen displayed (for up to 15 minutes), press the F1 key.
- To close the screen, press **ESC**.

Absolute and Relative positions usually have the same value, but the user can define a new relative position, e.g. if there is a specific point on the workpiece that the user wants to set as 'zero position'.



Setting a Relative Zero Position

To define a new relative zero position:

Move the cursor to the rel.: position, then press **Enter**. This position (in the current example, Y position 18.0 mm) is now set as the new relative zero.



When the cutting arm is moved, the relative position now shows the value in regard to the previously defined zero position.



To return to the previously defined zero position:

Move the cursor to the "Go to relative zero position" icon and press Enter.



Removing a Relative Zero Position

■ Move the cutting arm to *abs.:* zero position (0.0 mm).



■ With the cursor highlighting the *rel.*: position, press **Enter**.



**Reference** Positions

Magnutom will perform a search for reference positions for every 20th start up.

To perform an additional reference search:

- Switch Magnutom off.
- Press the Emergency STOP button then switch Magnutom on (a message will appear to inform that the Emergency STOP is activated).
- Release the Emergency STOP and a pop-up will appear. Press Enter to start the search for reference positions.

**OptiFeed function** 



During cutting, Magnutom continuously measures the load on the cutting motor and force in the cutting arm.

The factors that determine the load are the shape and properties of the workpiece.

The OptiFeed function:- Whenever the maximum force or motor load is reached (150%), Magnutom will automatically reduce the Feed speed. As soon as the force or load drops below the set limit, the speed will be increased to the original setting.

The figure below illustrates the increase in force as the cut-off wheel gets closer to the centre of a round workpiece.



Clamping Irregular Workpieces Irregular workpieces without plane clamping surfaces must be clamped using special clamping tools, as the workpieces must not move during cutting. This could result in damage to the cut-off wheel or to the sample itself. Use the T-slots to mount the special clamping tools. Struers offers a large selection of Clamping Tools (See Accessories).

will cut the smallest possible cross-section.

Removing the left cutting table (models with fixed cutting table)

If required, the left side of the cutting table can be removed to make room for very large or irregular workpieces.

Safety Features Doors	The doors are equipped with a safety switch to prevent the cut-off wheel from starting while the doors are open. Furthermore, a locking mechanism prevents opening of the doors before the cut-off wheel has come to a complete stop.
	The doors can only be opened when the machine is connected to a power supply and the main power switch is ON. To open the doors when the power is not connected, use the triangle key to release the safety lock. Remember to re-activate the safety lock release before operating Magnutom.
	The PETG safety glass in the windows are designed to withstand flying fragments of cut-off wheel or samples. A sticker on the window indicates when the glass should be replaced. Replace immediately if the window is damaged or cracked.
Hold-and-run operation	When moving the cut-off wheel whilst the doors are open, the Hold- to-run button on the front of Magnutom must be pressed and held immediately before the joystick is operated.
Motor overload	The motors of Magnutom are protected against overload. Should the motors overheat and/or overload, the motors will disengage until a normal temperature has been obtained.
Lockable Main Switch	When the main switch is turned Off, a hole in the lock's housing can be used to attach a padlock or a plastic strip. This is particularly useful during Servicing as it will prevent Magnutom from being switched on before the Service is complete.

#### **Optimising Cutting Results**

The following table shows possible answers to a number of common questions:

Optimising the Cutting Results		
Question	Answer	
How can I avoid discoloration or burning of the sample?	Use a lower Feed Speed.	
	Change the cut-off wheel as the hardness of the present cut-off wheel may be inappropriate for the hardness of the sample. *)	
How can I avoid burrs?	Use a softer cut-off wheel. *)	
	Clamp the workpiece securely at the right hand clamping device. Tighten the left hand clamping device just enough to prevent the workpiece from moving when being cut.	
How can I avoid the cut-off wheels wearing too quickly?	Use a lower Feed Speed, a different cutting mode or a harder cut-off wheel. *)	
How can I achieve faster cutting?	Position the workpiece so as to cut the smallest possible cross-section. Use a high Feed Speed.	

\*) Please refer to the Selection Guide in the <u>Struers Cut-off Wheels brochure</u>.

## 2. Struers Knowledge

Materialographic sectioning is where most microstructure analysis begins. A good understanding of the abrasive cutting process can help to select suitable clamping and cutting methods and thereby ensure the high quality cut. Minimizing cutting artefacts will help the remaining materialographic process and act as a good base for efficient and high quality preparation.



For further information, see the section on <u>*Cutting*</u> on the Struers website.

## 3. Accessories

Please refer to the *Magnutom brochure* for details of the range available.

**Clamping Tools** 

Please refer to the *Struers Clamping Tools brochure* for details of the range available.

## 4. Consumables

The use of Struers consumables is recommended. Other products (e.g. coolants) may contain aggressive solvents, which dissolve e.g. rubber seals. The warranty may not cover damaged machine parts (e.g. seals and tubes), where the damage can be directly related to the use of non-Struers consumables.

When Wheel selection mode is set to Intelligent, a cut-off wheel is automatically suggested based on the material (hardness) selected by the user.



Alternatively, please refer to the Selection Guide in the *Struers Cut*off *Wheels brochure*.

Specification	Cat no.
Corrozip	
Additive for Cooling Fluid To protect the machine from corrosion and to improve cutting and cooling qualities.	
11	49900045
5	49900046
Grease for maintenance/lubrication of the spindle	16080802
Oil for maintenance of cutting table	16080845

#### Other Consumables

Cut-off Wheels

## 5. Trouble-Shooting

Error	Explanation	Action
Machine Problems		
Water leaking.	Leak in recirculation water hose.	Check the hose and tighten the hose clamp.
	Water overflow in the Cooling fluid tank.	Remove the excess cooling fluid from the tank.
Samples or cutting chamber corroded.	Insufficient additive for cooling fluid.	Add Struers Additive for cooling fluid to the cooling fluid, using the correct concentration. Check with a refractometer. Follow the instructions in the Maintenance section.
	The machine is left with closed protection hood.	Leave the protection hood open to let the cutting chamber dry.
Quick-clamping device unable to hold the workpiece.	The quick-clamping device is not in balance.	Adjust the screw underneath the clamping column. Use a 3mm Allen key.
	Clamping heart worn.	Call a Struers Service Technician.
Door will not close	There is an obstruction in the cutting chamber.	Remove the obstruction.
Machine is locked	Incorrect Pass Code used.	Re-start machine using Main Switch.
		Enter the correct Pass Code.
		If error remains, contact a Struers Service Technician.

Error	Explanation	Action
Cutting Problems		
Discoloration or burning of the sample.	The hardness of the cut-off wheel is inappropriate for the hardness / dimensions of the sample.	See Consumables section, Cut-off Wheels.
	Inadequate cooling.	Check that there is enough water in the recirculation cooling unit.
		tray.
	Feed Speed or RPM too high.	Reduce the Feed Speed or reduce the RPM.
Unwanted burrs.	Wheel too hard.	Reduce RPM to make wheel softer, or change to different wheel: See Consumables section, Cut-off Wheels.
	Insufficient support of workpiece.	Add further support to the workpiece.
The cutting quality differs.	Cooling liquid hose clogged.	Clean the cooling liquid hose and the cooling tubes. Check the water flow by turning the cooling valve to cleaning position.
	Insufficient cooling liquid.	Refill tank with water. Remember to add Struers Additive.
The cut bends to one side.	Feed Speed is too high.	Reduce the Feed Speed.
The cut-off wheel breaks.	Incorrect mounting of the cut-off wheel.	Check that the centre-hole has the correct diameter.
		Check cardboard washer on both sides of the cut-off wheel. The nut must be tightened properly.
	Incorrect clamping of the workpiece.	Make sure that only one of the quick-clamping devices is tight. The other device should only press lightly.
		Use support tools if the geometry of the workpiece makes support necessary.
	Wheel too hard.	See Consumables section, Cut-off Wheels.
	Feed Speed too high.	Reduce the Feed Speed.
	Inadequate cooling.	Check that there is enough water in the recirculation cooling unit.
		Check the cooling liquid hoses.
I he cut-off wheel wears down too	Feed Speed too high.	Reduce the Feed Speed.
	Insufficient cooling.	Check that there is enough water in the recirculation cooling unit. Check the cooling liquid hoses.
	The cut-off wheel is too soft for the task.	See Consumables section, Cut-off Wheels.
	Magnutom vibrates (worn bearings).	Call a Struers Service Technician.

Error	Explanation	Action
The cut-off wheel does not cut through the sample.	Incorrect choice of cut-off wheel.	See Consumables section, Cut-off Wheels.
	Cut-off wheel worn.	Replace the cut-off wheel.
	The cut-off wheel gets caught in the workpiece.	Support the workpiece and clamp it on both sides of the cut-off wheel in such a way to allow the cut to stay open.
	Incorrect choice of Cutting Mode. AxioCut is intended for large workpieces.	See section on Operation, Cutting Mode.
The workpiece breaks when clamped.	The workpiece is brittle.	Place the workpiece between two polystyrene plates.
		NB! Always cut brittle workpieces very carefully.
The sample is corroded.	The sample is not resistant to water.	Use a neutral liquid as cooling fluid or cut without using cooling fluid at all. DO NOT USE AN INFLAMMABLE
	cutting chamber for too long.	when you leave the machine.
	Insufficient additive for cooling fluid.	Add Struers Additive for cooling fluid to the cooling liquid in the correct concentration. Check with a refractometer. See Maintenance section.
<i>AutoStop</i> does not stop the cutting action.	The cross-section of the workpiece is too small or irregular to detect a change in load.	Use the Stop position function.
AxioCut cutting mode is not working.	AxioCut does not function when AutoCut is selected.	Go to the <i>Process options</i> menu and set AutoCut to <b>Off.</b>

Error Messages	Error messages are divided into two classes: Messages and Errors
Messages	Messages are intended to inform the operator of the machine's progress and advise about minor operational errors.
_	

Errors

Errors must be rectified before cutting can be continued.

Message		Explanation	Action
Information #4 It is not possible to create new position. Required space must be equal or higher than wheel thickness plus 500 microns.	ŀ		Move X-table to increase the distance between the positions.
Information #5 Selected MULTICUT mode cannot be combined with actual cutting mode.	5		Change to another Cutting Mode.
Information       #40         Selected cut size exceeds available cutting capacity.       #4         Press Enter to set max. available cutting capacity.       Esc Cancel Ok	10	The selected cut size is not possible due to the current position of cut-off wheel.	Press <b>ENTER</b> to set cut size to the maximum available or press Esc to return to the method.
Information #41 Cutting process can not start because max. cutting position is reached. Please move the cutting arm from max. position. Ok	1	The cutting arm is very close to the end position, therefore the cutting process cannot be started.	Move the cutting arm from the end position and re-start the process.

Message		Explanation	Action
Information#47Selected batch job can't be executed, because the x-table displacement is fully used.Esc CancelDecrease no. of samples	#47	Multicut is activated but the current position of the X-table is not enough to make all the cuts.	Press <b>ENTER</b> to automatically decrease the number of cuts or press <b>Esc</b> to return to edit the method.
Error       # 14         Main motor overheated!         Please wait until motor temperature         decreases before continuing cutting.         Motor temperature:       0 °C         Press Enter to start cooling. Stops         automatically when temp. < 100°C.         Or press Esc. to cancel.         Esc. Cancel         Esc. Cancel	#14	Motor temperature is > 150°C during start of process or Motor temperature is > 170°C during cutting process.	Press <b>ENTER</b> and wait until the motor has cooled.
Error #16 Main motor rotation lost. Process halted! Please try to restart cutting process. Alarm register: Oh Error register: Oh Error register: Oh	#16	Cutting motor has stopped rotating. A fluctuation may have occurred in the mains power supply.	Check the mains power supply then Re-start. If error remains, contact a Struers Service Technician.
Error #17 Main motor overloaded! Process halted! Please decrease feed speed before continuing. © Ok	#17	Magnutoms OptiFeed feature automatically reduces feed speed if the load on the motor exceeds a set value. However, under some circumstances (e.g. initial feed speed too high or joystick movement activated during cutting) the Optifeed cannot decrease the feed speed fast enough and the motor current reaches a set limit.	Reduce the set feed speed before re-starting the cutting process.

Message		Explanation	Action
Warning #46 The specimen has not been detected. Specimen must be placed in the cutting trajectory or feed speed is too low.	#46	AutoCut is activated and the workpiece has not been detected during the cutting trajectory. Detection of workpiece is based on a change in the cutting motor current. The idle cutting motor current must increase by a set value upon contact with the workpiece. At low feed speeds, the difference between idle current and cutting current is very low so the workpiece is not detected	Increase the feed speed and re-start the cutting process.
Warning # 53 Pressure of cooling water is too low! Process halted! Please check the cooling system.	#53	Water flow is too low at the beginning of a cutting process.	Press <b>ENTER</b> to abort the cutting process. Then check the cooling system. If error remains, contact a Struers Service Technician.
Warning       #48         Pressure of cooling water is too low!       Please check the cooling system.         Please check the cooling system.       Esc Continue	#48	Water flow is too low at the beginning of a cutting process.	Press ENTER to abort the cutting <b>process</b> . Then check the cooling system. or The water sensor may be malfunctioning. Check if there is sufficient water flow, press Esc to continue with cutting process. If error remains, contact a Struers Service Technician.
Warning       # 102         Machine found a problem during         Power On Self Testing. Machine can         continue but specified accessory         will be disabled.         Reason:       Automatic X-table drive         offline!         Cacept	#102		Re-start. If error remains, contact a Struers Service Technician.

Message		Explanation	Action
Warning # 108	#108	Force sensor is not detected.	Re-start.
Force sensor signal is out of range. Collision with flange and detection of sample by cut-off wheel can not be provide. If problem persists contact Struers technical support. Ok	4	Magnutom will operate without this sensor but some features, e.g. detection of workpiece will not function.	If error remains, contact a Struers Service Technician.
Warning# 111Force sensor is disabled!Detection of specimen and limitation of max. cutting force can not be provide.	#111	The force sensor is disabled in the service menu. Magnutom will operate without this sensor but some features, e.g. detection of workpiece will not function.	Contact a Struers service technician to re-enable the force sensor.
C Ok			

Error		Explanation	Action
Error #23 Searching of reference position failed. Movement with cut-off wheel can be dangerous!	#23	A problem has occurred during search for reference positions and the procedure has been interrupted. Attention!:	Press <b>ENTER</b> to acknowledge this message.
Contact service technician if problem persists. Ok	ļ	to enable some movement with the joystick but a cutting process cannot be started. Searching of reference position is required again after pext	
		switching on.	
Error #24 Y-servo motor error detected during	#24	General Y-servo motor error detected during movement.	Press <b>ENTER</b> to clear this error inside servo motor to re- enable the movement.
movement.			Re-start.
Reason code: Oh Error status register: 1h Cok			If error remains, contact a Struers Service Technician. Please make a note of the <i>Reason code</i> and <i>Error</i> <i>status register</i> displayed on Magnutom.
Error #25 Z-servo motor error detected during	#25	General Z-servo motor error detected during movement.	Press <b>ENTER</b> to clear this error inside servo motor to re- enable the movement.
movement.			Re-start.
Reason code: Oh Error status register: 1h 🕑 Ok			If error remains, contact a Struers Service Technician. Please make a note of the <i>Reason code</i> and <i>Error</i> <i>status register</i> displayed on Magnutom.
Error #35 X-table actuator error detected during movement. Please check and remove all potential obstacles from X trajectory.	#35	An over current is detected during movement of the X-table. Usually due to an obstacle in the X-tables trajectory (e.g. sample is squeezed or the X-table is blocked).	Remove any obstacles then press <b>ENTER</b> to clear the error.

Error		Explanation	Action
Error #36 Excessive cutting force or collision with flange detected! Please decrease cutting speed or remove all potencial obstacles from the trajectory. Ok	#36	Excessive force is detected during a cutting process. This can happen if flange or cut- off wheel guard collides with a cutting sample.	Remove any obstacles from cutting direction then press <b>ENTER</b> to clear the error.
Error #51 Excessive temperature detected in Y-servo motor.	#51 #52	An over temperature (> 73°C) is detected during movement with Y or Z servo motor.	Press <b>ENTER</b> to clear this error inside servo motor to re- enable the movement. Turn Magnutom off and allow the servo motor to cool down.
Error status register: 1h COX Error #52 Excessive temperature detected in Z-servo motor. Error status register: 1h COX			

## 6. Service

#### **Service Information**

Struers recommends that a regular service check be carried out yearly (or after every 1,500 hours of use, if this is sooner). Struers offers a range of comprehensive maintenance plans to suit the requirements of our customers. This range of services is called **ServiceGuard.** 

The maintenance plans include equipment inspection, replacement of wear parts, adjustments/calibration for optimal operation, and a final functional test.

Information on total operation time and servicing of the machine is displayed on the screen at start-up:

<b>Magnutom 500</b>	Version 0.75 X	
<u>SERVICE INFO:</u> Total operating time: Time since last service: Time until next service: Memory module installed	0 h 0 h 1500 h	

A pop-up message will alert the user that the recommended service interval has been exceeded.

Lubricating Moving Parts Part of Struers ServiceGuard

Magnutom is equipped with an automatic lubrication system for the movement mechanisms.



Grease cartridge

The grease cartridge will need to be replaced after 1,500 hours of service.

Replacing the grease cartridge is part of the Struers **ServiceGuard** program.

## 7. Technical Data

Please refer to the <u>Magnutom brochure</u> for the general Technical Specifications.

Mains Cable Specification	Wire gauge	EU: Canada, USA:	10mm <sup>2</sup> AWG8		
	Cable diameter	Max. 20 mm			
External Short Circuit Protection	Magnutom must be Suggested fuse size	protected with e	xternal fuses.		
Residual Current Circuit Breaker	<u>Electrical Installations with RCCB</u> For Magnutom connected to electrical installations with RCCB, <b>a type B time delayed, 30 mA RCCB is REQUIRED</b> (ref. EN 50178 / 5.2.11.1).				
	<u>Electrical Installations without RCCB</u> The equipment must be protected by an insulation transformer (double-wound transformer).				
	<i>Important:</i> Local standards may overrule the recommendations for the main supply cable. If necessary, please contact a qualified electrician to verify which option is suitable for the local installation setup.				
Operating environment	Surrounding temperature	5-40°C/41-104°F			
	Humidity	0-95% RH non c	ondensing		



**CUTTING CAPACITY** 



The graph shows the projected cutting capacity under the following conditions:

- A new cut-off wheel.
- The workpiece is laid directly on the cutting table, with overhang where appropriate.
- Vertical clamping is used.
- Using ZY or AxioCut cutting mode

The actual cutting capacity depends on the sample material, cut-off wheel and clamping technique.

Clamping the Workpiece ■	Place the workpiece between the clamp and the back stop on the right hand cutting table. Push the clamp towards the workpiece and lock the quick- clamping device with the locking handle.
Starting the Cutting ■	Position the cut-off wheel. Close the protection hood. Press START $\Phi$ . The cut-off wheel starts rotating and the cooling water starts running.
Stopping the Cutting	
Automatic	Set a Cut length or set to AutoStop in the <i>Edit method</i> menu.
Manual Stop ■	Press STOP $\heartsuit$ and the cutting process stops. The cut-off wheel stops rotating and the cooling water stops.
Setting the Cutting Parameters	In the Cutting Display, set each of the cutting parameters using the TURN/PUSH knob.
Changing the Cut-off Wheel ■ ■ ■	Remove the nut with a fork spanner. Remove the flange and the cut-off wheel. Mount the new cut-off wheel. Mount the flange and nut. Tighten carefully and close the guard.
Cleaning the Cutting Chamber AxioWash ■	Press the AxioWash key to clean the cutting chamber.
Manual	Point the flushing gun towards the bottom of the cutting chamber. Turn on the water by pressing the Flush key on the Control Panel. Clean the cutting chamber thoroughly. Turn off the water by pressing the Flush key again. Return the flushing gun to its holder.



English	<b>Declaration of Conformity</b>	
Manufacturer	Struers ApS Pederstrupvej 84 DK-2750 Ballerup, Denmark Telephone +45 44 600 800	
Herewith declares that	Name:         Magnutom-400/-500           Cat. No.:         06276146, 06276246 / 06146146, 06146246           Function:         Cut-off machine           Type No.:         627 / 614	
fulfils all the relevant provision	s of the	
Machinery Directive 2006/42/EC	according to the following standard(s): EN ISO 12100:2010, EN ISO 13849-1:2015, EN ISO 13850:2015, EN 60204-1:2006/AC:2010, EN 61010-1:2	010.
and is in conformity with the:		
EMC Directive 2014/30/EU	according to the following standard(s): EN 61000-6-2:2005, EN 61000-6-3:2007/A1:2011, EN 61000-6-4:2007/A1:2011, EN 61000-3-3:2013, EN 61326-1:	2013.
RoHS Directive 2011/65/EU	according to the following standard(s): EN 50581:2012.	
Supplementary Information	The equipment complies with the following standards: NFPA70:2014, NFPA79:2012, FCC 47 CFR PART 15.	
The above has been declared a	ccording to the global approach, module A.	
Authorized to compile the	Technical File:	
Alevs Trense	_	
Klavs Tvenge Director of Business Development Struers ApS		
Pederstrupvej 84 DK-2750 Ballerup, Denmark		Date of Issue: 2017.10.05



# **Magnutom** Spare Parts and Diagrams



Manual No.: 16147001 Revision A

Date of Release 2018.01.16



## Always state *Serial No* and *Voltage/frequency* if you have technical questions or when ordering spare parts.

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations:

Instruction Manuals: Struers Instruction Manuals may only be used in connection with Struers equipment covered by the Instruction Manual.

Service Manuals: Struers Service Manuals may only be used by a trained technician authorised by Struers. The Service Manual may only be used in connection with Struers equipment covered by the Service Manual.

Struers assumes no responsibility for errors in the manual text/illustrations. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the present version of the equipment.

The contents of this manual are the property of Struers. Reproduction of any part of this manual without the written permission of Struers is not allowed.

All rights reserved. © Struers 2018.

Struers Pederstrupvej 84 DK-2750 Ballerup Denmark Telephone +45 44 600 800 Fax +45 44 600 801

## **Spare Parts and Diagrams**

#### Table of contents

Drawing

Drawings	
Magnutom-500 XT, complete	16140006J
Magnutom 500, complete (2 pages)	16140010M
Base with bottom covers (2 pages)	16140012I
Base with Cutting chamber (2 pages)	16140013L
Base w .Cutting chamber, Ass.	16140014L
Base with both movements	16140015J
Base with Y movement(2 pages)	16140018G
Base with legs and drain	16140019E
Ball Screw Complete, Y-movement	16140020D
Ball Screw Complete, Z-movement	16140021D
Servomotor Base Z-movement, assy	16140022G
Motor Base, assembled	16140023D
Lubrication Manifold, assembled	16140024C
Spindle, Assembled (2 pages)	16140030K
Wheel Guard, Assembled	16140031H
Cutting Table f. X-table (2 pages)	16140036E
Left Net, assembled	16140037B
Spindle lock, Assembled	16140039B
Side door, assembled	16140040G
Side door bottom labyrinth, assy	16140041A
Right Net, assembled	16140042B
Front door, assembled	16140045G
Chamber plates, assembly	16140048H
Chamber Bellows, assembled	16140050E
Upper Profile, Bellows, assemble	16140053A
Sample barrier Left, assembled	16140058D
Protection Shield 1,assembled	16140061B
Protection Shield 2,assembled	16140062C
Protection Cover, assembled	16140064G
HW box and hoses, assembled	16140068D
Front panel, assembly	16140071H
HW box and parts, assembly (2 pages)	16140075J
Power supply, assembled	16140076D
X-table Connection Box	16140077F
Plate for Valves, Assembled (2 pages)	16140095I
Water pistol	16140097D

Magnutom Spare Parts and Diagrams

Diagrams	
Block Diagram	16143050F
Circuit Diagram (9 pages)	
Main voltages, page 1	16143100H
Grounding, page 2	16143100H
YZ servomotors, page 3	16143100H
Safety controller, page 4	16143100H
Control board, page 5	16143100H
Accessories, page 6	16143100H
Door Guards, page 7	16143100H
X-Table box, page 8	16143100H
Safety related circuits, page 9	16143100H
Air diagram Magnutom 500	16141001B
Variant overview	16147600B

Some of the drawings may contain position numbers not used in connection with this manual.
The following is a list of the spare parts that may need replacement during the lifetime of the equipment.

To check the availability of other replacement parts, please contact your local Struers Service Technician. It may help identify the part by referral to its position number on the assembly drawings included in this manual.

	Spare	part list for Magnutom			
Drawing	Pos.	Spare Part	Ref.	Cat no:	
16140006		Magnutom-500 XT, complete			
	30	Seal, assembled		16140088	
	200	Spindle Housing Sealing		16140319	
	210	Protection Stop		16140384	
16140015		Base with both movements			
	110	MotorAM132M/A4B53x400VD/50Hz16	M1	2ME13211	
	200	Shaft Ring SD_45x55x4		2IS30045	
	210	Shaft Ring SD_50x62x5, 4 pcs		2IS30055	
16140018		Base with Y movement			
	80	Buffer, 2 pcs		16140389	
	90	Buffer, 4 pcs		16140413	
	120	Caged Ball LM Guide THK, 2 pcs		2BF41076	
	140	Coupling Rotex GS28,d22-d25,64		2JH02864	
	150	Gearbox i=5 Technoingranaggi		2MG40070	
	160	Servomotor JVL MAC800-D2,230V	M2	2MM51755	
	450	MAC00-FC41 Communication module		2PA00400	
16140019		Base with legs and drain			
	50	Cutting table plate long		16140352	
	60	Cutting table plate Short 22 pcs.		16140340	
	220	O-ring 72 NBR 872. 28.00-2.00		21020057	
	250	Nut M20 A2 DIN934, 6 pcs		2TA10200	
	350	GEKA coupling 3/4 int. thread		2NF60002	
16140010		Magnutom complete			
	210	Ball guide WSQ-16		2BF27416	
16140020		Ball Scrow Complete V movement			
10140020		Dail Screw Complete, 1-movement			

	Spare part list for Magnutom			
Drawing	Pos.	Spare Part	Ref.	Cat no:
	10	Spindle, Y-movement		16140150
16140021		Ball Screw Complete, Z-movement		
	10	Spindle, Z-movement		16140156
	50	Ball screw housing Hiwin BK-30		2BS00031
16140022		Servomotor Base Z-movement, assy		
	30	Gearbox APEX i=5PE090-005_JVL		2MG40901
	50	Coupling Rotex GS28,d22-d25,64		2JH02864
	80	Bearing Spherical SKF GE_30_C		2BK20030
	90	Servomotor JVL MAC800-D5,230V, (1 pc)	M3	2MM51756
	100	MAC00-FC41 Communication module		2PA00400
16140023		Motor Base, assembled		
	50	Bearing Spherical SKF GE_30_C, 5 pcs		2BK20030
	90	Nut Housing, 3 pcs		16140209
16140024		Lubrication Manifold, assembled		
	20	Grease Nipple M8x1 Zn, 6 pcs		2GN01011
16140030		Spindle, Assembled		
	80	Outer Flange		16140314
	90	Fixed Flange		16140316
	150	Inner Ring, machined		16140324
	160	RetainingRingForBoreJ85DIN472, 2 pcs		2ZL20850
	170	Rad. shaft Seal SKF CR45x85x10, 3 pcs		2IS10050
	190	Nut M24 DIN934 A2 Left hand		2TA70240
	220	O-ring 40x2 NBR 60, 4 pcs		21020074
	230	SphericalRollerBear.SKF22209E, 4 pcs		2BK20045
	240	Coupling Rotex 65,d38-d38,64Sh, 2 pcs		2JH06564
	270	V-ring VA65		2IV10065
	320	Oil-MOGUL LK32, 2 pcs (140 ml)		2LS00033
	360	SightWindowGN743.4-11-M16x1.5B		2NG17434
	380	Shaft Seal REINZ 81-26392-10		2IP22457
16140031		Wheel Guard, Assembled		
	40	Hose PVC w. steel coil ø12/ø18		2NU21218
	50	Hose PVC w. steel coil d18-d25		2NU21226

	Spare part list for Magnutom			
Drawing	Pos.	Spare Part	Ref.	Cat no:
	60	Snap-lock tube, assembly		16140033
	80	Spray nozzle TP14 120 316 PP.		2YD20004
	100	Ballofix valve 1/2-1/2in		2YH10603
16140033		Snap-lock tube, assembly		
	10	Hose SNAP-LOC. Complete		2NU90001
	20	14 guiding elements. SNAP-LOC, 2 pcs		2NU90002
16140036		Cutting Table f. X-table		
	30	Cutting table plate long		16140352
	190	Wire set for X-table		16143593
16140037 + +		Left Net, assembled + Right Net, assembled +		
16140042 +		Perforated sheet, assembled		
	20	Handle GN 565-20-100-SW		2GH22665
16140040		Side door, assembled		
	20	PETG Window Pane for side door		16140556
	50	Handle E+G M.643-200-B-M8		2GH00007
	120	Hinge GN 237-NI-50-50-A		2GG23750
16140041		Side door bottom labyrinth, assy		
	40	Sealing strip Black 3x15 4030, 2 pcs		2IP00302
	50	Sealing strip Black 3x15 4030		2IP00302
16140045		Front door, assembled		
	90	Support for window		16140573
	120	Roll with bearing ZABI R- 44		2BK50044
	210	Rubber stop		2GB00015
	330	Straight Actuator AZM161-B6, 2 pcs		2SS01616
16140050		Chamber Bellows, assembled		
	20	Bellows Hestego-Horizontal		16140535
	30	Bellows Hestego-Vertical Rear, 2 pcs		16140534
	40	Bellows Hestego-Vertical Front		16140533
	100	Roll with bearing ZABI R-33, 4 pcs		2BK50033
	130	Roll with bearing ZABI R-25, 4 pcs		2BK50026

	Spare part list for Magnutom				
Drawing	Pos.	Spare Part	Ref.	Cat no:	
16140051		Shield Capsule,Assembled			
	70	Roll with bearing ZABI R- 44		2BK50044	
	80	V-ring VA200, 4 pcs		2IV10200	
16140053		Upper Profile, Bellows, assemble			
16140013	70	Sealing rubber for Gap		16140736	
16140013	120	Seal. Rubber for Front Panel Gap		16140735	
16140012	170	Cam, Southco E5-6014		2GL35583	
16140012	190	Sliding bearing GFM-081017-15		2BG01715	
16140012	360	HT tube dia. 110		2NP00100	
16140058		Sample barrier Left, assembled			
	290	Wire set for X-table		16143593	
16140061 +		Protection Shield 1,assembled			
16140062		Protection Shield 2,assembled			
	40	Profile Gumex 0536013		2IP01065	
16140070		Control Panel, assembly			
	100	USB protective with link		2IB00105	
16140071		Front panel, assembly			
	30	PCB MultiPurpose,, Tested + Bootload.	A4	16013000	
	40	Turn - Push Button		15090600	
	50	Foil, Magnutom-500	SB2	16140725	
	70	PCB for SMM, Tested		15483004	
	110	Cap,PushButt.transp.5.5200.102		2SA30102	
	130	Emergency stop ø22 type RV		2SA10400	
	160	Turn/pushEncoderW.cable6inch24	BR1	2HR12413	
	180	Joystick 3 axis HFX-36P01	SQ1	2RQ00036	
	210	Push butt. Head RVAT DG Stainl.		2SA00400	
	220	Module holder. 3 elem. MHR-3, 2 pcs		2SA41603	
	250	Display 320x240 TFT-color, LED	A5	2HD50200	
	330	Display Flat Cable, 40p, 210mm		2WF01150	
16140075		HW box and parts, assembly			
	80	Hinge RT26-6		2GG02606	

	Spare part list for Magnutom			
Drawing	Pos.	Spare Part	Ref.	Cat no:
	90	Filter A1000-FIV3050-SE, (1 pc)	Z1	2MO93050
	100	Brake resistor CBV405-20R	R1	2RK05605
	110	AC reactor AX-RAI00360500-DE	L1	2MO60500
	120	Inverter V1000 VZA4015FAA	A1	2PU12050
	130	DC Link Reactor AX-RC01750430-DE	L2	2MO50430
	140	On/Off Switch 3x80A KG80	QS1	2SE20321
	150	Power Supply 88-264V, 24V/10A	A2	2PA90240
	160	Autotransformer ASUL120C	TM1	2MT12085
	180	Contactor J7KNA-AR-22-24D	K3+K4	2KM70900
	220	Standalone Safety Controller, (1 pc)	A3	2KS10010
	230	Omron J7KNG-14-01-24D Contactor	K5+K6	2KM71411
	260	Circuit breaker J7MN-3P-1E6, (1 pc)	FT1	2KM71805
	280	Omron J7KNG-10-10-24D contactor, (1pc)	K2	2KM71410
	290	Omron J7KNG-40-24D contactor	K1	2KM74010
	300	Solid state relay G3RV- SL500- D AC/DC24	KL1+KL2	2KL50024
	310	Lock Southco E3-55-85, 2 pcs		2GL35585
_	530	4.0 AT CICC fuse 10x38 CSA, (2pcs)	F1+F2	2FC11040
_	540	6.0 AT CICC fuse 10x38 CSA, (1pc)	F3	2FC11060
	620	Noise filter FN2412H-8-44	Z2	2MO12844
16140077		X-table Connection Box		
	190	3.5A 150V BYV28-150 FAST RECT., (3pcs)	VD1-VD3	2VD35028
16140088		Seal, assembled		
	20	Simmer ring SD 30x40x4		2IS30030
	30	V-ring V-50A		2IV10050
16140095		Plate for Valves, Assembled		
	40	Hose PVC w. steel coil d18-d25		2NU21226
	50	Press. hose. arm. PVC 3/8-ø10, 1 m		2NU29312
_	60	Hose PVC w. steel coil ø12/ø18, 1.8 m		2NU21218
	100	Solenoid valve 2/2, 24Vdc, 2 m	Y3	2YM10132
	230	Press. switch 0.5E BAR G¼in MS	BP1	2SP10011
	440	Valve ASCO 287	Y1	2YM10287
16140097		Water pistol		
	10	Cleaning Pistol		16140098
	30	Shower hose P04-HV-xxx-06		2NU99012

	Spare part list for Magnutom			
Drawing	Pos.	Spare Part	Ref.	Cat no:
		Loose parts, Magnutom-500		
		Wrench 36		2GR36410
		HT Elbow 110		2NG00110
		Key Southco E3-26-819-15		2GR81915
		Triangular key M5-100mm		2GR00887
		Wire Cut Control, Magnutom, item 16143585		
		Cable W18 Magnutom-500		16143518
		Cable W21 Magnutom-500		16143521
		Cable W22 Magnutom-500		16143522
		Cable W24 Magnutom-500		16143524
		Cable W27 Magnutom-500		16143527
		Cable W29 Magnutom-500		16143529
		Cable W38 Magnutom-500		16143538
		Earth wires Magnutom-500		16143554
		Wire Cut Magnutom, item 16143580		
		Cable W1 Magnutom-500		16143501
		Cable W2 Magnutom-500		16143502
		Cable W3 Magnutom-500		16143503
		Cable W5 Magnutom-500		16143505
		Cable W6 Magnutom-500		16143506
		Cable W7 Magnutom-500		16143507
		Cable W12 Magnutom-500		16143512
		Cable W13 Magnutom-500		16143513
		Cable W14 Magnutom-500		16143514
		Cable W15 Magnutom-500		16143515
		Cable W16 Magnutom-500		16143516
		Cable W25 Magnutom-500		16143525
		Cable W26 Magnutom-500		16143526
		Cable W32 Magnutom-500		16143532
		Cable W33 Magnutom-500		16143533
		Cable W34 Magnutom-500		16143534
		Cable W37 Magnutom-500		16143537
		Cable W39 Magnutom-500		16143539
		Cable W41 Magnutom-500		16143541
		Cable W42 Magnutom-500		16143542

Spare part list for Magnutom				
Drawing	Pos.	Spare Part	Ref.	Cat no:
		Wire Set, item 16143590		
		Cable W4 Magnutom-500		16143504
		Cable W9 Magnutom-500		16143509
		Cable W11 Magnutom-500		16143511
		Contact block 1 NC 1/2. MTO		2SB10071
		Contact block 1 NO 3/4. MTI		2SB10072
		Wire set for X-table, item 16143593		
		Actuator Linak_LA35	M4	2MK70001
		Wire set Control, Magnutom, item 16143595		
		Contact block 1 NC 1/2. MTO		2SB10071
		Contact block 1 NO 3/4. MTI		2SB10072
		Cable W10 Magnutom-500		16143510
		Cable W36 Magnutom-500		16143536
		Lamp HL1 w. plug	HL1	16143556
		Lamp HL2 w. plug	HL2	16143557
		Namur HQ1 w. plug	HQ1	16143558
		Strain Gauge YR1 w. plug	YR1	16143560
		Namur HQ2 w. plug	HQ2	16143562
		AZM 161SK-1212RKED-024, (2pcs)	YS1,YS2	2SS00120






























































	1	2	3			4	
A							
в	•	• •	•	**	(	20)	
с		10					
D		A		0			
Е					6X 30 Loctite 24	13	
	A 2012-11-30				2012 12 24	TDD	
	Rev Crea. date	Revision description		Draw Init	Appr. date	Appr Init	
	yyyy-mm-dd	→ ~ [Material:	Scolor	Eormat:	yyyy-mm-dd		
F	(€		1:5	A4	Veight ·	1350 4 a	
	Pederstrupvej 84	D: Description:	1.0		Treigint .	Rev:	
	DK-2750 Ballerup Copenhagen Denmark	16140053 Upper Prof	40053 Unner Profile Bellows assembled				
	Phone :+45 44600 800 Fax : +45 44600 804						

















































## **Overview, Variant Parts in Magnutom:**

Country Nom. voltage / freq.	Guard for Cut–off wheel (wheel diameter)	Table	Machine Item No .:	
	432 mm	Fixed	06276146	
3x380-480V		Automatic	06276246	
50-60Hz	508 mm	Fixed	06146146	
		Automatic	06146246	

