

# LaboForce-Mi Instruction Manual

**Original Instructions** 

CE

Doc. no.: 16387025-02\_B\_en Date of release: 2025.01.14

#### Copyright

The contents of this manual are the property of Struers ApS. Reproduction of any part of this manual without the written permission of Struers ApS is not allowed.

All rights reserved. © Struers ApS.

## **Table of Contents**

1	About this manual			
2	Safe	Safety		
	2.1	Intended use	5	
	2.2	LaboForce-Mi safety precautions	5	
		2.2.1 Read carefully before use	5	
	2.3	Safety messages	6	
	2.4	Safety messages in this manual	7	
3	Get	Get started		
	3.1	Device description	7	
	3.2	LaboForce-Mi	8	
4	Trai	nsport and storage	9	
	4.1	Long-term storage or shipping	9	
5	Inst	allation	9	
	5.1	Unpack the machine	9	
	5.2	Check the packing list	10	
	5.3	Installation - LaboForce-Mi	10	
	5.4	The specimen mover	11	
		5.4.1 Insert a specimen mover plate	11	
		5.4.2 Adjust the height of the specimen mover plate	12	
		5.4.3 Adjusting the horizontal position of the specimen mover plate	12	
		5.4.4 Complete the adjustment	12	
	5.5	Noise and vibration levels	13	
6	Оре	rate the device	13	
	6.1	Control panel functions	13	
7	Mai	Maintenance and service		
	7.1	General cleaning	14	
	7.2	Daily	14	
	7.3	Weekly	14	
	7.4	Spare parts	15	
	7.5	Service and repair	15	
	7.6	Disposal	15	
8		ubleshooting	15	
9		hnical data	16	
	9.1	Technical data	16	

Declaration of Incorporation of Partly Completed Machinery			21
10 Manufacturer			19
	9.3	Legal and regulatory information	19
	9.2	Diagrams	17

## **1** About this manual



CAUTION

Struers equipment must only be used in connection with and as described in the Instruction Manual supplied with the equipment.



Read the Instruction Manual carefully before use.



Note

Note

If you want to view specific information in detail, see the online version of this manual.

## 2 Safety

## 2.1 Intended use

The device is for professional semi-automatic materialographic preparation (grinding or polishing) of materials for further materialographic inspection.

LaboForce-Mi is intended to be used in combination with:

- LaboPol-30
- LaboPol-60

**Do not use the device for** Other equipment than stated in the section "Intended use". **the following** 

Model

LaboForce-Mi

## 2.2 LaboForce-Mi safety precautions

#### 2.2.1

#### Read carefully before use

Ignoring this information and mishandling of the equipment can lead to severe bodily injuries and material damage.

#### **General safety precautions**

- 1. The machine must be installed in compliance with local safety regulations. All functions on the machine and any connected equipment must be in working order.
- 2. This device must be mounted securely on the machine.

- 3. The operator must read the safety precautions and Instruction Manual, as well as relevant sections of the manuals for any connected equipment and accessories.
- 4. If you observe malfunctions or hear unusual noises, switch off the machine and call technical service.
- 5. In case of fire, alert bystanders and the fire brigade. Disconnect the electrical power supply. Use a powder fire extinguisher. Do not use water.
- 6. This machine must be operated and maintained only by skilled/trained personnel.
- 7. Always switch off the electrical power supply and remove the plug or power cable before dismantling the machine or installing additional components.
- 8. The machine must be disconnected from the electrical power supply before any service. Wait 5 minutes until residual potential on the capacitors is discharged.
- 9. Struers equipment must only be used in connection with and as described in the Instruction Manual supplied with the equipment.
- 10. If the equipment is subjected to misuse, incorrect installation, alteration, neglect, accident or incorrect repair, Struers will accept no responsibility for damage to the user or the equipment.
- 11. Dismantling of any part of the equipment, during service or repair, should always be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.)

### 2.3 Safety messages

Struers uses the following signs to indicate potential hazards.



#### **ELECTRICAL HAZARD**

This sign indicates an electrical hazard which, if not avoided, will result in death or serious injury.



#### DANGER

This sign indicates a hazard with a high level of risk which, if not avoided, will result in death or serious injury.



#### WARNING

This sign indicates a hazard with a medium level of risk which, if not avoided, could result in death or serious injury.



#### **CRUSHING HAZARD**

This sign indicates a crushing hazard which, if not avoided, could result in minor, moderate or serious injury.



#### HEAT HAZARD

This sign indicates a heat hazard which, if not avoided, can result in minor, moderate or serious injury.





Emergency stop Emergency stop

result in minor or moderate injury.

CAUTION

#### **General messages**

Note



This sign indicates that there is a risk of damage to property, or a need to proceed with special care.

This sign indicates a hazard with a low level of risk which, if not avoided, could



This sign indicates that additional information and hints are available.

### 2.4 Safety messages in this manual



#### WARNING

Struers equipment must only be used in connection with and as described in the Instruction Manual supplied with the equipment.



#### CAUTION

Keep clear of rotating parts during operation.



#### CAUTION

When working at machines with rotating parts, take care to prevent clothes and/or hair from being caught by the rotating parts.

## **3 Get started**

## 3.1 Device description

LaboForce-Mi is a control panel/specimen mover unit which can be mounted on the following machines for semi-automatic mineralographic preparation (grinding/polishing): LaboPol-30 and LaboPol-60.

The operator selects the grinding/polishing surface and the cooling fluid/abrasive suspension to be applied.

A holding device for the preparation of specimens must be used with the specimen mover. The holding device holds 8 specimens.

The specimen mover head has pressurized feet which can be adjusted individually to keep the specimens in place.

The operator decides which holding device should be used.

The operator sets the rotation speed of the preparation disc before starting the process.

The operator starts the machine by pressing the **Start** button on the control panel.

The operator applies other fluids manually or by using a separate dosing unit.

The operator stops the machine by pressing the Stop button on the control panel.

The operator cleans the specimens before the next preparation step or inspection.

The machine must always be used with the splash guard in place.

If the emergency stop is activated, the power to all moving parts is cut.


### 3.2 LaboForce-Mi

#### **Front view**




- A Control panel
- **B** Disc speed control
- **C** Specimen mover head
- **D** Control panel column

**Rear view** 



- **A** Rotation switch (Specimen mover head)
- B Steel pin
- **C** Locking handle
- **D** Name plate

#### The specimen mover



- A Mover cross for forced rotation
- **B** Force indicators
- **C** Force adjustment screw
- D Quick release button
- E Pressure feet
- F Specimen mover plate

## 4 Transport and storage

If, at any time after the installation, you have to move the unit or place it in storage, there is a number of guidelines we recommend that you follow.

- Package the unit securely before transportation. Insufficient packaging could cause damage to the unit and will void the warranty. Contact Struers Service.
- We recommend that you use the original packaging and fittings.

### 4.1 Long-term storage or shipping

For instructions on how to prepare the unit for long-term storage or shipping, see the specific machine manual.

Pay specific attention to the following:

- Disconnect the unit from the electrical power supply.
- Remove the control panel or specimen mover.
- Place the machine and accessories in their original packaging.

## 5 Installation

## 5.1 Unpack the machine

Note



We recommend that you keep all original packaging and fittings for future use.

- 1. Cut the packing tape on the top of the box.
- 2. Remove the loose parts.
- 3. Remove the unit from the box.

## 5.2 Check the packing list

Optional accessories may be included in the packing box.

The packing box contains the following items:

Pcs.	Description
1	LaboForce-Mi
1	Spacing disc
1	Instruction Manual set

## 5.3 Installation - LaboForce-Mi

Note



This device must be mounted securely on the machine.




**Note** Do not use the speed control knob on the control panel to move the specimen mover.


#### Procedure


Install the specimen mover in the connection hole of the machine.

- 1. Remove the plastic disc protecting the communication cable.
- 2. Connect the communication cable on the column to the connection port on the machine.




- 3. Slide the black V-ring down the column until it covers the connection hole.
- 4. Use the Allen key to tighten the two fastening screws. Do not tighten the screws completely.





## 5.4 The specimen mover



#### 5.4.1 Insert a specimen mover plate

- 1. Insert a specimen mover plate and push until the pin is aligned in the groove.
- 2. Make sure that the specimen mover plate is securely mounted.



A Mover cross for forced rotation


**B** Force adjustment screw

**D** Quick-release button

**F** Specimen mover plate

**C** Force indicator


**E** Pressure feet



- A Pin
- B Groove

#### 5.4.2 Adjust the height of the specimen mover plate

- 1. Use the locking handle on the left side to unlock the specimen mover and let it move into upright position.
- 2. Select the 'thickest' preparation surface to be used and place it on the preparation disc. Usually, this will be SiC Foil on an MD-Gekko disc, or SiC Paper on an MD-Fuga disc, or an MD-Alto.
- 3. Place the supplied spacing disc on the preparation surface.
- 4. Support the specimen mover head and loosen the 2 fastening screws that hold the column.
- 5. Lift and support the specimen mover.
- 6. Press the specimen mover head down as far as possible.
- 7. Use the locking handle to lock the specimen mover head in operating position.




- A Spacing disc
- 8. Lower the column until the specimen mover plate rests on the spacing disc.
- 9. Adjust the horizontal position of the specimen mover plate.

SeeAdjusting the horizontal position of the specimen mover plate ► 12.

#### 5.4.3 Adjusting the horizontal position of the specimen mover plate

- 1. Moves the specimen mover head to the right.
- Place the specimen mover plate in a position that allows the specimen to run 3 – 4 mm over the edge of the preparation disc.



A Disc edge

**B** Specimen mover plate

#### **5.4.4 Complete the adjustment**

- 1. Tighten the 2 fastening screws firmly. The specimen mover will now remain in position.
- 2. Cover the holes with the two caps.

Allen key and cover caps are included in the package.

### 5.5 Noise and vibration levels

For information about noise and vibration levels, see the specific machine manual.

## **6 Operate the device**

For instructions on how to operate the device, see the specific machine manual. See also the section "Intended use" for the specific machine manual.

## 6.1 Control panel functions



#### CAUTION

Keep clear of rotating parts during operation.



#### CAUTION

When working at machines with rotating parts, take care to prevent clothes and/or hair from being caught by the rotating parts.



A Disc speed control

#### **Disc rotation**



Starts rotation of the disc (**Spin** function).

#### Water


Manual override

- Press the button to apply water. Water is applied when no process is running.
- Press the button again to stop applying water.

#### Start

- $\diamondsuit$
- Starts the preparation process.

LaboForce-Mi



## 7 Maintenance and service

#### Technical questions and spare parts

If you have technical questions or when you order spare parts, state serial number and voltage/frequency. The serial number and the voltage are stated on the name plate of the machine.

### 7.1 General cleaning

To ensure a longer lifetime for your machine, we strongly recommend regular cleaning.



**Note** Do not use a dry cloth as the surfaces are not scratch resistant. Grease and oil can be removed with ethanol or isopropanol.



Note

Do not use acetone, benzol or similar solvents.

#### If the machine is not to be used for a longer period of time

• Clean the machine and all accessories thoroughly.

### 7.2 Daily

• Clean all accessible surfaces with a soft, damp cloth.

### 7.3 Weekly

• Clean the control panel with a soft damp cloth and common household detergents.

### 7.4 Spare parts

#### Technical questions and spare parts

If you have technical questions or when you order spare parts, state serial number and voltage/frequency. This information is stated on the name plate of the machine.

For further information, or to check the availability of spare parts, contact Struers Service. Contact information is available on <u>Struers.com</u>.

### 7.5 Service and repair

We recommend that a regular service check be carried out yearly or after every 1500 hours of use.

When the machine is started up, the display shows information about total operation time and the machines service information.

After 1500 hours of operation time, the display will show a message reminding the user that a service check should be scheduled.



Note

Service must only be performed by a qualified technician (electromechanical, electronic, mechanical, pneumatic, etc.). Contact Struers Service.

## 7.6 Disposal



Equipment marked with a WEEE symbol contains electrical and electronic components and must not be disposed of as general waste.

Contact your local authorities for information on the correct method of disposal in accordance with national legislation.

For disposal of consumables and recirculation fluid, follow local regulations.

## 8 Troubleshooting

For troubleshooting guidance, see the specific machine manual.

## 9 Technical data

## 9.1 Technical data


Specimen mover	Specimen holder	
	Force	2.5-20 N individual specimen
	Rotational speed	
	in process	3 rpm
	Rotational direction	Counter-clockwise, Clockwise
	Motor	15 W
	Torque	> 49 Nm
Safety standards/directives/legislation		See the Declaration of Conformity/Instruction manual
REACH		For information about REACH, contact your local Struers office
Operating environment	Surrounding temperature	_
	During operation	5-40 °C (41-104 °F)
	During transport	0-60 °C (32-104 °F)
	Humidity	< 85% RH non- condensing
Power supply	Voltage/frequency	1x24 V DC
	Power inlet	Connected directly to LaboPol
	Power	
	Nominal load	15 W
	Current	
	Nominal	0.6 A
	Max.	1.5 A
	Current, largest load	0.5 A

Width	18.3 cm (7.2")
Depth	43.5 cm (17.1")
Height	44.8 cm (17.6")
Height (with open cover)	measured from preparation disc
Weight	13 kg (29 lbs)
A-weighted sound emission pressure level at workstations	See the LaboPol technical data
	Depth Height Height (with open cover) Weight A-weighted sound emission pressure level at

## 9.2 Diagrams

Title	No.
LaboForce-50/LaboForce-Mi, Block diagram 1635305	

#### 16353050 B



## 9.3 Legal and regulatory information

#### FCC notice

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

## **10** Manufacturer

Struers ApS Pederstrupvej 84 DK-2750 Ballerup, Denmark Telephone: +45 44 600 800 Fax: +45 44 600 801 www.struers.com

#### **Responsibility of the manufacturer**

The following restrictions should be observed, as violation of the restrictions may cause cancellation of Struers legal obligations.

The manufacturer assumes no responsibility for errors in the text and/or illustrations in this manual. The information in this manual is subject to change without notice. The manual may mention accessories or parts not included in the supplied version of the equipment.

The manufacturer is to be considered responsible for effects on safety, reliability, and performance of the equipment only if the equipment is used, serviced, and maintained in accordance with the instructions for use.

## Declaration of Incorporation of Partly Completed Machinery

Manufacturer	Struers ApS • Pederstrupvej 84 • DK-2750 Ballerup • Denmark
Name	LaboForce-Mi
Model	N/A
Function	Specimen mover for LaboPol-30/60 Grinding/Polishing machine
Туре	638
Cat. no.	06386130

The above listed machinery is only intended to be used with Struers machinery and consumables, and shall not be put into service until the final machinery into which it is to be incorporated has been declared in conformity with this regulation, where appropriate.

Serial no.

CE

Struers Ensuring Certainty

#### Module H, according to global approach

We declare that the product mentioned is in conformity with the following legislation, directives and standards:

2006/42/EC	EN ISO 12100:2010, EN ISO 13849-1:2015, EN ISO 13849-2:2012, EN 60204-1:2018, EN 60204-1- 2018/Corr.:2020
2011/65/EU	EN 63000:2018
2014/30/EU	EN 61000-3-2:2014, EN 61000-3-3:2013, EN 61000-6-1:2007, EN 61000-6-3:2007, EN 61000-6-3- A1:2011, EN 61000-6-3-A1-AC:2012
Additional standards	NFPA 79, FCC 47 CFR Part 15 Subpart B

Authorized to compile technical file/ Authorized signatory Date: [Release date]

EU



- en For translations see
- bg За преводи вижте
- cs Překlady viz
- da Se oversættelser på
- de Übersetzungen finden Sie unter
- el Για μεταφράσεις, ανατρέξτε στη διεύθυνση
- es Para ver las traducciones consulte
- et Tõlked leiate aadressilt
- fi Katso käännökset osoitteesta
- fr Pour les traductions, voir
- hr Za prijevode idite na
- hu A fordítások itt érhetők el
- it Per le traduzioni consultare
- ja 翻訳については、
- lt Vertimai patalpinti
- lv Tulkojumus skatīt
- nl Voor vertalingen zie
- no For oversettelser se
- pl Aby znaleźć tłumaczenia, sprawdź
- pt Consulte as traduções disponíveis em
- ro Pentru traduceri, consultați
- se För översättningar besök
- sk Preklady sú dostupné na stránke
- sl Za prevode si oglejte
- tr Çeviriler için bkz
- zh 翻译见

www.struers.com/Library