



		LaboForce-Mi
Specimen mover	Specimen holder	
	Force	2.5-20 N individual specimen
	Rotational speed	
	in process	3 rpm
	Rotational direction	Counter-clockwise, Clockwise
	Motor	15 W
	Torque	> 49 Nm
Safety standards/directives/legislation		See the Declaration of Conformity/Instruction manual
REACH		For information about REACH, contact your local Struers office
Operating environment	Surrounding temperature	–
	During operation	5-40 °C (41-104 °F)
	During transport	0-60 °C (32-104 °F)
	Humidity	< 85% RH non-condensing
Power supply	Voltage/frequency	1x24 V DC
	Power inlet	Connected directly to LaboPol
	Power	
	Nominal load	15 W
	Current	
	Nominal	0.6 A
	Max.	1.5 A
	Current, largest load	0.5 A
Dimensions and weight	Width	18.3 cm (7.2")
	Depth	43.5 cm (17.1")
	Height	44.8 cm (17.6") measured from preparation disc
	Height (with open cover)	
Weight	13 kg (29 lbs)	
Noise level	A-weighted sound emission pressure level at workstations	See the LaboPol technical data